# splat CDWA3

## Cleaning indicator

For cleaning performance tests.

## Usage

Designed for routine control of cleaning processes and for monitoring automatic washing machines effectiveness as well as for monitoring ultrasonic washing machines which operate at frequencies equal to or greater than 35 kHz, and temperatures between 30 °C and 70 °C.

## Applicable Regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN ISO 13485:2016. ISO/TS 15883-5:2005.

## Characteristics

Challenge level: high.

Allows detection of all factors affecting the outcome of cleaning process.

Initial Color: blue.

Size: 70 mm x 30 mm (approximate).

Synthetic temperature-stable support.

Result: elimination of blue organic substance.

To interpret results visually, use the results reference guide on product's instructions for use.

Optionally, the readout of the cleaning indicator can be done using Trazanto Automatic Quality Control and Traceability System (see instructions for use in device's manual).

NOTE: If unfavorable results are obtained, check quality/quantity of detergent, cleaning time and temperature, positioning of the load, cleaning indicator location and general condition of the washing machine.

## Environmental conditions during manufacture

T = 15-30 °C, RH = 30-80 %.

## Storage conditions

T = 10-30 °C, RH = 30-80 %, store in dark place within its original package.

## Transportation conditions

Storage conditions should be strictly followed.

Transport in closed and strengthened boxes in order to avoid damages. Transportation of this product does not represent a risk for health.

## Shelf-life

5 years.

**Packing**

200 units per aluminum foil bag.

Packing information: product description, process for intended use, reference to regulation, instructions for use, storage conditions, manufacturer information and data on label.

**Labelling**

On product's label: product's code and description, product's presentation, batch number manufacture and expiration date, bar code and datamatrix code.

**Possible target markets**

Healthcare, Food, Pharmaceutical and Medical industries.

**Other important information**

Read product's instructions for use thoroughly before use.

**Description of identified and/or predictable risks which could represent a risk for the product**

Store the product within its original package, away from light, disinfectants, gases and humidity.

Do not re-use the product.

Do not touch indicator ink with your hands.

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