




**What's New?**

**02/17**

Stackable Step Stools - p. 22

Anti-Fatigue Step stool mat p.21 & 24

Stainless Steel Micro Brush - p. 49




**BATRIK**

●● Medical Manufacturing

**Catalogue 2017**



**Interactive PDF Tip:**

When on a Product page, click on these tabs  to go back to the Table of Contents at any time!

**Surgical Products**

Suture Tags	06
Clamp Jaw Covers	07
Vascular Silicone Ties	08
GOLFF - Anti-Fog Solution	09
Cautery Tip Cleaner	10
SKINS marker	11

**Surgical & General Accessories**

Pin Guards	14
Safe-Soak Mats & SPAC	15
Magnetic Mat	16
Tourniquets	17
Blue Cloud-Gel	18
Blue Cloud-Nimbus	20
<b>NEW!</b> Stackable Step Stools	22

**Endoscopy Products**

Bite Blocks	26
Safe-Soak Endo Sponge	28
Valve Covers	29
Channel Brushes	30
Valve/Control Brushes	33
External Cleaning Brushes	34

**Sterilization Products**

Silicone Pin Mat	36
Corner Pockets	37
Instrument Protectors	38
STICS & ID tape	40
Padlocks	41
Face shield	42
Scrub Pads	43

**Brushes**STERILIZATION  
BRUSHES


ARC Wire Brushes	46
Stainless Steel Micro Brush	49 <b>NEW!</b>
ARC Double Wire Brushes	50
Antimicrobial Twisted Wire Brushes	52

INSTRUMENT  
CLEANING BRUSHES

Flat Brushes	53
Precision Brush	54
Jar Brush	55
Acetabular Brush	56
All Purpose Brushes	57
Autoclave Cleaning Brush Kit	58









**BATRIK Suture Clamp Tags** are designed to protect and safely hold sutures without damage

- packed in a convenient foam base with easy peel-off (even with gloves) adhesive backing
- the foam base is affixed to a sterile area with the clamps presented in front of the Clamp Tags to allow for easy insertion
- available in **sterile** and **non sterile**
- available in 2, 3 or 5 pairs
- radio opaque
- non-toxic
- pre-counted
- latex free



**Easy Peel Corner**



### Suture Clamp Tags

Sterile item #	Non Sterile item #	Packaging	Size	Color
11-1002	911-1002	2 pairs / foam base	Standard $\varnothing$ 0.06" x 0.5" L	Yellow
11-1003	911-1003	3 pairs / foam base	Standard $\varnothing$ 0.06" x 0.5" L	Yellow
11-1005	911-1005	5 pairs / foam base	Standard $\varnothing$ 0.06" x 0.5" L	Yellow
11-1025	911-1025	5 pairs / foam base	Standard $\varnothing$ 0.06" x 0.5" L	Assorted
11-1045	911-1045	5 pairs / foam base	Standard $\varnothing$ 0.06" x 0.5" L	White
11-1065	911-1065	5 pairs / foam base	Standard $\varnothing$ 0.06" x 0.5" L	Blue
11-1085	911-1085	5 pairs / foam base	Standard $\varnothing$ 0.06" x 0.5" L	Red
	911-3002*	2 pairs / foam base	Large $\varnothing$ 0.09" x 0.75" L	Pink
	911-3003*	3 pairs / foam base	Large $\varnothing$ 0.09" x 0.75" L	Pink
	911-3005*	5 pairs / foam base	Large $\varnothing$ 0.09" x 0.75" L	Pink

**Our "Easy Peel Corner" is a feature unique to BATRIK**

**Packaging:** Sterile: 1 foam base/pouch – 5 pouches/box  
Non Sterile: sold per foam base



**\*Only available in a non sterile version. Please Call for more details.**

www.surgmed.com • info@surgmed.com • +1 514-631-7988 • Fax +1 514-631-9083

**BATRIK**  
 ●●● Medical Manufacturing

**BATRIK Clamp Jaw Covers** reduce trauma to vessels requiring occlusion, particularly when a hemostatic clamp is utilized. The thread-like indents allow for gentle, constant and non-slip tension.

- manufactured from Medical-Grade silicone
- gamma sterilized
- available in 3 colors and 3 widths to fit most clamp sizes
- convert almost any instrument to atraumatic
- radio opaque
- latex free



**Clamp Jaw Covers**

Sterile Item #	Non Sterile Item #	Dimensions – inches / cm	Color
12-1000	912-1000	0.1" x 5" / 0.25 x 12.7cm	Red
12-1001	912-1001	0.2" x 5" / 0.51 x 12.7cm	White
12-1002	912-1002	0.3" x 5" / 0.76 x 12.7cm	Blue

**Packaging:** Sterile: 2 clamps per sterile pouch – 10 pouches per box  
Non Sterile: 20 clamps per bag

**BATRIK Vascular Silicone Ties** permit safe and smooth retraction, occlusion and identification of vessels during surgical procedures.

Unlike clamps and retractors which can pinch the vessels flat, or even crack them, the BATRIK elliptical shape allows for the gentle handling of arteries, veins, nerves, tendons, spermatic cords, ureters, or whenever delicate retraction / occlusion is required.

- manufactured of Medical-Grade silicone
- do not rotate or twist on themselves when unpacked
- easy to separate — eliminating the need for cardboard inserts
- radio opaque
- available in 3 colors and 4 sizes
- latex free



### Vascular Silicone Ties

Item #	Dimensions – inches / cm	Color
10-1000	0.075" x 16" / 0.2cm x 40cm	Mini-Red
10-1001	0.075" x 16" / 0.2cm x 40cm	Mini-Blue
10-1002	0.075" x 16" / 0.2cm x 40cm	Mini-Yellow
10-1003	0.075" x 16" / 0.2cm x 40cm	Mini-White
10-1010	0.120" x 16" / 0.3cm x 40cm	Maxi-Red
10-1011	0.120" x 16" / 0.3cm x 40cm	Maxi-Blue
10-1012	0.120" x 16" / 0.3cm x 40cm	Maxi-Yellow
10-1013	0.120" x 16" / 0.3cm x 40cm	Maxi-White
10-1101	0.040" x 16" / 0.1cm x 40cm	Extra Mini-Blue
10-1111	0.200" x 22" / 0.5cm x 56cm	Extra Maxi-Blue
10-1113	0.200" x 22" / 0.5cm x 56cm	Extra Maxi-White

**Packaging:** 2 ties per sterile pouch — 10 pouches per box



**BATRIK GOLFF™ Anti Fog Solution** Laparoscopic surgical procedures often have problems with fogging. The BATRIK Sterile Anti-fog Solution, GOLFF™, is the solution that sets the standard for all other products.

- ensures clearer vision and avoids unnecessary interruptions during surgery
- prevents the occurrence of fog on lenses
- adhesive radiopaque sponge for secure placement
- easy squeeze 5ml bottle
- latex free
- OEM (available)



**GOLFF™**

Item #	Description	Packaging
GOLFF	Adhesive sponge and anti-fog per set	20 sterile sets / box

**BATRIK Cautery Tip Cleaners** assist in removing material from cautery tips for both mono-polar and bi-polar.

- foam pad with integrated abrasive material
- adheres to disposable and non-disposable drapes
- only available non-sterile
- radio opaque
- latex free




### Cautery Tip Cleaners

Item #	Description	Packaging
60-10000	Non Sterile — Cautery tip cleaners	2200 / box

**BATRIK SKINS™ – Surgical Skin Markers** Designed with a beveled tip to enable the creation of multiple tracing lines (fine or bold) just by holding the pen at different angles.

- made of lasting gentian violet ink that is non-toxic and non-tattooing
- available in non-sterile for pre-operative use, or sterile for use in the operating room
- available with or without 6 labels and a flexible 6" calibrated (in/mm) ruler
- surgical prep resistant
- latex free



**SKINS™ Surgical Skin Markers**

Item #	Description	Packaging
S-3000	STERILE – Multi-Tip fine marker with ruler and label set	50 / box
S-3001	STERILE – Multi-Tip fine marker	50 / box
S-30005	STERILE – Multi-Tip fine marker with ruler and label set (Volume Packaging)	500 / box
S-30015	STERILE – Multi-Tip fine marker (Volume Packaging)	500 / box
S-93000	NON STERILE – Multi-Tip fine marker with ruler and label set	1000 / box
S-93001	NON STERILE – Multi-Tip fine marker	1000 / box
S-94000	NON-STERILE – Ruler	100 / bag
S-94001	NON-STERILE – Ruler (Bulk)	20,000 / box
S-95000	NON-STERILE – Set of 3 labels	100 / bag

**SAVE** by ordering our SKINS MARKERS VOLUME PACKAGING – more cost effective and reduces waste





**BATRIK Pin Guards** You will no longer have to use tape or tubing as our BATRIK Pin Guards will safely cover the sharp ends of your K-Wire and Steinmann Pins.

- color coded by wire size which eliminates the selection of the wrong size
- a removable sleeve is provided to avoid injury and facilitate application (sterile only)
- packaging designed to reduce waste
- non sterile (100 units/ bag) also available
- latex free



### Pin Guards

Sterile item #	Non Sterile item #	Color	Dimensions inches / mm	Sterile Packaging	Non-Sterile Packaging**
14-1028	914-1028	Yellow	0.028" / 0.71mm	80 units/box*	100 units/box
14-1035	914-1035	Light Blue	0.035" / 0.89mm	80 units/box*	100 units/box
14-1045	914-1045	White	0.045" / 1.14mm	80 units/box*	100 units/box
14-1054	914-1054	Dark Red	0.054" / 1.40mm	80 units/box*	100 units/box
14-1062	914-1062	Green	0.062" / 1.57mm	80 units/box*	100 units/box
14-1078	914-1078	Pink	0.078" / 1.98mm	80 units/box*	100 units/box
14-1093	914-1093	Grey	0.093" / 2.36mm	80 units/box*	100 units/box
14-1109	914-1109	Beige	0.109" / 2.77mm	80 units/box*	100 units/box
14-1125	914-1125	Purple	0.125" / 3.17mm	80 units/box*	100 units/box
14-1028-1	N/A	Yellow	0.028" / 0.71mm	40 units/box	N/A
14-1035-1	N/A	Light Blue	0.035" / 0.89mm	40 units/box	N/A
14-1045-1	N/A	White	0.045" / 1.14mm	40 units/box	N/A
14-1062-1	N/A	Green	0.062" / 1.57mm	40 units/box	N/A
14-1141	914-1141	Brown	0.141" / 3.58mm	80 units/box*	100 units/box
14-1156	914-1156	Dark Blue	0.156" / 3.96mm	80 units/box*	100 units/box
14-1188	914-1188	Red	0.188" / 4.77mm	80 units/box*	100 units/box

\*Sterile Packaging: 2 units per sterile pouch; 40 pouches per box.

\*\* Non-Sterile: Does not include plastic sleeve

**BATRIK Safe-Soak™ Absorbent Floor Mats** are designed to absorb fluids from the operating room floor quickly and efficiently, to ensure that the floor beneath remains slip-free, safe and dry during surgical procedures. A skid-resistant backing is available with the Blue Absorbent Mat. All mats are single use and non-sterile.

- absorbs a variety of fluids
- provides quick and easy clean-up
- highly absorbent
- latex free, disposable, and non-sterile
- conforms to uneven surfaces



**Safe-Soak™ Absorbent Floor Mats**

Item #	Description	Size	Packaging	Absorbency
SS-20-1000	Safe-Soak™ / White extra absorbent mat	28" x 48" (71 cm x 122cm)	20 per box	4.5 litres / mat
SS-20-2000	Safe-Soak™ / Blue absorbent mat with skid-resistant backer	28" x 40" (71 cm x 102cm)	20 per box	2.5 litres / mat
SS-20-3000	Safe-Soak™ / Yellow absorbent mat	30" x 52" (76cm x 132cm)	15 per box	3.0 litres / mat

**BATRIK Safe-Soak™ SPAC** floor suction disk, quickly removes fluid waste from the operating room floor through the application of suction. These suction disks have a highly absorbent base material which allows quick fluid wicking and removal when suction is applied. Safe-Soak™ suction disks are latex free and non sterile.

- Easily maneuverable around operating room floor
- fluid is captured from bottom through a highly absorbent material
- opening permits fluid collection that come into contact with its surface
- highly absorbent and fast wicking
- latex free, disposable, and non-sterile
- highly durable and won't be damaged if stepped on
- the opening allows for easy movement using your foot
- the hose is securely affixed on its side to reduce risk of catching equipment



**Safe-Soak™ SPAC**

Item #	Description	Size	Packaging
SS-21-1000	Safe-Soak™ SPAC	8.625" / 22cm diameter	12 / box




**BATRIK Magnetic Mats** provide a **sterile, neutral field** for the safe placement and **hands-free transfer** of instruments.

- **raised magnets** hold instruments firmly in place while allowing easy grasping of instruments
- **round** magnets to avoid silicone rubber perforations
- **flexible** and **non-slip** backing drape conforms to any surface shape

**Reusable and Autoclavable (up to 134°C)**

**Latex free**



### Magnetic Mats

Item #	Dimensions – inches / cm	Packaging
MAG-02-1	12" x 15" / 30cm x 38cm	1 unit
MAG-02-5	12" x 15" / 30cm x 38cm	5 units



**BATRIK Tourniquets** are comfortable and easy to use on all patients, including; the elderly, undernourished, children and newborns.


- safeguards hospital staff and patients from latex allergens
- **new design:** Our material will feel like latex, but has a more comfortable texture for the patients and less stress on the skin
- available in different sizes
- compact sizes is convenient for most phlebotomy trays
- latex free

**Single Use Offers the Following Benefits :**

- most tourniquets are impossible to disinfect and re-usable units are often contaminated
- will reduce the risk of cross contamination (patient to patient or patient to healthcare provider)
- may help reduce the spread of nosocomial infections
- may improve your facility’s compliance with OSHA’s 1910-1030 (d)(1) “Blood-borne Pathogens Standard”
- 5-6% of the North American population have a latex allergy, particularly children with spina bifida



K-002 or K-004



K400PB



K-008 or K-009



K-010



K-400 or K-2-18RL

Facilitate single-use policy with our **25 units per roll**, which comes in a convenient dispenser box for easy access!



K-100 (or order by color)

**TOURNY™ + Tourniquets**


Item #	Description	Color	Dimensions	Packaging
K-002	TOURNY™ + Non Sterile	Blue	1" x 18" / 2.5cm x 45cm	1 unit / bag
K-004	TOURNY™ + Non Sterile	Blue	1" x 18" / 2.5cm x 45cm	100 units / box
K-008	TOURNY™ + Sterile	Blue	0.200" x 22" / 0.5cm x 56cm	20 units / box
K-009	TOURNY™ + Non Sterile	Blue	0.200" x 22" / 0.5cm x 56cm	20 units / box
K-010	TOURNY™ + Non Sterile	White	0.200" x 22" / 0.5cm x 56cm	20 units / box
K-400	TOURNY™ + Non Sterile	Blue	1" x 18" / 2.5cm x 45cm	16 rolls - 400 units/box
K-400PB	TOURNY™ + Non Sterile (Pebbled)	Light Blue	1" x 18" / 2.5cm x 45cm	16 rolls - 400 units/box
K-100	TOURNY™ + Non Sterile	Assorted	1" x 18" / 2.5cm x 45cm (3 rolls of each Blue, Orange, Green and Pink)	12 rolls - 300 units/box
K-254	TOURNY™ + Non Sterile (Pebbled)	Light Blue	1" x 18" / 2.5cm x 45cm	250 units/box
K-2-18RL	TOURNY™ + Non Sterile	Blue	1" x 18" / 2.5cm x 45cm	1 roll
K-2-18PBRL	TOURNY™ + Non Sterile (Pebbled)	Light Blue	1" x 18" / 2.5cm x 45cm	1 roll
K-3-18RL	TOURNY™ + Non Sterile	Green	1" x 18" / 2.5cm x 45cm	1 roll
K-6-18RL	TOURNY™ + Non Sterile	Pink	1" x 18" / 2.5cm x 45cm	1 roll
K-7-18RL	TOURNY™ + Non Sterile	Orange	1" x 18" / 2.5cm x 45cm	1 roll
K-2-29RL	Long TOURNY™ + Non Sterile	Blue	1" x 29.5" / 2.5cm x 75cm	1 roll



**BATRIK Medical Blue Cloud-Gel™ Anti-Fatigue Mat** is truly the medical industry standard. It has been scientifically engineered using patented gel technology combined with premium foam to provide the ultimate in comfort and support to surgeons, nurses and technicians.

Both the fabrication and material composition is specifically designed for use in all hospital applications, particularly where movement is restricted for prolonged durations. Our material composition will withstand the repeated exposure to a wide range of fluids.

Feel the **Blue Cloud-Gel™** Advantage  
designed for the O.R. !



Both top and bottom are stain resistant /  
anti-microbial & anti-bacterial

Integrated **Blue Cloud-Gel™**

Superior Quality Foam

Anti-Slip Backing



Latex-Free

Cart Washer Compatible (Contact Us  
for more details)

## COMPARE THE BENEFITS AND FEEL THE DIFFERENCE.

Integrated **Blue Cloud-Gel™** will not "Bottom-Out" like conventional foam mats.

Properly conforms and supports feet while promoting efficient circulation, reducing physical discomfort and fatigue.

Promotes proper posture while diminishing back, knee and foot discomfort.

Improves comfort and ergonomics while also enhancing energy.

Fabricated with the high standard for the operating room, but the anti-fatigue properties will benefit every department.


Studies have demonstrated that anti-fatigue mats improve performance, energy and post-operative discomfort.

[www.BlueCloudGel.com](http://www.BlueCloudGel.com)



**FEATURES INCLUDE:**

- Unique integrated **Blue Cloud-Gel™** technology provides superior long lasting anti-fatigue advantages over conventional foam based products.
- Incorporated high performance foam.
- Both the top and bottom are stain resistant and anti-microbial / anti-bacterial.
- The bottom material is Non-Slip Certified by the National Floor Safety Institute.
- Proprietary 3/4" Batrik **Blue Cloud-Gel™** Foam technology.
- Excellent friction properties when damp or dry surgical booties are worn.
- Water-resistant and non-absorbent.
- Easy to clean with most hospital quaternary cleaners.\*
- US standards: CAL 117 and NFPA 260 Compliant.
- Reusable
- 1 Year Warranty



**Blue Cloud-Gel Anti-Fatigue Mat**

Item #	Dimensions – inches / cm	Packaging
800-2401	18 x 24 / 46 x 61	1 / Box
800-2405	18 x 24 / 46 x 61	5 / Box
800-3201	20 x 32 / 51 x 81	1 / Box
800-3203	20 x 32 / 51 x 81	3 / Box
800-4801	20 x 48 / 51 x 122	1 / Box
800-4802	20 x 48 / 51 x 122	2 / Box
800-7201	30 x 72 / 76 x 183 (special order)	1 / Box


\*With the wide array of cleaners, we recommend that you test your solution beforehand.

**Did You Know?** We also offer the excellent anti-fatigue mat **"Nimbus"**, using only our superior foam based material (NOT for O.R. Application). You will be amazed with the benefits and our comfort technology.

[www.BlueCloudGel.com](http://www.BlueCloudGel.com)



**BATRIK Medical's Blue Cloud-Nimbus™ Anti-Fatigue Mat** proprietary high-density "energy-return" polyurethane materials will provide unmatched comfort, relief and ergonomic support.



### COMPARE THE BENEFITS AND SEE THE DIFFERENCE

- Our Energy-Return properties will help reduce leg, foot, lower back discomfort and work-related injuries.
- Feel the difference compared to any other foam or rubber based anti-fatigue mat on the market as our bounce-back memory restores energy.
- Longer lasting support as 3/4" (1.9 cm) proprietary foam blend will not bottom out like most other foam based anti-fatigue mats.
- Eco-Friendly materials are used, so the **Blue Cloud-Nimbus™** is not just good for you, but also good for the environment.
- If within 30 days you are not 100% satisfied with the added benefits, contact your dealer for a refund.\*



[www.BlueCloudGel.com](http://www.BlueCloudGel.com)

## FEATURES INCLUDE

- Anti-Microbial properties blended into the mat to provide added protection.
- Highly Water-Resistant Surface will not absorb water or dirt. (Not commonly found on other typical foam based mats.)
- One piece construction ensures longer useful life.
- Textured design on surface and high-traction base ensures greater safety. (US Standard NFSI 101-C)
- Durable design is ideal for commercial applications, yet still ergonomically & aesthetically pleasing for home or office environments.
- Flammability/Burn resistant per CFR 16 Parts 1630/1631
- Non-Toxic and Latex Free
- Light weight and easy to clean.\*\*
- 5 Year Warranty.\*\*\*



### Blue Cloud Nimbus Anti-Fatigue Mat

Item #	Dimensions – inches / cm	Packaging	Color	Item #	Dimensions – inches / cm	Packaging	Color
809-3201	20 x 32 / 51 x 81	1 / Box	Black	810-3201	20 x 32 / 51 x 81	1 / Box	Taupe
809-3204	20 x 32 / 51 x 81	4 / Box	Black	810-3204	20 x 32 / 51 x 81	4 / Box	Taupe
809-4801	20 x 48 / 51 x 122	1 / Box	Black	810-4801	20 x 48 / 51 x 122	1 / Box	Taupe
809-4804	20 x 48 / 51 x 122	4 / Box	Black	810-4804	20 x 48 / 51 x 122	4 / Box	Taupe
809-3601	24 x 36 / 61 x 91	1 / Box	Black	810-3601	24 x 36 / 61 x 91	1 / Box	Taupe
809-7201	20 x 72 / 51 x 183	1 / Box	Black	810-7201	20 x 72 / 51 x 183	1 / Box	Taupe
809-6001	36 x 60 / 91 x 152	1 / Box	Black	810-6001	36 x 60 / 91 x 152	1 / Box	Taupe

### NEW PRODUCTS

809-1705 Reusable Anti-Fatigue Mat for our step stool OR-36363 (p.24)<sup>1</sup>

899-1710 Single-Use Anti-Fatigue Mat for our step stool OR-36363 (p.24)<sup>1</sup>

<sup>1</sup>With adhesive strip to secure it in place.

\* Refund applies only to the purchase price, excludes any transport, handling or taxes.

\*\* We recommend using most soap based products. Do not submerge the mat in water.

\*\*\* 5 Year Warranty for Non-Industrial usage; 4 Year Warranty for Industrial usage.

## Stacking Step Stools by IMPERIAL SURGICAL

**BATRIK** is proud to be distributing on behalf of IMPERIAL SURGICAL the unique, antimicrobial Stackable Step Stool. Made of Aluminium, it has truly been designed with Infection Prevention as a prime focus. It is Cart Washable and can even be sterilized.\*

The IMPERIAL SURGICAL Stackable Stepstools are designed for the O.R., but are equally useful hospital wide. Made of light weight aluminum, our Stackable Stepstools are protected with **an anti-microbial powder coating**. Furthermore, they can support more weight than most of our competitors (see Figure 1).

- **Improved model**
- **Anti-slip surface** for better stability
- **No rubber corners:** this was a source of previous fluid collection (see the picture on next page, circled in red)
- **No Welds** or large Hollow opening to trap debris or fluids
- **Transport Cart** Available (Figure 2)
- **NEW:** You may now purchase a fully autoclavable stepstool without the mat! See last page for order details.

**NEW!**



\* The OR-36364 can be sterilized. See last page for details.

Can support the weight of half a Jeep Wrangler!



Figure 1



Figure 2

**Transport Cart**  
(see last page for order details)

[www.StackingStepStool.com](http://www.StackingStepStool.com)

## CAN THE COMPETITION COMPARE?

### **Anti-Microbial Coating**


Patent pending technology. These stools are often in contact with moisture, fluid and blood; and will benefit from anti-microbial protection. The benefits include:

- Reduced risks of infections<sup>1</sup>
- Prevents the spread of dangerous micro-organisms and suppresses the growth of mold, bacteria, yeast and mildew
- Provides superior durability and protection for essentially the life of the product.<sup>2</sup>

### **No Corner Rubbers**

Did you know that the corner rubbers found on some of the competition stackable stepstools can become a source for bacteria? With time, they erode and create deep cavity channels that can collect fluids and thus, create an environment for bacteria to propagate.

IMPERIAL SURGICAL has removed these rubber corners for your safety.



### **Compatible**

If you have older models or similar looking competitor stepstools, these newly designed models will stack and inter-lock perfectly.

### **High density foot pegs**

Our foot pegs are made of denser material than the competition, and will not break or come loose easily. And if they do; we will replace them free of charge! Only Transport fees will apply (see Warranty).

### **Now Available**

Reusable and Single-use (O.R.)\* Anti-Fatigue Mats. They provide comfort, support and reduce leg & back fatigue when standing for extended period.

\*Single-Use Replacement Mats to reduce cross contamination when needed.

### **7 years warranty**

We offer one of the best warranty in the market. Our 7 years warranty covers all parts (step stool only).

### **World-Class Customer Service**

Your satisfaction is our utmost priority.

## ORDERING INFORMATION

**Material:** Cast Aluminium


**Weight:** 8.3 lbs (3.8 kg)

**Weight Capacity:** 300 lbs (136 kg)

**Dimensions:** 17.4" x 13" x 4.7" (L x W x H)

**Packaging (1):** 19" x 14" x 6" (10 lbs)

**Packaging (4):** 26.5" x 20" x 15" (45 lbs)



OR-36363

### Stacking Step Stools

Item #	Description	Color
OR-36363	Stackable Stepstool with mat	Black
OR-36363-02	Stackable Stepstool with mat	Dark Blue
OR-36363-021	Stackable Stepstool with mat	Light Blue
OR-36363-03	Stackable Stepstool with mat	Green
OR-36363-07	Stackable Stepstool with mat	Orange
OR-36363-11	Stackable Stepstool with mat	Dark Grey
OR-36364	Stackable Stepstool <b>without</b> mat	N/A

**Autoclavable**



OR-3693

### Optional Accessory

Item #	Description	Packaging
OR-3693	Transport Cart	1
809-1705	Reusable Anti-Fatigue Mat*	5 units / box
899-1710	Single-Use Anti-Fatigue Mat*	10 units / box

\*With adhesive strips to secure it in place.



899-1710

### Replacement Parts

Item #	Description
OR-36363-1A	Replacement Pad for OR-36363
OR-36363-1B	Replacement Foot Pegs

Visit us at [www.StackingStepStool.com](http://www.StackingStepStool.com) or  
Write us at [info@surgmed.com](mailto:info@surgmed.com)  
for more product details





**BATRIK Endoscopy BO2™ Bite Blocks with Nasal Oxygen Delivery**

significantly enhances patient safety by providing supplemental oxygenation during gastrointestinal endoscopy preventing hypoxia in sedated patients, while ensuring superior patient comfort.

- lower procedural costs by obviating the need for separate nasal prongs and oxygen masks, the BO2 Bite Block helps save time and money
- ensure patient comfort and safety as with other Bite Blocks from BATRIK Medical, the BO2 Bite Block is designed to be free of sharp edges to ensure patient comfort and safety. Additionally, our bite blocks are guaranteed to be latex-free, lead-free and DEHP-free.
- choose from disposable and reusable options. The BO2 Bite Block is available in single use and reusable options. Use our flexible ARC Wire Brush #41-32505, which will allow you to properly clean the nasal prongs and inlet of the reusable BO2.
- ensure secure fit. An optional retention strap is available to ensure a secure fit. These straps, like the bite block, are latex free, lead free and DEHP free.



**BO2™ Bite Block**

Item #	Description	Packaging
BO2	18mm / 54FR - Single use, includes strap	100 / Box
BO2-N	18mm / 54FR - Single use, without strap	100 / Box
BO2-RN	18mm / 54FR - Reusable, without strap	25 / Box
BO2-R	18mm / 54FR - Reusable, includes strap	25 / Box

**BATRIK Endoscopy Bitestrong™ Bite Blocks**

are designed to provide a safe passage through the patient's mouth during TEE thus protecting the patient's teeth as well as the endoscope. Strict quality control procedures and the careful selection of raw materials used in the manufacture of our Bite Blocks is your guarantee that they are latex free, lead free and DEHP free.

**Single-Use**



**Bite Blocks**

Item #	Description	Packaging
B4	20mm / 60FR - Single use	100 / Box
B4T	20mm / 60FR - Single use, includes strap	100 / Box
B4T-N	20mm / 60FR - Single use, without strap	100 / Box
B1-25	16mm / 48FR - <b>Pediatric</b> single use, includes strap	25 / Box
B1	16mm / 48FR - <b>Pediatric</b> single use, includes strap	100 / Box
B1-N	16mm / 48FR - <b>Pediatric</b> single use, without strap	100 / Box

**BEST SELLER**

**BATRIK Endoscopy Bitestrong™ Single Use Bite Blocks**

- our molding process ensures a smooth texture finish with no sharp edges (flashing) that can often irritate a patient’s mouth
- our bite blocks are individually wrapped in sealed, easy-to-open PE bags
- available without dental retention
- additional straps are also available, if required



**Bitestrong™ Single Use Adult Bite Block**

Item #	Description	Packaging
B3	20mm / 60FR - Single use, includes strap	100 / Box
B3-N	20mm / 60FR - Single use, without strap	100 / Box
B3-NDR	Same as B3, no dental retention rim	100 / Box
B3-NDR-N	Same as B3-N, no dental retention rim, without strap	100 / Box

**BATRIK Endoscopy Bitestrong™ Reusable Bite Blocks**


- our reusable Bite Blocks are designed for repeated usage and will withstand commonly practiced chemical and steam sterilization cycles (137°C / 279°F)<sup>(1)</sup>
- our molding process ensures a smooth textured finish with no sharp edges (flashing) that can often irritate a patient’s mouth
- our Bite Blocks are individually wrapped in sealed, easy-to-open PE bags
- available with or without an adjustable retention strap
- additional straps are also available, if required



**Bitestrong™ Reusable Adult Bite Block and Straps**

Item #	Description	Packaging
B5	20mm / 60FR Reusable, includes strap	100 / Box
B5-N	20mm / 60FR Reusable, without strap	100 / Box
B5-25	20mm / 60FR Reusable, includes strap	25 / Box
B5-N-25	20mm / 60FR Reusable, without strap	25 / Box
B6	20mm / 60FR Reusable	100 / Box
B6T	20mm / 60FR Reusable, includes strap	100 / Box
B6T-N	20mm / 60FR Reusable, without strap	100 / Box
B6-25	20mm / 60FR Reusable	25 / Box
B6T-25	20mm / 60FR Reusable, includes strap	25 / Box
B6T-N-25	20mm / 60FR Reusable, without strap	25 / Box
BS	55cm x 20cm / 21.7" x 7.9" Single use Bite Block straps	100 / Box

Order your single use Bite Block sample box today! **Item # B-SAMPLE-S**  
The BATRIK BITESTRONG BiteBlock Single Use Sample Box includes:



Order your reusable Bite Block sample box today! **Item # B-SAMPLE-R**  
The BATRIK BITESTRONG BiteBlock Reusable Sample Box includes:



**BATRIK Safe-Soak Endo Sponge™** are specifically fabricated for use with all endoscopes. Unlike other colored sponges, our sponges are near white as they will better evidence the effectiveness of the cleaning process. Dark colors easily hide the contaminants, where with a white sponge you can visually evaluate the pre-clean efficacy.

**Available in two formats, Flat or Tubular. Both provide the identical advantages:**

- 100% Medical Grade High Density Foam
- will not flake or scratch expensive equipment
- high absorbency as it retains more solution than conventional sponges or paper wipes
- compatible with Enzymatic or Disinfectant solutions
- green friendly



**Safe-Soak Endo Sponge™**

Item #	Description
17-100	Tubular Sponge / White Dispenser Box
17-101	Tubular / Bulk
17-200	Flat / White Dispenser Box
17-201	Flat / Bulk


Packaging
25 Units / Box – 4 Boxes / Case
560 Units / Case
50 Units / Box – 4 Boxes / Case
1000 Units / Case

Note: The white sponge may at times appear slightly yellowed. This is only due to UV exposure, not a reflection to its age.

**BATRIK Couvrette™ Valve Covers** eliminate the need for manual cleaning and reprocessing of reusable valves. Our Single-Use Valve Covers effectively reduce cross contamination between patients and also reduce fluid leakage from the parts during procedures.

- easy and secure fit to all valves
- Reduce the potential for errors
- single-use
- non sterile
- latex free

**Not yet available in the USA**



**Reusable valve covers are difficult to clean. Save time and reduce risk.**

**Biopsy Valve Covers**

Item #	Description	Color	Packaging
16-0002	Biopsy Valve Cover - OLYMPUS	Blue	100 / Box
16-0004	Biopsy Valve Cover - PENTAX	Red	100 / Box

## BATRIK ENDOSCOPY Reusable Channel Brushes DAYBRUSH™


*If a facility still prefers a Reusable brush, then the DayBrush™ is the perfect solution.*

BATRIK has developed a disposable channel brush that can be used for cleaning endoscopes for an entire day. You no longer need to worry if your present reusable brush is cleaning effectively. After processing ten (10) endoscopes<sup>(1)</sup>, simply discard the brush.

### The unique structure of the Reusable DayBrush™ enables the following

- soft tip ensures easy manipulation within the channel with no risk of damage
- the back-end is sealed and softened to prevent fluid from getting trapped inside the tube. The softening process also prevents damage if the cleaning brush is used in reverse
- can be steam sterilized
- the volume of bristles on the DayBrush™ is 20mm in length — almost double most other brushes. A large brush provides more surface contact during the cleaning process
- the BATRIK Reusable DayBrush™ has been tested to effectively clean up to 10 scopes<sup>(1)</sup>
- significant per use savings when compared to reusable or disposable channel cleaning brushes

<sup>(1)</sup> Discard after cleaning 10 scopes or sooner if any brush degradation is noticed. After ten applications the brush MAY have the appearance that it can still effectively clean, but we discourage prolonged use due to risk of cross contamination. Actual number of uses will vary based on handling and processing.



#### ADD THE SUFFIX...

F= receive a kit with our Flat Sponge (Packaging: Box of 40). Applicable on certain models only, call for details.

K= receive a kit with the single use Valve / Control Cleaning Brush 5/11 Double Brushes (item # 27-57787). Applicable on all models.

FK= includes both sponge and valve cleaning brush.



### 5mm Cylindrical DayBrush™

Item #	Length inches / cm	Brush Length inches / mm	Brush (Ø) Inches / mm	Sheath (Ø) Inches / mm	Brush Ends
30-18050-T	70" / 180cm	0.79" / 20mm	0.20" / 5mm	0.07" / 1.75mm	Single
30-23050-T	90" / 230cm	0.79" / 20mm	0.20" / 5mm	0.07" / 1.75mm	Single
30-23050-2T	90" / 230cm	0.79" / 20mm	0.20" / 5mm	0.07" / 1.75mm	Double

**Packaging:** 50 / Box


## BATRIK ENDOSCOPY Single Use ARC™ Channel Cleaning Brushes

BATRIK offers one of the largest selection of Endoscopy Channel Cleaning Brushes in the industry. Our brushes are available in a variety of lengths, material, and styles to accommodate the widest range of endoscopic devices used in the healthcare industry today.

No matter what your requirements, you'll be able to find the brush that suits your specific needs. If you don't see the brush you require, just call us. We have well over 50 different channel brushes available and if you can't find the exact one that you need, no problem — we can custom fabricate to your specifications.

### Why Choose a BATRIK ARC Channel Brush?


- one of the industry's longest brush heads on the market: 15mm to 24 mm in length. More bristles translate into more cleaning contact within the channels
- bristle fan tip on the head of the brush protects and cleans the channel. No sharp end or hard metal tips that can scratch or cause damage
- latex free
- available in two different materials in order to accommodate different customer preferences:
  - POLYPROPYLENE is the industry standard as it is a low cost material and offers greater flexibility.
  - DELRIN is slightly more expensive, but due to its more rigid feeling, it is easier to insert and manipulate. This in turn reduces time and labor which can offset the added cost. ALL BATRIK small diameter brushes are supplied only with DELRIN, otherwise they would be too fragile



*F= receive a kit with our Flat Sponge (Packaging: Box of 40). Applicable on certain models only; call for details.*

*K= receive a kit with the single use Valve / Control Cleaning Brush 5/11 Double Brushes (item # 27-57787). Applicable on all models.*

*FK= includes both sponge and valve cleaning brush.*



### 5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.85 MM – PUNCHY PINK POLYPROPYLENE

Item #	Length inches / cm	Brush Length inches / mm	Brush (Ø) Inches / mm	Sheath (Ø) Inches / mm	Brush Ends
311-23050	90" / 230cm	0.79" / 20mm	0.20" / 5mm	0.073" / 1.85mm	Single
311-23050-2	90" / 230cm	0.79" / 20mm	0.20" / 5mm	0.073" / 1.85mm	Double
311-18050	70" / 180cm	0.79" / 20mm	0.20" / 5mm	0.073" / 1.85mm	Single
311-10050	39" / 100cm	0.79" / 20mm	0.20" / 5mm	0.073" / 1.85mm	Single

**Packaging:** 50 / Box

*Please note that you may receive the brush in a RED color instead of PINK. It is perfectly normal as we are slowly phasing out the RED color. The quality of the product remains identical (31 Series).*

### 5/6 MM TAPERED BRUSHES – SHEATH Ø 1.85 MM – PUNCHY PINK POLYPROPYLENE

Item #	Length inches / cm	Brush Length inches / mm	Brush (Ø) Inches / mm	Sheath (Ø) Inches / mm	Brush Ends
333-23056	90" / 230cm	0.79" / 20mm	0.20, 0.24" / 5, 6mm	0.073" / 1.85mm	Single
333-23056-2	90" / 230cm	0.79" / 20mm	0.20, 0.24" / 5, 6mm	0.073" / 1.85mm	Double

**Packaging:** 50 / Box

*Please note that you may receive the brush in a RED color instead of PINK. It is perfectly normal as we are slowly phasing out the RED color. The quality of the product remains identical (33 Series).*

**BATRIK ENDOSCOPY Single Use ARC™ Channel Cleaning Brushes****5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.75 MM – METALLIC GREY DELRIN**

Item #	Length inches / cm	Brush Length inches / mm	Brush (Ø) Inches / mm	Sheath (Ø) Inches / mm	Brush Ends
300-18050	70" / 180cm	0.79" / 20mm	0.20" / 5mm	0.070" / 1.75mm	Single
300-23050	90" / 230cm	0.79" / 20mm	0.20" / 5mm	0.070" / 1.75mm	Single
300-23050-2	90" / 230cm	0.79" / 20mm	0.20" / 5mm	0.070" / 1.75mm	Double

**Packaging:** 50 / Box

Please note that you may receive the brush in a **WHITE** color instead of **GREY**. It is perfectly normal as we are slowly phasing out the **WHITE** color. The quality of the product remains identical (30 Series)

**5/6 MM TAPERED BRUSHES – SHEATH Ø 1.75 MM – METALLIC GREY DELRIN**

Item #	Length inches / cm	Brush Length inches / mm	Brush (Ø) Inches / mm	Sheath (Ø) Inches / mm	Brush Ends
322-23056	90" / 230cm	0.79" / 20mm	0.20, 0.24" / 5, 6mm	0.070" / 1.75mm	Single
322-23056-2	90" / 230cm	0.79" / 20mm	0.20, 0.24" / 5, 6mm	0.070" / 1.75mm	Double

**Packaging:** 50 / Box

Please note that you may receive the brush in a **WHITE** color instead of **GREY**. It is perfectly normal as we are slowly phasing out the **WHITE** color. The quality of the product remains identical (32 Series)

**3 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.75 MM – METALLIC PURPLE DELRIN**

Item #	Length inches / cm	Brush Length inches / mm	Brush (Ø) Inches / mm	Sheath (Ø) Inches / mm	Brush Ends
344-23030	90" / 230cm	0.79" / 20mm	0.12" / 3mm	0.060" / 1.55mm	Single
344-18030	70" / 180cm	0.79" / 20mm	0.12" / 3mm	0.060" / 1.55mm	Single
344-10030	39" / 100cm	0.79" / 20mm	0.12" / 3mm	0.060" / 1.55mm	Single

**Packaging:** 50 / Box

Please note that you may receive the brush in a **GREEN** color instead of **PURPLE**. It is perfectly normal as we are slowly phasing out the **GREEN** color. The quality of the product remains identical (34 Series)

**2.5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.00 MM – YELLOW DELRIN**

Item #	Length inches / cm	Brush Length inches / mm	Brush (Ø) Inches / mm	Sheath (Ø) Inches / mm	Brush Ends
35-23025	90" / 230cm	0.59" / 15mm	0.098" / 2.5mm	0.039" / 1.00mm	Single

**Packaging:** 50 / Box**2 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.00 MM – YELLOW DELRIN**

Item #	Length inches / cm	Brush Length inches / mm	Brush (Ø) Inches / mm	Sheath (Ø) Inches / mm	Brush Ends
35-23020	90" / 230cm	0.59" / 15mm	0.079" / 2.0mm	0.039" / 1.00mm	Single

**Packaging:** 50 / Box

\* We are pleased to announce that as we cycle through the existing inventory our brush length will increase to 22 mm +/- 2mm



## BATRIK Valve / Control Cleaning Brushes



### VALVE / CONTROL CLEANING BRUSH – 5/7

Item #	Length inches / cm	Brush Length inches / mm	Brush 1 (Ø) inches / mm	Brush 2 (Ø) inches / mm	Description
27-05197	5.00" / 13 cm	1.00" / 25 mm	0.20" / 5 mm	0.28" / 7 mm	Dble. Trumpet

**Packaging:** (Single-use recommended) Bag of 3 - **ITEM #** 27-05197-50 / Bag of 50



### VALVE / CONTROL CLEANING BRUSH – 5/11

Item #	Length inches / cm	Brush 1 Length inches / mm	Brush 2 Length inches / mm	Brush 1 (Ø) inches / mm	Brush 2 (Ø) inches / mm	Description
27-57787	6.00" / 15cm	1.18" / 30 mm	0.79" / 20 mm	0.43" / 11 mm	0.20" / 5 mm	Dble. Universal

**Packaging:** (Single-use recommended) Bag of 3 - **ITEM #** 27-57787-50 / Bag of 50



### VALVE / CONTROL CLEANING REUSABLE BRUSH – 5/11 – EXTRA LONG

Item #	Length inches / cm	Brush 1 Length inches / mm	Brush 2 Length inches / mm	Brush 1 (Ø) inches / mm	Brush 2 (Ø) inches / mm	Description
40-18511	8.07" / 20.5cm	1.57" / 40 mm	0.79" / 20 mm	0.43" / 11 mm	0.20" / 5 mm	Dble. Universal

**Packaging:** (Re-usable and can be steam sterilized) Bag of 3



or



### VALVE / CONTROL CLEANING BRUSH – SINGLE USE BRUSH

Item #	Length inches / cm	Brush Length inches / mm	Brush (Ø) inches / mm	Description
36-10010*	4.3" / 11cm	1.57" / 40mm	0.39" / 10mm	Single end large

**Packaging:** (Single-use recommended) Bag of 3 - **ITEM #** 36-10010-50 / Bag of 50

\*Actual brush handle may vary, but brush head is unchanged.

**BATRIK ENDOSCOPY V-Brush - Reusable External Cleaning Brush** External cleaning brush for endoscopes / optical instruments. With its unique shape and innovative design, the BATRIK Reusable External Cleaning Brush is the ultimate tool for cleaning the exterior of your endoscopic instruments without the risk of damaging expensive equipment.

- easy-grip handle ensures comfort and control
- numerous dense bristles provide additional surface contact
- superior endoscopic scope cleaning with a brush specifically designed for this purpose

**Handle Length / Diameter**  
120mm / 12mm



**Brush Length / Diameter**  
90mm / 25mm

**Total Length**  
210mm

- polyamide (TYNEX®) bristles
- superior bristle durability
- steamable thermoplastic handle
- ergonomic handle
- steam / chemical sterilizable
- latex free

**Reusable External Cleaning Brush**

Item #	Description	Packaging
36-20010	Reusable External Cleaning Brush	3 units / Bag

**BATRIK ENDOSCOPY V-Brush Single Use External Cleaning Brush** External cleaning brush for endoscopes / optical instruments. An ideal fit for most endoscopic instruments, the BATRIK Single Use External Cleaning Brush will provide excellent cleaning capabilities. It's U-shape bristle design provides continuous and effective cleaning. There are no sharp or metallic ends that can scratch or cause damage to valuable equipment.

**Many of the same features as the reusable brush**

- polypropylene bristles
- vinyl caps
- ergonomic handle
- latex free

**Handle Length / Diameter**  
80mm / 3mm




**Brush Length / Diameter**  
70mm / 25mm

**Total Length**  
150mm

**Single Use External Cleaning Brush**

Item #	Description	Packaging
36-20005	Single Use External Cleaning Brush	50 units / Box





**BATRIK Silicone Pin Mat** cushions, protects and organizes delicate surgical instruments for sterilization. Residual moisture collection below the silicone mat is significantly reduced due to our unique **elevated "Nub"** design. These "Nubs" allow for greater air-flow ventilation and drying.

**New Features**

- now with **more holes** to facilitate better venting and reduce water collection
- unique **elevated nubs** also provide superior venting and a reduction in water collection
- **Custom sizes** available on request (call for details)

**Reusable and Autoclavable**  
**Latex free**



**Silicone Pin Mat**

Item #	Dimensions – inches / cm	Packaging
13-2001	13" x 18" / 33cm x 45cm	1 unit
13-2005	13" x 18" / 33cm x 45cm	5 units


www.surgmed.com • info@surgmed.com • +1 514-631-7988 • Fax +1 514-631-9083



**BATRIK Corner Pockets™** Avoid frustration and costs from having to reprocess your surgical trays due to ripped sterile wrapping. Corner Pockets™ will prevent tearing associated with sharp corners or feet on surgical trays. They will fit all types of corners (even the round ones), and the bright orange-red color makes them easily identifiable and will reduce the risk of being discarded.

We are also pleased to introduce one of the only non-foam single-use corner protectors (see below).

- no moisture retention
- reusable or single-use options
- non sterile
- latex free



**The Reusable corner pockets can withstand repeated sterilization up to 134°C; after 10 cycles, no product degradation was observed.**

**Corner Pockets™**

Item #	Description	Packaging
18-CP07	Corner Pockets™ — Reusable	80 pcs / Bag
18-CORNER	Corner Pockets™ — Reusable	4 Bags / Box (320 pcs)
18-EVAL7	Corner Pockets™ — Reusable	24 pcs / Bag
18-CP02*	Corner Pockets™ — Single-Use	240 pcs / Bag
18-EVAL2*	Corner Pockets™ — Single-Use	24 pcs / Bag

*\* Although the product appears similar to the reusable device, please be advised that it is blended differently to reduce costs and has only been validated as single use. Re-using the single use Corner Pocket may compromise sterility.*




**\*NEW\* AN INDUSTRY FIRST – SINGLE USE CORNER POCKETS *NOT* MADE OF FOAM!**

Our Single-Use corner pockets are competitive in price to foam protectors, but have the following advantages:

- lower transport costs
- requires less storage space
- save time; easier to use and position
- no moisture retention
- recyclable
- better value

**COMING SOON**  
**\*NEW\* FOAM CORNER PROTECTORS**



**BATRIK Instrument Protectors** will protect instrument tips and edges from damage during handling and sterilization.


- available in a wide array of sizes for almost every application
- sizes are color coded
- opaque and tinted protectors are radiopaque
- steam and gas permeable
- non sterile
- latex free



**Tinted**

**Variations Include**

- perforated protectors (holes) which reduce vacuum seal and moisture build-up
- osteotome protectors
- double protectors
- retractor protectors
  - "Weitlaner" Perforated "Square" Style
  - Single-Tip Perforated Clear "Tube" Style
- validated for pre-vacuum sterilization



**Opaque – Solid Color**

Item # (regular)	Item # (perforated)
<b>1</b> RS-01	PS-01
<b>2</b> RS-02	PS-02
<b>3</b> RS-03	PS-03
<b>4</b> RS-04	PS-04
<b>5</b> RS-05	PS-05
<b>6</b> RS-06	PS-06
<b>7</b> RS-07	PS-07
<b>8</b> RS-08	PS-08
<b>9</b> RS-09	PS-09
RS-10	PS-10

**Packaging: 100 units/Bag**

**Tinted Color**

Item # (regular)	Item # (perforated)	Color
RT-01	PT-01	White
RT-02	PT-02	Blue
RT-03	PT-03	Green
RT-04	PT-04	Red
RT-05	PT-05	Yellow
RT-06	PT-06	Maroon
RT-07	PT-07	Orange
RT-08	PT-08	Brown
RT-09	PT-09	Black
RT-10	PT-10	Assorted

**Packaging: 100 units/Bag**

**Clear - No Color**

Item # (regular)	Item # (perforated)
<b>1</b> RC-01	PC-01
<b>2</b> RC-02	PC-02
<b>3</b> RC-03	PC-03
<b>4</b> RC-04	PC-04
<b>5</b> RC-05	PC-05
<b>6</b> RC-06	PC-06
<b>7</b> RC-07	PC-07
<b>8</b> RC-08	PC-08
<b>9</b> RC-09	PC-09
RC-10	PC-10

**Packaging: 100 units/Bag**



**Clear**


**Solid**

### Double Protectors

Item #	Dimensions – inches / mm	Color
DG-01	0.0625 x 0.75" / 1.59 x 19.05mm	White
DG-02	0.078 x 0.75" / 1.98 x 19.05mm	Blue
DG-03	0.109 x 0.75" / 2.78 x 19.05mm	Green
DG-04	0.125 x 1" / 3.18 x 25.4mm	Red
DG-05	0.1875 x 1" / 4.76 x 25.4mm	Yellow
DG-07	0.0625 x 0.375 x 1" / 1.59 x 9.53 x 25.4mm	Orange
DG-10	Assorted	Assorted



**Packaging: 50 units/Bag - Perforated only**



### Osteotome Protectors

Item #	Dimensions – inches / mm
OP-11	0.25 x 0.75" / 6.35 x 19.05mm
OP-12	0.375 x 1" / 9.53 x 25.4mm
OP-13	0.5 x 1" / 12.7 x 25.4mm
OP-14	0.75 x 1" / 19.05 x 25.4mm
OP-15	1 x 1" / 25.4 x 25.4mm
OP-16	1.25 x 1.25" / 31.75 x 31.75mm
OP-17	1.5 x 1.5" / 38.1 x 38.1mm
OP-18	1.75 x 1.25" / 44.45 x 31.75mm
OP-20	Assorted

**Packaging: 50 units/Bag - Perforated only**



### Retractor Protectors "Weitlaner"

Item #	Size	Color
RGM-01	Small	Pink
RGM-02	Small	Clear
RGM-03	Medium	Green
RGM-04	Medium	Clear

**Packaging: 50 units/Bag**

### Retractor Protectors Single-Tip

Item #	Dimensions – inches / cm	Packaging
RGS-01	1.5" / 3.8cm	55 units / Bag
RGS-02	2" / 5.1cm	45 units / Bag
RGS-03	3" / 7.62cm	40 units / Bag



## BATRIK S.T.I.C.S.™ a Simple Tracking Identification Color System

- quickly identifies instruments by department or procedure
- easy to apply
- long lasting
- tape adheres securely
- a free protective cover is provide with every sheet
- latex free
- non sterile

**Each 8.5 x 11" (216 x 279mm) sheet contains the following adhesive strips:**

- 1" (25.4mm) x 220
- 1.5" (38.1mm) x 44
- 2" (50.8mm) x 44
- total tape per sheet 374" (9500mm)

### Solid Colors

Item #	Color	
T-850-2-01	White	
T-850-2-02	Royal Blue	
T-850-2-03	Kelly Green	
T-850-2-04	Red	
T-850-2-05	Yellow	
T-850-2-06	Purple	
T-850-2-07	Orange	
T-850-2-08	Brown	
T-850-2-09	Black	
T-850-2-10	Gray	
T-850-2-11	Lavender	
T-850-2-12	Beige	
T-850-2-13	Light Blue	
T-850-2-14	Lime Green	
T-850-2-15	Pink	
T-850-2-17	Plum	
T-850-2-18	Mint Green	

### Black Stripes

Item #	Color	
T-850-3-01	Black / White	
T-850-3-03	Black / Kelly Green	
T-850-3-04	Black / Red	
T-850-3-05	Black / Yellow	
T-850-3-07	Black / Orange	
T-850-3-10	Black / Gray	
T-850-3-11	Black / Lavender	
T-850-3-12	Black / Beige	
T-850-3-13	Black / Light Blue	
T-850-3-14	Black / Lime Green	
T-850-3-15	Black / Pink	
T-850-3-17	Black / Plum	

### White Stripes

Item #	Color	
T-850-4-02	White / Royal Blue	
T-850-4-03	White / Kelly Green	
T-850-4-04	White / Red	
T-850-4-05	White / Yellow	
T-850-4-07	White / Orange	
T-850-4-08	White / Brown	
T-850-4-11	White / Lavender	
T-850-4-13	White / Light Blue	
T-850-4-14	White / Lime Green	
T-850-4-15	White / Pink	
T-850-4-17	White / Plum	

### Tri-Color Stripe

Item #	Color	
T-850-5-51	Red / White / Blue	
T-850-5-54	Red / Black / Green	

### Specialty Stripes

Item #	Color	
T-850-7-04	Yellow / Red	

### Four Color

Item #	Color
T-850-2-31	Red / Yellow / Purple / Orange
T-850-2-32	Blue / White / Black / Green

**BATRIK Identification Tapes** are a simple and inexpensive identification solution for marking surgical instruments.

- quickly identifies instruments by department or procedure
- easy to apply
- withstands repeated exposure to sterilization
- has more tape per roll (300" / 762cm)
- are available in 9 colors and 2 widths (0.25" / 0.635cm and 0.125" / 0.3175cm)
- packaged in a convenient clear box
- non sterile

\*sterilization reports are available. All samples tested validated the efficiency of sterilization, even underneath the tape

### Identification Tapes

Item #	Item #	Color
0.125" / 0.3175 cm	0.25" / 0.635 cm	
T-125-01	T-250-01	White
T-125-02	T-250-02	Blue
T-125-03	T-250-03	Green
T-125-04	T-250-04	Red
T-125-05	T-250-05	Yellow
T-125-06	T-250-06	Purple
T-125-07	T-250-07	Orange
T-125-08	T-250-08	Brown
T-125-09	T-250-09	Black





**BATRIK Color Coding System Padlock Security Tags** plastic seals for instrument containers and general hospital applications are specifically designed for the healthcare market. Our material selection, fabrication process and engineering ensures consistent quality and high temperature compliance with the various steam sterilization temperatures in world wide hospital applications.

- larger base permits easier handling and application
- EasyTwist breaking design
- fully open design (45°) Facilitates insertion into instrument containers
- will not melt or soften with any established steam sterilization protocol
- specifically manufactured for the healthcare market
- made of polypropylene
- latex and DEHP free
- custom colors and OEM available on request



**EasyTwist Padlocks**

**Regular Padlocks**

**With Indicator Dot**

Color	100 Units	200 Units	500 Units	100 Units	200 Units	500 Units
White	15-701	15-711	15-710	15-701-IND	15-711-IND	15-710-IND
Blue	15-702	15-712	15-720	15-702-IND	15-712-IND	15-720-IND
Green	15-703	15-713	15-730	15-703-IND	15-713-IND	15-730-IND
Red	15-704	15-714	15-740	15-704-IND	15-714-IND	15-740-IND
Yellow	15-705	15-715	15-750	15-705-IND	15-715-IND	15-750-IND
Pink	15-706	15-716	15-760	15-706-IND	15-716-IND	15-760-IND
Orange	15-707	15-717	15-770	15-707-IND	15-717-IND	15-770-IND
Black	15-709	15-719	15-790	15-709-IND	15-719-IND	15-790-IND

**BATRIK Printed Padlock Security Tags** allow you to quickly identify instruments or containers that require specific attention



**Printed Padlock Security Tags**

Item #	Color and Description	Packaging
15-701C	White "CLEAN"	100 / Package
15-702ST	Blue "STERILE"	100 / Package
15-704R	Red "REPAIR"	100 / Package
15-705S	Yellow "SHARPEN"	100 / Package
15-705L	Yellow "LOANER"	100 / Package
15-709D	Black "DISCARD"	100 / Package



**BATRIK Superior Face Protection SAVEFACE™ Shields** use superior material (7.5 mil thickness PET) to ensure excellent visibility and extended usage due to the design and material quality.

**Features**

- flat viewing area that reduces distortion and glare
- self-shaping foam head piece molds comfortably to forehead
- high quality 7.5 mil thickness PET
- shields can be used for longer periods due to comfort and improved materials
- superior anti-fog properties on both sides of the shield
- wide elastic band 1" (25mm) provides additional comfort
- meets OSHA Blood Borne Pathogen Regulation
- ideal safety addition to all pandemic kits
- FDA 510K registered
- latex free

**REDUCE YOUR SHIELD USAGE DUE TO OUR IMPROVED DESIGN, QUALITY AND COMFORT!**



**SaveFace™ Shield**

Item #	Description	Packaging**
SF-100	SAVEFACE™ Shields – 9" (23cm)	100 / Box



**SafeFace™ Shield with Drape**


Item #	Description	Packaging
SF-180	SAVEFACE™ 15" (38cm) Shield with Drape – 6" L (15cm)	100 / Box
SF-180L	SAVEFACE™ 21" (53cm) Shield with Extra Long Drape – 12" L x 8" W (30cm x 20cm)	100 / Box
SF-180W	SAVEFACE™ 20" (51cm) Shield with Extra Wide Drape – 11" L x 13" W (28cm x 33cm)	100 / Box

*Fluid drape barrier includes closing system below the chin for even greater protection*

*\*\* More units per box reduces transport costs and waste.*

**BATRIK Scrub Pads** are ideal for pre-cleaning your surgical instruments without leaving any metal residue that often occurs when using steel wool or metal sponges. The abrasiveness is perfectly balanced to provide you with efficient cleaning without damaging your instruments. Its fine fibers allow the users to **polish**, **clean** and **debur** all type of material surfaces.

- Abrasive Material: Aluminum Oxide
- Grit: 400
- Style: Flexible
- Pad material: Non-Woven
- Color: Maroon
- Grade: Very Fine
- Type: Hand Pad



### Scrub Pads

Item #	Description	Size	Packaging
60-46050	Scrub Pads for Instruments	4" x 6" (10 cm x 15 cm)	3 u. / Bag - 150 units / box
60-46003	Scrub Pads for Instruments	4" x 6" (10 cm x 15 cm)	3 u. / bag





## BATRIK STERILIZATION Single Use and Reusable Wire Style Brushes - ARC™

### Brush Selection Tips:

#### PART 1: Instrument Specifications

##### A. What is the maximum length of the instrument(s) to be cleaned?

Add at least 2" (5cm) to that length, so you have an adequate length to work with.

##### B. What is the inner diameter of the instrument?

Add at most 0.125" (3mm) to determine the diameter of the bristle to ensure effective cleaning. The bristle diameter has to slightly exceed the inner tube diameter. Bristles that are too large will bend too much and not clean effectively.

#### PART 2: Brush Selection

Now that you have the basic information for your required brush, simply find the section in our catalogue that best corresponds to the overall length required as determined in PART 1. Now select the brush that best corresponds to the diameter determined in PART 2.

As a rule, if your present brush is 8" (20.32cm) in length and we only offer 10" (25.4cm) length with your required diameter, that is entirely acceptable.

In many cases, the BATRIK Brush will have more bristles on the brush. As an example, many 5mm (0.19") diameter brushes are only 1" (2.5cm) in length. At BATRIK, we have strived to provide the greatest amount of bristles to ensure the most effective cleaning.


#### At Last... Customizable Wire Brushes!

A unique aspect of single-use ARC™ Wire Brushes is that they can be customized to any length:

- up to 2000mm for Wire Brushes on plastic semi-rigid sheaths, rigid brushes from 5 to 20mm in diameter
- no limit in length for Wire Brushes on flexible plastic sheaths, flexible brushes from 2-5mm in diameter

#### BATRIK Arc Brushes

ARC™ Wire Brushes are guaranteed to last longer than twisted wire and preserve the integrity of your valuable instruments.



**RIGID**  
BRUSH & SHEATH

### 12" / 30cm Single-End Wire Brush

	Single Use Item #	Steamable Item #	Color	Total Length Inches / cm	Brush (Ø) Inches / mm	Brush Length Inches / mm	Sheath (Ø) Inches / mm
<b>FLEXIBLE</b> BRUSH & SHEATH	41-31502-S	41-31502	White	12.401" / 31.5cm	0.078" / 2mm	0.59" / 15mm	0.39" / 1mm
	41-32503-S	41-32503	Yellow	12.795" / 32.5cm	0.118" / 3mm	0.98" / 25mm	0.39" / 1mm
	41-32504-S	41-32504	Green	12.795" / 32.5cm	0.157" / 4mm	0.98" / 25mm	0.61" / 1.55mm
	41-32505-S	41-32505	White	12.795" / 32.5cm	0.19" / 5mm	0.98" / 25mm	0.68" / 1.75mm
	41-35005-S	41-35005	White	13.779" / 35cm	0.19" / 5mm	1.968" / 50mm	0.118" / 3mm
<b>RIGID</b> BRUSH & SHEATH	41-35006-S	41-35006	Green	13.779" / 35cm	0.23" / 6mm	1.968" / 50mm	0.118" / 3mm
	41-35008-S	41-35008	Blue	13.779" / 35cm	0.31" / 8mm	1.968" / 50mm	0.157" / 4mm
	41-35010-S	41-35010	Red	13.779" / 35cm	0.39" / 10mm	1.968" / 50mm	0.157" / 4mm
	41-35012-S	41-35012	Yellow	13.779" / 35cm	0.49" / 12.5mm	1.968" / 50mm	0.157" / 4mm
	41-35015-S	41-35015	Green	13.779" / 35cm	0.59" / 15mm	1.968" / 50mm	0.228" / 5.8mm
	41-35020-S	41-35020	Blue	13.779" / 35cm	0.78" / 20mm	1.968" / 50mm	0.228" / 5.8mm

**Packaging:** Single Use: Bag of 20 / Steamable: Bag of 3

## 18" / 45cm Single-End Wire Brush

FLEXIBLE

BRUSH &amp; SHEATH

Single Use Item #	Steamable Item #	Color	Total Length Inches / cm	Brush (Ø) Inches / mm	Brush Length Inches / mm	Sheath (Ø) Inches / mm
41-46502-S	41-46502	White	18.307" / 46.5cm	0.078" / 2mm	0.59" / 15mm	0.39" / 1mm
41-47503-S	41-47503	Yellow	18.7" / 47.5cm	0.118" / 3mm	0.98" / 25mm	0.39" / 1mm
41-47504-S	41-47504	Green	18.7" / 47.5cm	0.157" / 4mm	0.98" / 25mm	0.61" / 1.55mm
41-47505-S	41-47505	White	18.7" / 47.5cm	0.19" / 5mm	0.98" / 25mm	0.68" / 1.75mm
41-50005-S	41-50005	White	19.685" / 50cm	0.19" / 5mm	1.968" / 50mm	0.118" / 3mm
41-50006-S	41-50006	Green	19.685" / 50cm	0.23" / 6mm	1.968" / 50mm	0.118" / 3mm
41-50008-S	41-50008	Blue	19.685" / 50cm	0.31" / 8mm	1.968" / 50mm	0.157" / 4mm
41-50010-S	41-50010	Red	19.685" / 50cm	0.39" / 10mm	1.968" / 50mm	0.157" / 4mm
41-50012-S	41-50012	Yellow	19.685" / 50cm	0.49" / 12.5mm	1.968" / 50mm	0.157" / 4mm
41-50015-S	41-50015	Green	19.685" / 50cm	0.59" / 15mm	1.968" / 50mm	0.228" / 5.8mm
41-50020-S	41-50020	Blue	19.685" / 50cm	0.78" / 20mm	1.968" / 50mm	0.228" / 5.8mm

Packaging: Single Use: Bag of 20 / Steamable: Bag of 3

## 24" / 60cm Single-End Wire Brush

FLEXIBLE

BRUSH &amp; SHEATH

Single Use Item #	Steamable Item #	Color	Total Length Inches / cm	Brush (Ø) Inches / mm	Brush Length Inches / mm	Sheath (Ø) Inches / mm
41-61502-S	41-61502	White	24.212" / 61.5cm	0.078" / 2mm	0.59" / 15mm	0.39" / 1mm
41-62503-S	41-62503	Yellow	24.6" / 62.5cm	0.118" / 3mm	0.98" / 25mm	0.39" / 1mm
41-62504-S	41-62504	Green	24.6" / 62.5cm	0.157" / 4mm	0.98" / 25mm	0.61" / 1.55mm
41-62505-S	41-62505	White	24.6" / 62.5cm	0.19" / 5mm	0.98" / 25mm	0.68" / 1.75mm
41-65005-S	41-65005	White	25.59" / 65cm	0.19" / 5mm	1.968" / 50mm	0.118" / 3mm
41-65006-S	41-65006	Green	25.59" / 65cm	0.23" / 6mm	1.968" / 50mm	0.118" / 3mm
41-65008-S	41-65008	Blue	25.59" / 65cm	0.31" / 8mm	1.968" / 50mm	0.157" / 4mm
41-65010-S	41-65010	Red	25.59" / 65cm	0.39" / 10mm	1.968" / 50mm	0.157" / 4mm
41-65012-S	41-65012	Yellow	25.59" / 65cm	0.49" / 12.5mm	1.968" / 50mm	0.157" / 4mm
41-65015-S	41-65015	Green	25.59" / 65cm	0.59" / 15mm	1.968" / 50mm	0.228" / 5.8mm
41-65020-S	41-65020	Blue	25.59" / 65cm	0.78" / 20mm	1.968" / 50mm	0.228" / 5.8mm

Packaging: Single Use: Bag of 20 / Steamable: Bag of 3

## 30" / 75cm Single-End Wire Brush


FLEXIBLE

BRUSH &amp; SHEATH


Single Use Item #	Steamable Item #	Color	Total Length Inches / cm	Brush (Ø) Inches / mm	Brush Length Inches / mm	Sheath (Ø) Inches / mm
41-76502-S	41-76502	White	30.118" / 76.5cm	0.078" / 2mm	0.59" / 15mm	0.39" / 1mm
41-77503-S	41-77503	Yellow	30.511" / 77.5cm	0.118" / 3mm	0.98" / 25mm	0.39" / 1mm
41-77504-S	41-77504	Green	30.511" / 77.5cm	0.157" / 4mm	0.98" / 25mm	0.61" / 1.55mm
41-77505-S	41-77505	White	30.511" / 77.5cm	0.19" / 5mm	0.98" / 25mm	0.68" / 1.75mm
41-80005-S	41-80005	White	31.496" / 80cm	0.19" / 5mm	1.968" / 50mm	0.118" / 3mm
41-80006-S	41-80006	Green	31.496" / 80cm	0.23" / 6mm	1.968" / 50mm	0.118" / 3mm
41-80008-S	41-80008	Blue	31.496" / 80cm	0.31" / 8mm	1.968" / 50mm	0.157" / 4mm
41-80010-S	41-80010	Red	31.496" / 80cm	0.39" / 10mm	1.968" / 50mm	0.157" / 4mm
41-80012-S	41-80012	Yellow	31.496" / 80cm	0.49" / 12.5mm	1.968" / 50mm	0.157" / 4mm
41-80015-S	41-80015	Green	31.496" / 80cm	0.59" / 15mm	1.968" / 50mm	0.228" / 5.8mm
41-80020-S	41-80020	Blue	31.496" / 80cm	0.78" / 20mm	1.968" / 50mm	0.228" / 5.8mm

Packaging: Single Use: Bag of 20 / Steamable: Bag of 3

**FLEXIBLE**



**RIGID**





## BATRIK STERILIZATION Reusable Stainless Steel ARC™ Micro Brushes

### Our Unique 1.5 mm Stainless Steel Brush is Now Available!


**NEW!**  
1.5 mm Brush

One of the smallest brushes in the industry: They are designed to clean your channels with an inner diameter ranging from 0.5 mm to 1.40 mm.

It is highly recommended using the **pull-through technique** (insert the blunt, non-brush end first).

#### Features

- Flexible stainless steel sheath - will not bend like the regular metal twisted wire brushes
- Stainless steel sheath
- Reusable\*
- Packaging: 6 units / Tube



#### 1.5 mm Brush

Item #	Total Length Inches / cm	Brush Size (Ø) Inches / mm	Brush Length Inches / mm	Sheath (Ø) Inches / mm
35-03015	12" / 30 cm	0.06" / 1.5 mm	0.6" / 15 mm	0.02" / 0.40 mm
35-04515	18" / 45 cm	0.06" / 1.5 mm	0.6" / 15 mm	0.02" / 0.40 mm

**\* Due to the delicate nature of these brushes, although all the material used in the fabrication can be re-processed repeatedly (as with any other wire brush), any re-use is conditional as to the handling and application of this brush.**

**BATRIK STERILIZATION Reusable Asymmetric Double ARC™ Wire Brushes****14" / 35cm Asymmetric Double-End Wire Brush**

Item #	Color	Length Inches / cm	Brush 1 (Ø) Inches / mm	Brush 2 (Ø) Inches / mm	Brush Length Inches / mm	Sheath (Ø) Inches / mm
40-35056	White	13.779" / 35cm	0.19" / 5 mm	0.23" / 6mm	1.968" / 50mm	0.118" / 3mm
40-35810	Blue	13.779" / 35cm	0.31" / 8mm	0.39" / 10mm	1.968" / 50mm	0.157" / 4mm
40-35112	Yellow	13.779" / 35cm	0.39" / 10mm	0.49" / 12.5mm	1.968" / 50mm	0.157" / 4mm
40-35152	Green	13.779" / 35cm	0.59" / 15mm	0.78" / 20mm	1.968" / 50mm	0.228" / 5.8mm

**Packaging:** Bag of 3**20" / 50cm Asymmetric Double-End Wire Brush**


Item #	Color	Length Inches / cm	Brush 1 (Ø) Inches / mm	Brush 2 (Ø) Inches / mm	Brush Length Inches / mm	Sheath (Ø) Inches / mm
40-50056	White	19.685" / 50cm	0.19" / 5 mm	0.23" / 6mm	1.968" / 50mm	0.118" / 3mm
40-50810	Blue	19.685" / 50cm	0.31" / 8mm	0.39" / 10mm	1.968" / 50mm	0.157" / 4mm
40-50112	Yellow	19.685" / 50cm	0.39" / 10mm	0.49" / 12.5mm	1.968" / 50mm	0.157" / 4mm
40-50152	Green	19.685" / 50cm	0.59" / 15mm	0.78" / 20mm	1.968" / 50mm	0.228" / 5.8mm

**Packaging:** Bag of 3**26" / 65cm Asymmetric Double-End Wire Brush**

Item #	Color	Length Inches / cm	Brush 1 (Ø) Inches / mm	Brush 2 (Ø) Inches / mm	Brush Length Inches / mm	Sheath (Ø) Inches / mm
40-65056	White	25.59" / 65cm	0.19" / 5 mm	0.23" / 6mm	1.968" / 50mm	0.118" / 3mm
40-65810	Blue	25.59" / 65cm	0.31" / 8mm	0.39" / 10mm	1.968" / 50mm	0.157" / 4mm
40-65112	Yellow	25.59" / 65cm	0.39" / 10mm	0.49" / 12.5mm	1.968" / 50mm	0.157" / 4mm
40-65152	Green	25.59" / 65cm	0.59" / 15mm	0.78" / 20mm	1.968" / 50mm	0.228" / 5.8mm

**Packaging:** Bag of 3**31" / 80cm Asymmetric Double-End Wire Brush**

Item #	Color	Length Inches / cm	Brush 1 (Ø) Inches / mm	Brush 2 (Ø) Inches / mm	Brush Length Inches / mm	Sheath (Ø) Inches / mm
40-80056	White	31.496" / 80cm	0.19" / 5 mm	0.23" / 6mm	1.968" / 50mm	0.118" / 3mm
40-80810	Blue	31.496" / 80cm	0.31" / 8mm	0.39" / 10mm	1.968" / 50mm	0.157" / 4mm
40-80112	Yellow	31.496" / 80cm	0.39" / 10mm	0.49" / 12.5mm	1.968" / 50mm	0.157" / 4mm
40-80152	Green	31.496" / 80cm	0.59" / 15mm	0.78" / 20mm	1.968" / 50mm	0.228" / 5.8mm

**Packaging:** Bag of 3

## BATRIK STERILIZATION Reusable Symmetric Double ARC™ Wire Brushes

### 14" / 35cm Symmetric Double-End Wire Brushes

Item #	Color	Length Inches / cm	Brush 1 (Ø) Inches / mm	Brush Length Inches / mm	Sheath (Ø) Inches / mm
40-35005	White	13.779" / 35cm	0.19" / 5 mm	1.968" / 50mm	0.118" / 3mm
40-35006	Green	13.779" / 35cm	0.23" / 6mm	1.968" / 50mm	0.118" / 3mm
40-35008	Blue	13.779" / 35cm	0.31" / 8mm	1.968" / 50mm	0.157" / 4mm
40-35010	Red	13.779" / 35cm	0.39" / 10mm	1.968" / 50mm	0.157" / 4mm
40-35012	Yellow	13.779" / 35cm	0.49" / 12.5mm	1.968" / 50mm	0.157" / 4mm
40-35015	Green	13.779" / 35cm	0.59" / 15mm	1.968" / 50mm	0.228" / 5.8mm
40-35020	Blue	13.779" / 35cm	0.78" / 20mm	1.968" / 50mm	0.228" / 5.8mm

**Packaging:** Bag of 3

### 20" / 50cm Symmetric Double-End Wire Brushes

Item #	Color	Length Inches / cm	Brush 1 (Ø) Inches / mm	Brush Length Inches / mm	Sheath (Ø) Inches / mm
40-50005	White	19.685 / 50cm	0.19" / 5 mm	1.968" / 50mm	0.118" / 3mm
40-50006	Green	19.685 / 50cm	0.23" / 6mm	1.968" / 50mm	0.118" / 3mm
40-50008	Blue	19.685 / 50cm	0.31" / 8mm	1.968" / 50mm	0.157" / 4mm
40-50010	Red	19.685 / 50cm	0.39" / 10mm	1.968" / 50mm	0.157" / 4mm
40-50012	Yellow	19.685 / 50cm	0.49" / 12.5mm	1.968" / 50mm	0.157" / 4mm
40-50015	Green	19.685 / 50cm	0.59" / 15mm	1.968" / 50mm	0.228" / 5.8mm
40-50020	Blue	19.685 / 50cm	0.78" / 20mm	1.968" / 50mm	0.228" / 5.8mm

**Packaging:** Bag of 3

### 26" / 65cm Symmetric Double-End Wire Brushes

Item #	Color	Length Inches / cm	Brush 1 (Ø) Inches / mm	Brush Length Inches / mm	Sheath (Ø) Inches / mm
40-65005	White	25.59" / 65cm	0.19" / 5 mm	1.968" / 50mm	0.118" / 3mm
40-65006	Green	25.59" / 65cm	0.23" / 6mm	1.968" / 50mm	0.118" / 3mm
40-65008	Blue	25.59" / 65cm	0.31" / 8mm	1.968" / 50mm	0.157" / 4mm
40-65010	Red	25.59" / 65cm	0.39" / 10mm	1.968" / 50mm	0.157" / 4mm
40-65012	Yellow	25.59" / 65cm	0.49" / 12.5mm	1.968" / 50mm	0.157" / 4mm
40-65015	Green	25.59" / 65cm	0.59" / 15mm	1.968" / 50mm	0.228" / 5.8mm
40-65020	Blue	25.59" / 65cm	0.78" / 20mm	1.968" / 50mm	0.228" / 5.8mm

**Packaging:** Bag of 3

### 31" / 80cm Symmetric Double-End Wire Brushes

Item #	Color	Length Inches / cm	Brush 1 (Ø) Inches / mm	Brush Length Inches / mm	Sheath (Ø) Inches / mm
40-80005	White	31.496" / 80cm	0.19" / 5 mm	1.968" / 50mm	0.118" / 3mm
40-80006	Green	31.496" / 80cm	0.23" / 6mm	1.968" / 50mm	0.118" / 3mm
40-80008	Blue	31.496" / 80cm	0.31" / 8mm	1.968" / 50mm	0.157" / 4mm
40-80010	Red	31.496" / 80cm	0.39" / 10mm	1.968" / 50mm	0.157" / 4mm
40-80012	Yellow	31.496" / 80cm	0.49" / 12.5mm	1.968" / 50mm	0.157" / 4mm
40-80015	Green	31.496" / 80cm	0.59" / 15mm	1.968" / 50mm	0.228" / 5.8mm
40-80020	Blue	31.496" / 80cm	0.78" / 20mm	1.968" / 50mm	0.228" / 5.8mm

**Packaging:** Bag of 3

## BATRIK STERILIZATION In-Brush™ **ANTIMICROBIAL** Twisted Wire Brushes

### Twisted Wire Brushes

Designed for use with suction tubes, our twisted wires have a fine / medium diameter yet are rigid enough to ensure excellent brush cleaning capabilities

### Nylon Twisted Wire Brushes - with **Antimicrobial** Bristles

We offer a wide selection of brushes for most of your instrument cleaning applications. Ideal for large bore, uniquely shaped, and endoscopic instruments.



#### Frazier / Baron Suction Tubes 10" / 250mm

Item #	Size	Length Inches / mm	Brush Length Inches / mm	Brush (Ø) Inches / mm
20-10079	6FR	10" / 250mm	1" / 25.4mm	0.079" / 2mm
20-10090	7FR	10" / 250mm	1" / 25.4mm	0.090" / 2.3mm
20-10106	8FR	10" / 250mm	1" / 25.4mm	0.106" / 2.7mm
20-10118	9FR	10" / 250mm	1" / 25.4mm	0.118" / 3mm
20-10129	10FR	10" / 250mm	1.5" / 38mm	0.129" / 3.3mm
20-10146	11FR	10" / 250mm	1.5" / 38mm	0.146" / 3.7mm
20-10157	12FR	10" / 250mm	1.5" / 38mm	0.157" / 4.0mm

**Packaging:** Bag of 3

#### Nylon Twisted Wire Brushes 12" / 305mm

Item #	Length Inches / mm	Brush Length Inches / mm	Brush (Ø) Inches / mm
22-12157	12" / 305mm	3" / 76mm	0.157" / 4mm
22-12197	12" / 305mm	3" / 76mm	0.197" / 5mm
22-12236	12" / 305mm	3" / 76mm	0.236" / 6mm
22-12393	12" / 305mm	3" / 76mm	0.393" / 10mm

**Packaging:** Bag of 3



#### Nylon Twisted Wire Brushes 16" / 406mm

Item #	Length Inches / mm	Brush Length Inches / mm	Brush (Ø) Inches / mm
23-16197	16" / 406mm	2.5" / 63.5mm	0.197" / 5mm
23-16236	16" / 406mm	2.5" / 63.5mm	0.236" / 6mm
23-16276	16" / 406mm	3" / 76mm	0.276" / 7mm
23-16394	16" / 406mm	3" / 76mm	0.394" / 10mm
23-16472	16" / 406mm	3" / 76mm	0.472" / 12mm
23-16591	16" / 406mm	3" / 76mm	0.591" / 15mm
23-16787	16" / 406mm	3" / 76mm	0.787" / 20mm

**Packaging:** Bag of 3

#### Nylon Twisted Wire Brushes 24" / 610mm

Item #	Length Inches / mm	Brush Length Inches / mm	Brush (Ø) Inches / mm
24-24098	24" / 610mm	2.5" / 63.5mm	0.118" / 3mm
24-24197	24" / 610mm	2.5" / 63.5mm	0.197" / 5mm
24-24276	24" / 610mm	3.25" / 82.6mm	0.276" / 7mm
24-24394	24" / 610mm	3.25" / 82.6mm	0.394" / 10mm
24-24472	24" / 610mm	3.25" / 82.6mm	0.472" / 12mm

**Packaging:** Bag of 3



#### Broncho Brushes

Item #	Description	Length Inches / mm	Brush Length Inches / mm
26-06175	Small	6" / 152mm	1.75" / 44.5mm
26-08325	Large	8" / 203mm	3.25" / 82.6mm

**Packaging:** Bag of 3

#### Tracheostomy Brushes

Item #	Description	Length Inches / mm	Brush Length Inches / mm	Brush (Ø) Inches / mm
25-08500	Large 6-12.3 <sup>3</sup> / <sub>4</sub> "	8" / 203mm	3.75" / 95mm	0.5" / 12.7mm
25-37550	Small 000-5.1 <sup>7</sup> / <sub>8</sub> "	4" / 93mm	1.88" / 48mm	0.5" / 12.7mm

**Packaging:** Bag of 3

**BATRIK Instrument Cleaning Flat Brushes** are specifically designed for hospital use in the pre-clean and preparation of instruments.

- long, ergonomic designed handle permits ease-of-use
- Extra long length of bristles (70 mm)
- the brushes fibers are transparent, allowing for quick verification of its cleanliness
- the plastic composition used to fabricate the handle will not release any toxic components

BATRIK manufactures four types of flat brushes:

- **EXTRA RIGID** instrument safe polyamide fibers offer steel-like rigidity and allows deep cleaning of soiled instruments, particularly in grooved or uneven surfaces.
- **RIGID** the rigidity of this brush, along with the length of the bristles and handle, improves and simplifies instrument cleaning.
- **FLEXIBLE** This brush is designed for large instruments that require a more gentle cleaning application such as optics, plastics or hard to reach spots which are easier to clean using a flexible bristle.
- **EXTRA FLEXIBLE** The very long bristles eliminate any stress during sanitation, and makes this the ideal brush for fragile instruments(optics, etc.) as well as both pliable and rigid endoscopes. The length of the bristles assist in cleaning areas that would otherwise be difficult to reach with shorter bristles, such as the accessory channel or tip deflection of pliable endoscopes



#### Extra Rigid Reusable / Autoclavable Flat Brush

Item #	Description	Length – mm	Brush Size – mm
50-40917	Reusable Flat brush (extra rigid)	220mm	77L x 20W x10H

**Packaging:** Bag of 3



#### Flexible Reusable / Autoclavable Flat Brush

Item #	Description	Length – mm	Brush Size – mm
50-40918	Reusable Flat brush (rigid)	220mm	77L x 20W x15H

**Packaging:** Bag of 3



#### Rigid Reusable / Autoclavable Flat Brush

Item #	Description	Length – mm	Brush Size – mm
50-40919	Reusable Flat brush (flexible)	220mm	77L x 20W x15H

**Packaging:** Bag of 3



#### Extra Flexible Reusable / Autoclavable Flat Brush

Item #	Description	Length – mm	Brush Size – mm
50-40920	Reusable Flat brush (extra flexible)	220mm	77L x 20W x25H

**Packaging:** Bag of 3




#### Rigid and Flexible Single Use / Autoclavable Flat Brush

Item #	Description	Length – mm	Brush Size – mm
50-40921	Single use flat brush (rigid)	220mm	60L x 20W x15H
50-40922	Single use flat brush (flexible)	220mm	60L x 20W x 15H

**Packaging:** Bag of 20

#### POLYAMIDE FIBERS



**BATRIK Instrument Cleaning Reusable Precision Brush** External cleaning brush is designed for cleaning instrument jaws. It is the ideal brush for the quick removal of difficult debris, the BATRIK Reusable Precision Brush can also easily clean deep into instrument jaws. Although the firm and durable bristles feel like stainless steel, there is absolutely no risk of instrument damage or personal injury.

- polyamide (TYNEX®) bristles
- bristles have an extra firm density and superior durability
- steamable thermoplastic handle
- ergonomic handle
- steam / chemical sterilizable
- latex free

No more  
Steel brushes



**ELIMINATES THE NEED FOR THE STAINLESS STEEL TOOTHBRUSH-STYLE BRUSH, WHICH CAN DAMAGE VALUABLE INSTRUMENTS!**




“**THROW AWAY YOUR WIRE BRUSHES, THIS IS THE BEST OPTION TO EXTEND YOUR INSTRUMENT USEFUL LIFE**”

#### Reusable Precision Brush

Item #	Description	Packaging
50-40800	Reusable Precision External Cleaning Brush Overall length is 6.29" / 160mm	3 units / Bag

**BATRIK Instrument Cleaning Reusable Jar / Bottle Brush** with its long and flexible bristles will maximize the effectiveness of large-surface cleaning. The easy-grip handle also provides excellent maneuverability. Designed specifically for hospital applications.

- high quality polyamide clear bristles
- easy to identify debris
- NEW ergonomic handle
- autoclavable
- latex free



**Reusable Jar / Bottle Brush**

**Item #**  
50-10295

**Description**  
Reusable Jar / Bottle Brush – Autoclavable  
Overall length is 13.19" / 335mm


**Packaging**  
3 units / Bag

**BATRIK Instrument Cleaning ANTIMICROBIAL Acetabular Reamer Brush - Durable and Reusable Cleaning Brush.** This double-ended cleaning brush is manufactured exclusively for cleaning acetabular bone reamers. The "Pig Tail" rounded style circular head brush is designed to effectively clean inside the reamer, while the rigid straight brush on the reverse side is used for cleaning individual cutting holes. Special fusion of the brush heads to the hand piece along with using high quality **antimicrobial** medical grade bristles ensures a durable and reusable effective cleaning tool.

#### Product Specifications\*

- overall length: 8" (200mm)
- round brush head diameter: 1.75" (45mm)
- straight brush diameter: 5/16" (8mm) Tapered
- straight brush length: 1 3/16" (30mm)

\*specifications subject to small variances



Brush ±  
1.75" / 45mm

Brush Length  
8" / 200mm

Tapered Brush ±  
0.3125" / 8mm

#### Acetabular Reamer Cleaning Brush

Item #	Description	Packaging
50-ACETAB	Acetabular Reamer Cleaning Brush	3 units / Bag



## BATRIK INSTRUMENT CLEANING All Purpose Brushes

### Plastic Handle Brushes

We offer a wide selection of brushes for most of your instrument cleaning applications. These brushes are ideal for large bore, uniquely shaped and endoscopic instruments.

### In-Brush™ Plastic Handle

Choose either the regular toothbrush style (nylon or stainless) or try our double-ended brush.

\*order in volume and save even more



### Plastic Single End Nylon – Handle Brush

Item #	Description	Length Inches / mm	Brush Length Inches / mm	Brush (Ø) Inches / mm
28-72550	Single End Nylon	7.25" / 184mm	1.5" / 38mm	0.5" / 12.7mm

**Packaging:** Bag of 3 - ITEM # 28-72550-50 / Bag of 50



### Plastic Single End Stainless Steel – Handle Brush

Item #	Description	Length Inches / mm	Brush Length Inches / mm	Brush (Ø) Inches / mm
29-72550	Single End Stainless Steel	7.25" / 184mm	1.5" / 38mm	0.5" / 12.7mm

**Packaging:** Bag of 3 - ITEM # 29-72550-50 / Bag of 50



### Plastic Double End Nylon – Handle Brush

Item #	Description	Length Inches / mm	Brush Length Inches / mm	Brush (Ø) Inches / mm
28-69394	Large Brush-end	6.9" / 175mm	1.38" / 35mm	0.39" / 10mm
	Small Brush-end	6.9" / 175mm	1" / 25mm	0.2" / 5mm

**Packaging:** Bag of 3 - ITEM # 28-69394-50 / Bag of 50



### Hand Scrub Nylon Brush – Non-Autoclavable

Item #	Length Inches / mm	Brush Length Inches / mm	Brush (Ø) Inches / mm
50-40911	4.09" / 104mm	3.74" / 95mm	1.1" / 28mm

**Packaging:** Bag of 3



## BATRIK Instrument Autoclave Cleaning Brush Kits

### Regular Autoclave Cleaning Brush Solution Kit includes:

- **1 Hand Brush**  
Length: 4.5 x 9" / 11.5 x 23cm
- **1 Anodized Aluminum Extension Pole**  
Size: 2-4ft / 60-122cm

\* Brushes are made with FDA approved polyester bristles

- **1 Rounded Corner Brush**  
Head size: 5.5 x 10" / 14 x 25.4cm



**Our Round Brush now has 3 holes underneath for pole placement!**

#### Our new Hand Brush now has an ergonomic handle!

- Better hand comfort
- Longer bristles
- Handle is easier to clean (no cavity)

- **More flexibility** in regards to the brush placement (parallel or perpendicular)
- Hard-to-reach corners will no longer be an issue thanks to the **new brush placement possibilities** and to the **longer brush length** (now 10" instead of 8.75").

\*Some brushes may still have only one center hole.



### Autoclave Regular Brush Cleaning Kit

Item #	Description
50-CLAVREG	Autoclave REGULAR Brush Cleaning Kit

### Superior Autoclave Cleaning Brush Solution Kit includes:

- **1 Hand Brush**  
Length: 4.5 x 9" / 11.5 x 23cm
- **1 Anodized Aluminum Extension Pole**  
Size: 2-4ft / 60-122cm

- **1 Rounded Corner Brush**  
Head size: 5.5 x 8.75" / 14 x 22.225cm
- **2 BATRIK Full Face Shields**

\* Brushes are made with FDA approved polyester bristles

### Autoclave Superior Brush Cleaning Kit

Item #	Description
50-CLAVSUP	Autoclave SUPERIOR Brush Cleaning Kit

### Replacement Components



**Hand Brush**  
50-35900  
**Packaging:** 1



**Brush Head**  
50-55800  
**Packaging:** 1



**Extension Pole**  
50-EXTPOLE  
**Packaging:** 1



**BATRIK Face Shields**  
SF-100  
**Packaging:** 100 / Box



**BATRIK MEDICAL MANUFACTURING**

850 Halpern Ave., Dorval, QC H9P 1G6

P.: 514-631-7988 • T.-F.: 1-800-661-5432 • F.: 514-631-9083 • E.: [info@surgmed.com](mailto:info@surgmed.com) • W.: [www.surgmed.com](http://www.surgmed.com)