



BATRIK Silicone Pin Mat cushions, protects and organizes delicate surgical instruments for sterilization. Residual moisture collection below the silicone mat is significantly reduced due to our unique **elevated “Nub”** design. These “Nubs” allow for greater air-flow ventilation and drying.

New Features

- now with **more holes** to facilitate better venting and reduce water collection
- unique **elevated nubs** also provide superior venting and a reduction in water collection
- **Custom sizes** available on request (call for details)

Reusable and Autoclavable
Latex free



Silicone Pin Mat

Item #	Dimensions – inches / cm	Packaging
13-2001	13" x 18" / 33cm x 45cm	1 unit
13-2005	13" x 18" / 33cm x 45cm	5 units

www.surgmed.com • info@surgmed.com • +1 514-631-7988 • Fax +1 514-631-9083