


# BATRIK CORNER POCKETS™


With three models to choose from, our Corner Pockets™ portfolio is a truly comprehensive line of products designed to prevent tearing associated with sharp corners or feet on surgical trays.



COLOUR	BLUE	RED	PINK
USAGE	SINGLE-USE	REUSABLE <i>withstands repeated sterilization up to 134°C</i>	REUSABLE <i>withstands repeated sterilization up to 134°C</i>
DIMENSIONS	2.8" x 4" (7.3 cm x 10.3 cm)	2.8" x 4" (7.3 cm x 10.3 cm)	3.4" x 4.7" (8.5 cm x 12 cm)
LATEX-FREE	✓	✓	✓
NO MOISTURE RETENTION	✓	✓	✓
FABRICATION	Flexible, silicone-like plastic	Flexible, silicone-like plastic	Our most flexible plastic for easiest placement and most secure protection

# SOFT POCKETS™

Our Soft Pockets™ Foam Corner Protectors are made from high density foam to effectively cushion and protect the corners of surgical trays to prevent tearing. They are especially useful with heavier instrument sets.



## PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING
<b>18-CP02**</b>	Corner Pockets™ <i>Single-Use, Blue</i>	240 / Bag
<b>18-EVAL2**</b>	Corner Pockets™ <i>Single-Use, Blue</i>	24 / Bag
<b>18-CP07</b>	Corner Pockets™ <i>Reusable, Red</i>	80 / Bag
<b>18-CORNER</b>	Corner Pockets™ <i>Reusable, Red</i>	320 / Box (4 Bags of 80 pcs)
<b>18-EVAL7</b>	Corner Pockets™ <i>Reusable, Red</i>	24 / Bag
<b>18-CP06</b>	Corner Pockets™ <i>Reusable, Pink</i>	80 / Bag
<b>18-CPF-02</b>	Soft Pockets™ <i>Single-Use Foam Corner Protectors</i>	80 / Box (24 Boxes / Pallet)

\*\* Although the single use Corner Pockets™ appear similar to the reusable models they are blended differently and have only been validated as single-use. Re-using them may compromise sterility.