

**BATRIK Instrument Cleaning Reusable Precision Brush** External cleaning brush is designed for cleaning instrument jaws. It is the ideal brush for the quick removal of difficult debris, the BATRIK Reusable Precision Brush can also easily clean deep into instrument jaws. Although the firm and durable bristles feel like stainless steel, there is absolutely no risk of instrument damage or personal injury.

- polyamide (TYNEX®) bristles
- bristles have an extra firm density and superior durability
- steamable thermoplastic handle
- ergonomic handle
- steam / chemical sterilizable
- latex free

No more  
Steel brushes



**ELIMINATES THE NEED FOR THE STAINLESS STEEL TOOTHBRUSH-STYLE BRUSH, WHICH CAN DAMAGE VALUABLE INSTRUMENTS!**



“**THROW AWAY YOUR WIRE BRUSHES, THIS IS THE BEST OPTION TO EXTEND YOUR INSTRUMENT USEFUL LIFE**”

#### Reusable Precision Brush

Item #	Description	Packaging
50-40800	Reusable Precision External Cleaning Brush Overall length is 6.29" / 160mm	3 units / Bag