


BATRIK Instrument Cleaning ANTIMICROBIAL Acetabular Reamer Brush - Durable and Reusable Cleaning Brush. This double-ended cleaning brush is manufactured exclusively for cleaning acetabular bone reamers. The "Pig Tail" rounded style circular head brush is designed to effectively clean inside the reamer, while the rigid straight brush on the reverse side is used for cleaning individual cutting holes. Special fusion of the brush heads to the hand piece along with using high quality **antimicrobial** medical grade bristles ensures a durable and reusable effective cleaning tool.

Product Specifications*

- overall length: 8" (200mm)
- round brush head diameter: 1.75" (45mm)
- straight brush diameter: 5/16" (8mm) Tapered
- straight brush length: 1 3/16" (30mm)

*specifications subject to small variances



Brush ±
1.75" / 45mm

Brush Length
8" / 200mm

Tapered Brush ±
0.3125" / 8mm

Acetabular Reamer Cleaning Brush

Item #	Description	Packaging
50-ACETAB	Acetabular Reamer Cleaning Brush	3 units / Bag