

BATRIK INSTRUMENT CLEANING All Purpose Brushes

Plastic Handle Brushes

We offer a wide selection of brushes for most of your instrument cleaning applications. These brushes are ideal for large bore, uniquely shaped and endoscopic instruments.

In-Brush™ Plastic Handle

Choose either the regular toothbrush style (nylon or stainless) or try our double-ended brush.

*order in volume and save even more



Plastic Single End Nylon – Handle Brush

| Item # | Description | Length Inches / mm | Brush Length Inches / mm | Brush (Ø) Inches / mm |
|----------|------------------|-----------------------|-----------------------------|--------------------------|
| 28-72550 | Single End Nylon | 7.25" / 184mm | 1.5" / 38mm | 0.5" / 12.7mm |

Packaging: Bag of 3 - ITEM # 28-72550-50 / Bag of 50



Plastic Single End Stainless Steel – Handle Brush

| Item # | Description | Length Inches / mm | Brush Length Inches / mm | Brush (Ø) Inches / mm |
|----------|----------------------------|-----------------------|-----------------------------|--------------------------|
| 29-72550 | Single End Stainless Steel | 7.25" / 184mm | 1.5" / 38mm | 0.5" / 12.7mm |

Packaging: Bag of 3 - ITEM # 29-72550-50 / Bag of 50



Plastic Double End Nylon – Handle Brush

| Item # | Description | Length Inches / mm | Brush Length Inches / mm | Brush (Ø) Inches / mm |
|----------|-----------------|-----------------------|-----------------------------|--------------------------|
| 28-69394 | Large Brush-end | 6.9" / 175mm | 1.38" / 35mm | 0.39" / 10mm |
| | Small Brush-end | 6.9" / 175mm | 1" / 25mm | 0.2" / 5mm |

Packaging: Bag of 3 - ITEM # 28-69394-50 / Bag of 50



Hand Scrub Nylon Brush – Non-Autoclavable

| Item # | Length Inches / mm | Brush Length Inches / mm | Brush (Ø) Inches / mm | Packaging |
|-------------|-----------------------|-----------------------------|--------------------------|-----------|
| 50-40911 | 4.09" / 104mm | 3.74" / 95mm | 1.1" / 28mm | Bag of 3 |
| 50-40911-50 | 4.09" / 104mm | 3.74" / 95mm | 1.1" / 28mm | Bag of 50 |