

PAD SIZE

Our pads are smaller than competitor pads, allowing for easier placement.

CONDUCTOR PLATE

Our aggressive hydrogel reduces the need for a large foam boarder allowing for a larger conductor plate than the leading competitors.

INSTRUCTIONS FOR USE

Checklist is featured on the opaque release liner.

Instructions for use:

1. Remove electrodes from peel pouch.
2. Ensure skin site is clean and dry.
3. Remove electrodes from release liner.
4. Apply electrodes firmly to patient's skin.
5. Connect electrodes to defibrillator or patient cable.
6. Connect ECG leads for demand pacing.

OPAQUE LINER

Prevents slipping accidents.

60" LEAD WIRES

Longest in the market. Provides more flexibility when placing the defibrillator.

COLOR CODED PADS

Anterior electrode is white foam and the posterior electrode is blue foam. This helps train staff to look for colors to aid in pad placement.

COLOR CODED PACKAGING

Vermed Defib Pads have a 36 month shelf life (3 years).
Approximately 50% longer than the leading competitors.

LEADS OUT OPTION

Allows pads to be pre-connected prior to opening pouch.

Latex-free



PHYSIO CONTROL



PHILIPS MEDICAL



CARDIAC SCIENCE



ZOLL MEDICAL

1.800.669.6905 | Vermed.com

Patient Monitoring & Medical Surgical Consumables

Defibrillation Electrodes

Defib Pads for all departments of the hospital

Vermed's Defib Pads were designed with your needs in mind. They are not just different - they are a step above the other brands.

PHILIPS MEDICAL

Product #	For Defibrillator Models	Leads	Pad & Leadwires
D311	XLT, MRX, HeartStart XL, MRX ALS, HP Codemaster XL, HeartStart FR2+, HeartStart Forerunner	In	Radiotransparent - Invisible
D309	XLT, MRX, HeartStart XL, MRX ALS, HP Codemaster XL, HeartStart FR2+, HeartStart Forerunner	In	Radiolucent - Slight Shadow
D317	XLT, MRX, HeartStart XL, MRX ALS, HP Codemaster XL	Out	Radiolucent - Slight Shadow Leads Out Package
D321	XLT, MRX, HeartStart XL, MRX ALS, HP Codemaster XL	In	Pediatric Radiolucent - Slight Shadow

PHYSIO CONTROL

Product #	For Defibrillator Models	Leads	Pad & Leadwires
D304	Lifepak 9P, 15, 10C, 20, 12, 500, 500 DPS, 1000	In	Radiotransparent - Invisible
D313	Lifepak 9P, 15, 10C, 20, 12, 500, 500 DPS, 1000	In	Radiolucent - Slight Shadow
D315	Lifepak 9P, 15, 10C, 20, 12, 500, 500 DPS, 1000	Out	Radiolucent - Slight Shadow Leads Out Package
D319	Lifepak 9P, 15, 10C, 20, 12	In	Pediatric Radiolucent - Slight Shadow


ZOLL MEDICAL

Product #	For Defibrillator Models	Leads	Pad & Leadwires
D303	M-Series, E-Series, R-Series, X-Series	In	Radiotransparent - Invisible
D312	M-Series, E-Series, R-Series, X-Series	In	Radiolucent - Slight Shadow
D314	M-Series, E-Series, R-Series, X-Series	Out	Radiolucent - Slight Shadow Leads Out Package
D318	M-Series, E-Series, R-Series, X-Series	In	Pediatric Radiolucent - Slight Shadow


CARDIAC SCIENCE

Product #	For Defibrillator Models	Leads	Pad & Leadwires
D326	PowerHeart AED G3 Plus, PowerHeart AED G3 Pro	Out	Can be used in place of Cardiac Science Part #9131-001


PHILIPS MEDICAL
Connector Style



PHYSIO CONTROL
Connector Style



ZOLL MEDICAL
Connector Style



CARDIAC SCIENCE
Connector Style

