




## HEARTSYNC

Defibrillation • Cardioversion • Pacing • ECG Monitoring

- 
**36 MONTH SHELF LIFE**  
 Extended shelf life reduces waste and saves money
- 
**60 INCH LEAD WIRES**  
 Increased flexibility with equipment placement & clinician mobility
- 
**ADVANCED HYDROGEL**  
 Unique gel provides maximum adhesion
- 
**SMALLER PAD SIZE**  
 Allows for easier placement without crowding
- 
**LARGER CONDUCTOR PLATE**  
 Larger conductive area disperses energy more effectively
- 
**IMAGING VISIBILITY**  
 Radiolucent & radiotransparent options to suit all your needs
- 
**VISIBLE LINER**  
 Prevents slip hazard to clinician

### WEAR TIME

24 hours

### INSTRUCTIONS FOR USE

User checklist on the release liner

### LEADS OUT OPTION

Pads can be pre-connected prior to opening

### COLOUR CODED PADS

Anterior electrode is blue foam and the posterior electrode is white foam. This helps to train staff look for colours, therefore aiding in pad placement.

### COLOUR CODED PACKAGING

The correct product, for your device.


Latex-free




**PHYSIO CONTROL**



**PHILIPS MEDICAL**



**CARDIAC SCIENCE**



**ZOLL MEDICAL**

**+44 1803 860100 | Vermed.co.uk**


Patient Monitoring & Medical Surgical Consumables

# Defibrillation Electrodes

## Defib pads for all departments of the hospital



Vermed's HeartSync Defibrillation Electrodes are designed with your needs in mind. Featuring long leadwires, advanced hydrogel and ease of placement, they will meet the needs of your clinical staff.



### PHYSIO CONTROL

Part #	Product Reference	For Defibrillator Models	Leads	Pad & Leadwires
6600017H	C100AC-PHYSIO / D304	Lifepak 9P, 15, 10C, 20, 12, 500, 500 DPS, 1000	In	
6600093H	T100AC-PHYSIO / D313	Lifepak 9P, 15, 10C, 20, 12, 500, 500 DPS, 1000	In	
6600100H	T100LOAC-PHYSIO / D315	Lifepak 9P, 15, 10C, 20, 12, 500, 500 DPS, 1000	Out	
6600068H	PEDIATRIC-PHYSIO / D319	Lifepak 9P, 15, 10C, 20, 12	In	

### PHILIPS MEDICAL

Part #	Product Reference	For Defibrillator Models	Leads	Pad & Leadwires
6600015H	C100AC-PHILIPS / D311	XLT, MRX, HeartStart XL, MRX ALS, HP Codemaster XL, HeartStart FR2+, HeartStart Forerunner	In	
6600091H	T100AC-PHILIPS / D309	XLT, MRX, HeartStart XL, MRX ALS, HP Codemaster XL, HeartStart FR2+, HeartStart Forerunner	In	
6600099H	T100LOAC-PHILIPS / D317	XLT, MRX, HeartStart XL, MRX ALS, HP Codemaster XL	Out	
6600066H	PEDIATRIC-PHILIPS / D321	XLT, MRX, HeartStart XL, MRX ALS, HP Codemaster XL	In	

### CARDIAC SCIENCE


Part #	Product Reference	For Defibrillator Models	Leads	Pad & Leadwires
6600088H	T100-CS / D326	PowerHeart AED G3 Plus, PowerHeart AED G3 Pro	Out	 Can be used in place of Cardiac Science Part #9131-001

### ZOLL MEDICAL


Part #	Product Reference	For Defibrillator Models	Leads	Pad & Leadwires
6600019H	C100AC-ZOLL / D303	M-Series, E-Series, R-Series, X-Series	In	
6600095H	T100AC-ZOLL / D312	M-Series, E-Series, R-Series, X-Series	In	
6600102H	T100LOAC-ZOLL / D314	M-Series, E-Series, R-Series, X-Series	Out	
6600070H	PEDIATRIC-ZOLL / D318	M-Series, E-Series, R-Series, X-Series	In	

For additional information, defibrillator compatibility and usage visit: [vermed.com/defib](http://vermed.com/defib)


**PHYSIO CONTROL**  
Connector Style



**PHILIPS MEDICAL**  
Connector Style



**CARDIAC SCIENCE**  
Connector Style



**ZOLL MEDICAL**  
Connector Style

