



VALVE SETS



Valves and adapters in a practical disposable set

- ✓ Olympus and Fujifilm
- ✓ High safety
- ✓ Constant performance






Valve sets for endoscope

In a practical disposable set



OLYMPUS



Features

Biopsy valve, air/water valve, suction valve

- ☐ Suitable for Olympus OES 40, EVIS 140, 145, 160, 165, 180, 185, 190, 240, 260, 290
- ☐ Minimises cross-contamination risks
- ☐ Supports constant performance

Jet channel adapter

- ☐ Used together with the rinsing hose in connection with a rinsing pump for flushing with rinsing fluids
- ☐ With integrated non-return valve
- ☐ Individually sterile packed

Set contents

Pieces per box

Manufacturer Product Code

Biopsy valve, air/water valve, suction valve & jet channel adapter	25	BP-41650
Biopsy valve, air/water valve, suction valve	25	BP-41651
Air/water valve, suction valve	25	BP-41652
Biopsy valve, air/water valve, suction valve, air/water cleaning valve	25	BP-41653
Biopsy valve, air/water valve, suction valve, air/water cleaning valve & jet channel adapter	25	BP-41654
Air/water valve, suction valve, air/water cleaning valve & jet channel adapter	25	BP-41655
Air/water cleaning valve	100	BP-41656

FUJIFILM



Features

- ☐ Suitable for Fujifilm series
- ☐ Minimises cross-contamination risks
- ☐ Supports constant performance

Set-Inhalt

Stück je Box

Manufacturer Product Code

Series 700		
Biopsy valve, air/water valve, suction valve	25	? AWV-KIT-01-FJ
Air/water valve, suction valve	25	? AWV-KIT-02-FJ
Air/water cleaning valve	100	? AWV-KIT-07-FJ
Series 500/600		
Biopsy valve, air/water valve, suction valve	25	BP-40060050
Air/water valve, suction valve	25	BP-40060051

