

## LED x4 MAGNIFIER



## LED 4x MAGNIFIER

Visual inspection of instruments can be difficult, the LED 4x magnifier is fitted with a 72mm diameter lens for up-close visual inspection of potential damage and residual dirt after cleaning.

The lightweight 85x85x60mm magnifier comes with four white LED lights featuring three levels of brightness, 4x the magnification for easy and thorough inspection of robotic instrumentation, as well as an on/off switch and a Micro-USB charging cable.

The magnifier is available for individual purchase.

Code	Description
MAG-002	LED 4x Magnifier

**MAGIC TOUCH MAGNIFIER**



**MAGIC TOUCH MAGNIFIER**

Designed for detailed inspection, the 9.5 x 9.5 x 4.5 cm Magic Touch has a magnification of 3X and weighs 204grams. Equipped with three LED lights, the Magic Touch magnifier is a touch activated unit that features an on/off switch for uninterrupted lighting. The optional stand is a perfect accessory to keep your hands free.

Code	Description	Dimensions	Quantity / Pack
RNDMAG-1	Magic touch magnifier	ø 9.53 x 4.57h	1
RNDMAGSTND-1	Stand for Magic Touch Magnifier	/	

**OPTIVISOR**



**OPTIVISOR**

Soft pliable wide headband for ultimate comfort. Easily adjust to any head size. The self-locking design remains in place. Glass lens ensures high quality image. Magnification is 3.5X. Working distance is 10 cms and weight 272 grams.

Code	Description	Quantity / Pack
26106	Optivisor	1