



CLEANING BRUSHES



Safe and gentle cleaning

- ✓ Rounded brush tips
- ✓ Resistant to bending
- ✓ Latex-free





For flexible Endoscopes



Cleaning brushes

Safe and gentle cleaning



Features

- Ball heads on the brush ends and rounded brush tips for gentle cleaning
- Brushes made of soft, durable nylon
- Ideal brush head size, with optimal density
- Brush-head diameter allows for maximum mechanical brushing of the endoscope channel
- Bend resistant thanks for first-class catheter rigidity

SINGLE/DUAL-ENDED CLEANING BRUSHES

Colour	Catheter		Brush 1		Brush 2		Pieces per box	Manufacturer Product Code
	Ø	Length	Ø	Length	Ø	Length		
For endoscopes with 2.8 – 4.2 mm working channels								
Yellow	1.7 mm	230 cm	5 mm	20 mm	-	-	100 / 400	BP-41450
Yellow	1.7 mm	180 cm	5 mm	20 mm	5 mm	20 mm	100 / 400	BP-41454
Yellow	1.7 mm	230 cm	6 mm	20 mm	-	-	100	BP-41469
Blue	1.7 mm	230 cm	6 mm	20 mm	6 mm	20 mm	100	BP-41455
White	1.7 mm	230 cm	5 mm	20 mm	5 mm	20 mm	100 / 400	BP-40392
Yellow	1.8 mm	230 cm	5 mm	20 mm	5 mm	20 mm	100	BP-41456
For bronoscopes with working channels from 2.0 mm								
Green	1.5 mm	120 cm	3 mm	20 mm	3 mm	20 mm	100	BP-41470
Green	1.5 mm	180 cm	3 mm	20 mm	3 mm	20 mm	100	BP-41471
Green	1.5 mm	230 cm	3 mm	20 mm	3 mm	20 mm	100 / 400	BP-41453
For children's endoscopes with working channels from 2.0 mm								
Green	1.5 mm	230 cm	3 mm	20 mm	3 mm	20 mm	100	DCBS-18-20-03-230?
For children's bronoscopes with working channels from 1.6 mm								
Yellow	1.0 mm	120 cm	2 mm	15 mm	2 mm	15 mm	100	BP-41473
For enteroscopes with working channels from 2.8 mm								
Blue	1.7 mm	300 cm	6 mm	20 mm	6 mm	20 mm	100	BP-41472
For bronoscopes with working channels from 1.2 mm								
Yellow	1.0 mm	100 cm	2 mm	15 mm	-	-	100 / 400	BP-41449
Blue	1.0 mm	130 cm	1.5 mm	10 mm	-	-	100 / 400	BP-40399
Blue	1.0 mm	130 cm	1.5 mm	10 mm	1.5 mm	10 mm	100 / 400	BP-41451
Blue	1.0 mm	230 cm	2 mm	15 mm	2 mm	15 mm	100 / 400	BP-41452
For endoscopes with working channels from 1.0 mm								
Blue	0.8 mm	120 cm	1.5 mm	10 mm	1.5 mm	10 mm	100	BP-40030420
Blue	0.8 mm	180 cm	1.5 mm	10 mm	1.5 mm	10 mm	100	BP-40030440
Blue	0.8 mm	230 cm	1.5 mm	10 mm	1.5 mm	10 mm	100	BP-40030421

Features

- Efficient cleaning thanks to two brushes on one shaft





For flexible endoscopes



Cleaning brushes


Safe and gentle cleaning



Features

- Ball heads on the brush ends and rounded brush tips for gentle cleaning
- Brushes made of soft, durable nylon
- Ideal brush head size and density
- Brush-head diameter allows for maximum mechanical brushing of the endoscope channel
- Bend resistant thanks for first-class catheter rigidity

MULTIPOINT DUAL-ENDED CLEANING BRUSH



Features

- Combines the advantages of a dual-ended cleaning brush with an efficient sweeper

Colour	Catheter		Brush 1		Brush 2		Brush 3		Pieces per box	Manufacturer Product Code
	Ø	Length	Ø	Length	Ø	Length	Ø	Length		
For endoscopes with 2.8 – 4.2 mm working channels										
Blue	1.7 mm	230 cm	6 mm	20 mm	6 mm	20 mm	6 mm	20 mm	100	BP-40393
Blue	1.7 mm	260 cm	6 mm	20 mm	6 mm	20 mm	6 mm	20 mm	100	MCBS-23-20-06-260

