




# COMPLETE PRODUCT CATALOGUE

Established in 1997. Proudly Canadian. 🇨🇦  
**Sold in over 80 countries worldwide.**



Proudly Canadian and established in 1997, we specialize in the manufacturing of enzymatic solutions offered in different formats as well as the fabrication of silicone and plastic products. Our products are sold in over 80 countries worldwide and used primarily in surgery, endoscopy, and central sterilization areas.

Our in-house controlled environment production facility in Montreal, Canada adheres to strictly enforced quality and regulatory control systems. We are registered as a device manufacturing establishment with the United States FDA and we are in full compliance with CE and Health Canada Regulations. We continually invest in upgrading our plant, machinery, and processes to ensure high quality products at high volumes. Most products are available off-the-shelf.

At Batrik, we are committed to innovation and providing our customers with the finest quality products backed by world-class customer care.

## SURGICAL PRODUCTS & GENERAL ACCESSORIES






- 4 SUTURE CLAMP TAGS
- 5 CLAMP JAW COVERS
- 6 VASCULAR SILICONE TIES
- 7 CAUTERY TIP CLEANER
- 8 SKINS™ [Surgical Skin Multi-Tip Markers](#)
- 9 PIN GUARDS
- 10 SAFE-SOAK™ SPAC & MAT
- 11 MAGNETIC MAT
- 12 TOURNY™+ TOURNIQUETS
- 14 SAVEFACE™ SHOE COVERS
- 15 STACKABLE STEP STOOLS
- 16 ANTI-FATIGUE MATS  
[Blue Cloud-Gel™ / Blue Cloud-Nimbus™](#)

## ENDOSCOPY PRODUCTS

- 18 GOLFF™ ANTI-FOG SOLUTION
- 19 BITE BLOCKS  
[Endoscopy BO2™ / Bitestrong™](#)
- 21 COUVRETTE™ VALVE COVERS
- 22 DAYBRUSH™ CHANNEL BRUSHES
- 23 ARC™ CHANNEL BRUSHES
- 26 VALVE / CONTROL BRUSHESW
- 27 V-BRUSHES

## INFECTION CONTROL & STERILIZATION PRODUCTS

- 28 SILICONE PIN MAT
- 29 CORNER POCKETS™ & SOFT POCKETS™
- 30 INSTRUMENT PROTECTORS
- 33 S.T.I.C.S.
- 34 ID TAPE
- 35 TAPE-B-GONE™ TAPE REMOVER
- 36 PADLOCKS SECURITY TAGS
- 40 FACE SHIELDS
- 42 SCRUB PAD
- 43 ENZYMATIC DETERGENTS

- 44 INSTRUMENT LUBRICANT
- 45 **Your First** ™ ENVIRO-SACHET ENZYMATIC GEL
- 46 **Your First** ™ FOAMING SPRAY
- 47 **Your First** ™ F1RST WIPES
- 48 **Your First** ™ BEDSIDE SPONGE KITS
- 51 **Your First** ™ TRANSPORT BAGS

## BRUSHES

- 52 ARC™ WIRE STYLE BRUSHES
- 56 DOUBLE ARC™ WIRE BRUSHES  
ASYMMETRIC / SYMMETRIC
- 59 IN-BRUSH™ TWISTED WIRE BRUSHES
- 61 STAINLESS STEEL ARC™ MICRO BRUSHES
- 62 INSTRUMENT CLEANING FLAT BRUSHES
- 63 ORTHOPEDIC RASPS CLEANING BRUSH
- 64 PRECISION BRUSH
- 65 JAR /BOTTLE & DRAIN BRUSHES
- 66 ACETABULAR REAMER BRUSH
- 67 ALL PURPOSE BRUSHES
- 68 AUTOCLAVE CLEANING BRUSH KITS

## WOMEN'S HEALTH



- 69 HYSTEROSONOGRAPHY CATHETER

# BATRIK SUTURE CLAMP TAGS

Non-toxic, radiopaque, and latex free, our Suture Clamp Tags are designed to protect and safely hold sutures without any damage.

Additionally, its foam base, affixed to a sterile area with the clamps, allows for easy insertion.

- Packed in a convenient foam base with easy peel-off adhesive backing, even while wearing gloves
- Available in sterile and non-sterile
- Available in 2, 3, or 5 pairs
- Pre-counted



OUR "EASY PEEL CORNER" IS A FEATURE UNIQUE TO BATRIK

## PRODUCT SELECTION

ITEM # STERILE	ITEM# NON-STERILE	COLOR	DIMENSIONS	PACKAGING
<b>11-1002</b>	<b>911-1002</b>	Yellow	Standard ø (0.06" x 0.5" / 0.15 cm x 1.27 cm)	2 pairs / foam base
<b>11-1003</b>	<b>911-1003</b>	Yellow	Standard ø (0.06" x 0.5" L / 0.15 cm x 1.27 cm)	3 pairs / foam base
<b>11-1005</b>	<b>911-1005</b>	Yellow	Standard ø (0.06" x 0.5" L / 0.15 cm x 1.27 cm)	5 pairs / foam base
<b>11-1025</b>	<b>911-1025</b>	Assorted	Standard ø (0.06" x 0.5" L / 0.15 cm x 1.27 cm)	5 pairs / foam base
<b>11-1045</b>	<b>911-1045</b>	White	Standard ø (0.06" x 0.5" L / 0.15 cm x 1.27 cm)	5 pairs / foam base
<b>11-1065</b>	<b>911-1065</b>	Blue	Standard ø (0.06" x 0.5" L / 0.15 cm x 1.27 cm)	5 pairs / foam base
<b>11-1085</b>	<b>911-1085</b>	Red	Standard ø (0.06" x 0.5" L / 0.15 cm x 1.27 cm)	5 pairs / foam base
<b>11-3002</b>	<b>911-3002</b>	Pink	Large ø (0.09" x 0.75" L / 0.23 cm x 1.91 cm)	2 pairs / foam base
<b>11-3003</b>	<b>911-3003</b>	Pink	Large ø (0.09" x 0.75" L / 0.23 cm x 1.91 cm)	3 pairs / foam base
<b>11-3005</b>	<b>911-3005</b>	Pink	Large ø (0.09" x 0.75" L / 0.23 cm x 1.91 cm)	5 pairs / foam base

STERILE PACKAGING: 1 Foam Base / Pouch – 5 Pouches / Box

NON-STERILE PACKAGING: Sold per foam base



## BATRIK CLAMP JAW COVERS

Manufactured from medical grade silicone, our Clamp Jaw Covers reduce trauma to vessels requiring occlusion, particularly when a hemostatic clamp is utilized. The thread-like indents allow for gentle, constant and non-slip tension.

GOOD TO KNOW: OUR CLAMP JAW COVERS CONVERT ALMOST ANY INSTRUMENT TO ATRAUMATIC



- Gamma sterilized

- Available in 3 colors and 3 widths to fit most clamp sizes

- Radiopaque and latex free

### PRODUCT SELECTION

ITEM # STERILE	ITEM # NON-STERILE	DIMENSIONS	COLOR
<b>12-1000</b>	<b>912-1000</b>	0.1" X 5" / 0.25 cm X 12.7 cm	Red
<b>12-1001</b>	<b>912-1001</b>	0.2" X 5" / 0.51 cm X 12.7 cm	White
<b>12-1002</b>	<b>912-1002</b>	0.3" X 5" / 0.76 cm X 12.7 cm	Blue

STERILE PACKAGING: 2 clamps per sterile pouch - 10 pouches per box

NON-STERILE PACKAGING: 20 clamps per bag (also available 100 units / bag)

# BATRIK VASCULAR SILICONE TIES

AVAILABLE IN 4 SIZES  
AND 4 COLORS



Our Vascular Silicone Ties permit safe and smooth retraction, occlusion and identification of vessels during surgical procedures.

Unlike clamps and retractors which can pinch the vessels flat, or even crack them, the BATRIK elliptical shape allows for the gentle handling of arteries, veins, nerves, tendons, spermatic cords, ureters, or whenever delicate retraction / occlusion is required.

## PRODUCT SELECTION

ITEM #	DIMENSIONS	COLOR
<b>10-1000</b>	0.075" X 16" 0.2 cm x 40 cm	Mini-Red
<b>10-1001</b>	0.075" X 16" 0.2 cm x 40 cm	Mini-Blue
<b>10-1002</b>	0.075" X 16" 0.2 cm x 40 cm	Mini-Yellow
<b>10-1003</b>	0.075" X 16" 0.2 cm x 40 cm	Mini-White
<b>10-1010</b>	0.120" x 16" 0.3 cm x 40 cm	Maxi-Red
<b>10-1011</b>	0.120" x 16" 0.3 cm x 40 cm	Maxi-Blue
<b>10-1012</b>	0.120" x 16" 0.3 cm x 40 cm	Maxi-Yellow
<b>10-1013</b>	0.120" x 16" 0.3 cm x 40 cm	Maxi-White
<b>10-1101</b>	0.040" x 16" 0.1 cm x 40 cm	Extra-Mini-Blue
<b>10-1111</b>	0.200" x 22" 0.5 cm x 56 cm	Extra-Maxi-Blue
<b>10-1113</b>	0.200" x 22" 0.5 cm x 56 cm	Extra-Maxi-White

PACKAGING: 2 ties per sterile pouch – 10 pouches per box



- Manufactured of medical grade silicone



- Do not twist on themselves when unpacked



- Easy to separate, eliminating the need the need for cardboard inserts



- Sterile, radiopaque, and latex free

**BATRIK**  
**CAUTERY TIP CLEANER**

Our Cautery Tip Cleaner assists in removing material from cautery tips for both mono-polar and bi-polar electrodes.

- Foam pad with integrated abrasive material
- Adheres to disposable and non-disposable drapes
- Radiopaque
- Available exclusively in Non-Sterile
- Latex free

**PRODUCT SELECTION**

ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
<b>60-10000</b>	Cautery Tip Cleaner Non-Sterile	1.8" x 1.8" / 4.5 cm x 4.5 cm	2,200 /Box

# BATRIK SKINS™ SURGICAL MULTI-TIP SKIN MARKERS

Our SKINS™ Surgical Multi-Tip Skin Markers are designed with a beveled tip to enable the creation of fine or bold tracing lines by holding the pen at different angles.

- Available with or without 6 labels and a flexible 6" calibrated (in/mm) ruler
- Available in non-sterile for pre-operative use, or sterile for use in the operating room
- Surgical prep resistant
- Made of lasting gentian violet ink that is non-toxic and non-tattooing
- Latex free



## PRODUCT SELECTION


ITEM #	TYPE	DESCRIPTION	PACKAGING
<b>S-3000</b>	Sterile	Multi-Tip fine marker with ruler and label set	50 / Box
<b>S-3001</b>	Sterile	Multi-Tip fine marker	50 / Box
<b>S-30005</b>	Sterile	Multi-Tip fine marker with ruler and label set (Volume Packaging)	500 / Box
<b>S-30015</b>	Sterile	Multi-Tip fine marker (Volume Packaging)	500 / Box
<b>S-93000</b>	Non-Sterile	Multi-Tip fine marker with ruler and label set	1000 / Box
<b>S-93001</b>	Non-Sterile	Multi-Tip fine marker	1000 / Box
<b>S-94000</b>	Non-Sterile	Ruler	100 / Bag
<b>S-94001</b>	Non-Sterile	Ruler (Bulk)	20,000 / Box
<b>S-95000</b>	Non-Sterile	Set of 3 labels	100 / Bag

SAVE by ordering our SKINS™ MARKERS in volume packaging. Cost effective and reduces waste.

# BATRIK PIN GUARDS



- Safely cover the sharp ends of your K-Wire and Steinmann Pins
- Color coded by wire size which eliminates the risk of selecting of the wrong size
- A removable sleeve is provided to avoid injury and facilitate application (sterile only)
- Non-sterile model is available with or without plastic sleeves
- Packaging designed to reduce waste
- Latex free



## PRODUCT SELECTION

ITEM # PIN GUARDS + SLEEVE	ITEM # PIN GUARDS + SLEEVE	ITEM # PIN GUARDS (NO SLEEVE)	COLOR	DIMENSIONS
<b>14-1028</b>	<b>914-1028 - SLV</b>	<b>914-1028</b>	Yellow	0.028" / 0.71 mm
<b>14-1035</b>	<b>914-1035 - SLV</b>	<b>914-1035</b>	Light Blue	0.035" / 0.89 mm
<b>14-1045</b>	<b>914-1045 - SLV</b>	<b>914-1045</b>	White	0.045" / 1.14 mm
<b>14-1054</b>	<b>914-1054 - SLV</b>	<b>914-1054</b>	Dark Red	0.054" / 1.40 mm
<b>14-1062</b>	<b>914-1062 - SLV</b>	<b>914-1062</b>	Green	0.062" / 1.57 mm
<b>14-1078</b>	<b>914-1078 - SLV</b>	<b>914-1078</b>	Pink	0.078" / 1.98 mm
<b>14-1093</b>	<b>914-1093 - SLV</b>	<b>914-1093</b>	Grey	0.093" / 2.36 mm
<b>14-1109</b>	<b>914-1109 - SLV</b>	<b>914-1109</b>	Beige	0.109" / 2.77 mm
<b>14-1125</b>	<b>914-1125 - SLV</b>	<b>914-1125</b>	Purple	0.125" / 3.17 mm
<b>14-1141</b>	<b>914-1141 - SLV</b>	<b>914-1141</b>	Brown	0.141" / 3.58 mm
<b>14-1156</b>	<b>914-1156 - SLV</b>	<b>914-1156</b>	Dark Blue	0.156" / 3.96 mm
<b>14-1188</b>	<b>914-1188 - SLV</b>	<b>914-1188</b>	Red	0.188" / 4.77 mm

STERILE PACKAGING: 80 units/box (2 units per sterile pouch - 40 pouches per box)

NON-STERILE PACKAGING: 100 units/box



## BATRIK SAFE-SOAK™ SPAC

Our Safe-Soak™ SPAC floor suction disk quickly removes fluid waste from the operating room floor through the application of suction.



- Highly durable: it won't be damaged if stepped on
- Easily maneuverable by foot around operating room floor
- Latex free, disposable, and non-sterile
- The bottom opening permits fluid collection through a highly absorbent material that comes into contact with its surface
- The hose is securely affixed on its side to reduce risk of catching equipment
- Its highly absorbent base material allows for quick fluid wicking and removal when suction is applied

## BATRIK SAFE-SOAK™ ABSORBENT FLOOR MATS

Our Safe-Soak™ Absorbent Floor Mats are designed to absorb a variety of fluids from the operating room floor quickly and efficiently, ensuring the floor beneath remains slip-free, safe, and dry during surgical procedures.

- Highly absorbent
- Facilitates quick and easy clean up
- Latex free, single-use, disposable, and non-sterile
- Conforms to uneven surfaces
- Skid-resistant backing (blue absorbent mat only)




### PRODUCT SELECTION

ITEM #	DESCRIPTION	ABSORBENCY	SIZE	PACKAGING
<b>SS-20-1000</b>	Safe-Soak™ / White extra absorbent mat	4.5 L / mat	28" x 48" (71 cm x 122 cm)	20 / Box
<b>SS-20-2000</b>	Safe-Soak™ / Blue absorbent mat with skid-resistant backer	2.5 L / mat	28" x 40" (71 cm x 102 cm)	20 / Box
<b>SS-20-3000</b>	Safe-Soak™ / Yellow absorbent mat	3 L / mat	30" x 52" (76 cm x 132 cm)	15 / Box
<b>SS-21-1000</b>	Safe-Soak™ SPAC	N/A	8.625" (22 cm diameter)	12 / Box



# BATRIK MAGNETIC MAT

Our Magnetic Mat provides a sterile, neutral field for the safe placement and hands-free transfer of instruments. It is latex free, reusable, and autoclavable (up to 134° C).




## PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
<b>MAG-02-1</b>	Magnetic Mat	12" x 15" / 30 cm x 38 cm	1 Unit
<b>MAG-02-5</b>	Magnetic Mat (Volume Packaging)	12" x 15" / 30 cm x 38 cm	5 Units

BATRIK

# TOURNY™ + TOURNIQUETS

Our Tourniquets are comfortable and easy to use on all patients, including the elderly, undernourished, children, and newborns.



## AVAILABLE COLORS



- Available in different sizes including compact sizes, which are convenient for phlebotomy trays
- Innovative latex free material has a more comfortable texture for patients, puts less stress on skin, and safeguards hospital staff and patients from latex allergens

## SINGLE-USE BENEFITS

- Reduced risk of cross-contamination (patient to patient or patient to healthcare provider)
- Reduced risk of spread of nosocomial infections
- May improve your facility's compliance with OSHA's 1910-1030 (d)(1) "Blood-borne Pathogens Standard"

FOR YOUR CONVENIENCE:  
 25 UNITS PER ROLL  
 COME IN A DISPENSER BOX  
 FOR EASY ACCESS



## PRODUCT SELECTION

ITEM #	DESCRIPTION	COLOR	DIMENSIONS	PACKAGING
<b>K-002</b>	TOURNY™ + Non-Sterile	Blue	1" x 18" / 2.5 cm x 45 cm	1 / Bag
<b>K-004</b>	TOURNY™ + Non-Sterile	Blue	1" x 18" / 2.5 cm x 45 cm	100 / Box
<b>K-009</b>	TOURNY™ + Non-Sterile	Blue	0.200" x 22" / 0.5 cm x 56 cm	20 / Box
<b>K-010</b>	TOURNY™ + Non-Sterile	White	0.200" x 22" / 0.5 cm x 56 cm	20 / Box
<b>K-400</b>	TOURNY™ + Non-Sterile	Blue	1" x 18" / 2.5 cm x 45 cm	16 rolls 400 / Box
<b>K-400PB</b>	TOURNY™ + Non-Sterile (Pebbled)	Light Blue	1" x 18" / 2.5 cm x 45 cm	16 rolls 400 / Box
<b>K-100</b>	TOURNY™ + Non-Sterile	Assorted (3 rolls of each) Blue, Orange, Green, and Pink	1" x 18" / 2.5 cm x 45 cm	12 rolls 300 / Box
<b>K-254</b>	TOURNY™ + Non-Sterile (Pebbled)	Light Blue	1" x 18" / 2.5 cm x 45 cm	250 / Box
<b>K-2-18RL</b>	TOURNY™ + Non-Sterile	Blue	1" x 18" / 2.5 cm x 45 cm	1 Roll
<b>K-2-18PBRL</b>	TOURNY™ + Non-Sterile (Pebbled)	Light Blue	1" x 18" / 2.5 cm x 45 cm	1 Roll
<b>K-3-18RL</b>	TOURNY™ + Non-Sterile	Green	1" x 18" / 2.5 cm x 45 cm	1 Roll
<b>K-6-18RL</b>	TOURNY™ + Non-Sterile	Pink	1" x 18" / 2.5 cm x 45 cm	1 Roll
<b>K-7-18RL</b>	TOURNY™ + Non-Sterile	Orange	1" x 18" / 2.5 cm x 45 cm	1 Roll
<b>K-2-29RL</b>	Long TOURNY™ + Non-Sterile	Blue	1" x 29.5" / 2.5 cm x 75 cm	1 Roll

BATRIK

# SAVEFACE™ SHOE COVERS

Our SaveFace™ superior medical shoe covers are durable, safe, and easy to use. Made of polypropylene fabric with robust welds, our shoe covers offer the ultimate in durability. The triple-elastic run and anti-skid print allows for a smooth fit around the foot and reduces the risk of slipping.

- Particle count tested and complies with IEST-RP-CC003.4
- 100% latex and silicone free



MADE IN CANADA 

## PRODUCT SELECTION

ITEM #	SIZE	PACKAGING
<b>SF-SC02UV</b>	Universal Up to shoe size 11	Regular Packaging 200 / dispenser box - 4 dispenser boxes / case
<b>SF-SC02UV-8</b>	Universal Up to shoe size 11	Volume Packaging 100 / bag - 8 bags / case
<b>SF-SC02XL</b>	Extra-Large Size 11 - 15	Regular Packaging 200 / dispenser box - 4 dispenser boxes / case
<b>SF-SC02XL-8</b>	Extra-Large Size 11 - 15	Volume Packaging 100 / bag - 8 bags / case



IMPERIAL SURGICAL LTD.

# STACKABLE STEP STOOLS

Made of lightweight cast aluminium and protected with an antimicrobial powder coating, Imperial Surgical Ltd.'s patented Stackable Step Stools have been designed with infection prevention as a prime focus.



## PRODUCT SELECTION

ITEM #	DESCRIPTION	COLOR
<b>OR-36363</b>	Stackable Step Stool with mat	Black
<b>OR-36363-02</b>	Stackable Step Stool with mat	Dark Blue
<b>OR-36363-03</b>	Stackable Step Stool with mat	Green
<b>OR-36363-07</b>	Stackable Step Stool with mat	Orange
<b>OR-36364</b>	Stackable Step Stool without mat	N/A

PACKAGING: Sold by Unit (Dealers: 4 / Box)

## OPTIONAL ACCESSORIES

ITEM #	DESCRIPTION	PACKAGING
<b>OR-3693</b>	Transport cart	1 Unit
<b>809-1705</b>	Reusable Anti-Fatigue Mat with adhesive strips to secure it in place	5 / Box
<b>899-1710</b>	Single-Use Anti-Fatigue Mat with adhesive strips to secure it in place	10 / Box

## REPLACEMENT PARTS

ITEM #	DESCRIPTION	PACKAGING
<b>OR-36363-1A</b>	Replacement mat for OR-36363	1 Unit
<b>OR-36363-1B</b>	Replacement Foot Pegs	1 Unit



**PATENTED**

- Antimicrobial coating prevents the spread of dangerous micro-organisms and suppresses the growth of mold, bacteria, yeast, and mildew
- Cart washable. Model OR-36364 (without mat) can be sterilized
- Anti-slip surface for enhanced stability
- Rubber-free corners avoids fluid collection
- No welds or large hollow openings to trap debris or fluids
- Transport cart available
- 300 lb (136 kg) weight capacity
- 7 year warranty covers all parts (step stool only)

TIP: COLOR-CODE YOUR STEP STOOLS BY DEPARTMENT



**LEARN MORE at**  
[www.stackingstepstool.com](http://www.stackingstepstool.com)

# BATRIK BLUE CLOUD-GEL™ ANTI-FATIGUE MATS


Our Blue Cloud-Gel™ Anti-Fatigue Mats have been scientifically engineered for use in the O.R. using patented gel technology and premium foam to provide the ultimate comfort and support for surgeons, nurses, and technicians particularly when movement is restricted for prolonged durations.

- Proprietary 3/4" Blue Cloud-Gel™ foam technology will not "bottom-out" like conventional foam mats
- Conforms to and supports feet promoting efficient circulation, reducing physical discomfort and fatigue
- Promotes proper posture while diminishing back, knee, and foot discomfort
- Latex free, water resistant, and nonabsorbent
- Reusable, easy to clean, and Cart Washer compatible (contact us for more details)
- Excellent friction properties when damp or dry surgical booties are worn
- Completely sealed with an RF weld for excellent protection against moisture absorption
- US Standards: CALL 117 and NFPA 260 Compliant
- 1 Year Warranty

DID YOU KNOW? STUDIES HAVE DEMONSTRATED THAT ANTI-FATIGUE MATS IMPROVE PERFORMANCE, ENERGY, AND POST-OPERATIVE DISCOMFORT.

## PRODUCT SELECTION

ITEM #	DIMENSIONS	PACKAGING
<b>800-2401</b>	18" X 24" / 46 cm x 61 cm	1 / Box
<b>800-2405</b>	18" X 24" / 46 cm x 61 cm	5 / Box
<b>800-3201</b>	20" X 32" / 51 cm x 81 cm	1 / Box
<b>800-3203</b>	20" X 32" / 51 cm x 81 cm	3 / Box
<b>800-4801</b>	20" X 48" / 51 cm x 122 cm	1 / Box
<b>800-4802</b>	20" X 48" / 51 cm x 122 cm	2 / Box
<b>800-7201</b>	30" X 72" / 76 cm x 183 cm (Special Order)	1 / Box






# BATRIK BLUE-CLOUD NIMBUS™ ANTI-FATIGUE MATS

Our Blue-Cloud Nimbus™ Anti-Fatigue Mats are designed for general hospital use. Made with proprietary high-density “energy-return” polyurethane materials, they provide unmatched comfort, relief, and ergonomic support.

- Our Energy-Return properties will help reduce leg, foot, lower back discomfort, and work-related injuries
- Longer lasting support as 3/4" proprietary foam blend will not bottom out like most other foam based anti-fatigue mats
- Antimicrobial properties blended into the mat provides added protection
- Highly water resistant surface will not absorb water or dirt
- Lightweight and easy to clean
- Textured design on surface and high-traction base ensures greater safety. (US Standard NFSI 101-C)
- Beveled edging reduces tripping hazards and anti-slip backing enhances safety
- Durable design is ideal for commercial applications, yet still ergonomically and aesthetically pleasing for home or office environments
- Flammability/Burn resistant per CFR 16 Parts 1630/1631
- Made with eco-friendly materials, is non-toxic, and latex free
- 5 Year Warranty (non-industrial usage), 4 Year Warranty (industrial usage)

GOOD TO KNOW: IF WITHIN  
30 DAYS YOU ARE NOT 100%  
SATISFIED, CONTACT YOUR  
DEALER FOR A REFUND  
(EXCLUDES TRANSPORT,  
HANDLING, AND TAXES)



Available in  
Taupe or Black

## PRODUCT SELECTION

ITEM #	DIMENSIONS INCHES / CM	PACKAGING
<b>BLACK COLOR</b>		
<b>809-3201</b>	20" X 32" / 51 cm x 81 cm	1 / Box
<b>809-3204</b>	20" X 32" / 51 cm x 81 cm	4 / Box
<b>809-4801</b>	20" X 48" / 51 cm x 122 cm	1 / Box
<b>809-4804</b>	20" X 48" / 51 cm x 122 cm	4 / Box
<b>809-3601</b>	24" X 36" / 61 cm x 91 cm	1 / Box
<b>809-7201</b>	20" X 72" / 51 cm x 183 cm	1 / Box
<b>809-6001</b>	36" X 60" / 91 cm x 152 cm	1 / Box
<b>TAUPE COLOR</b>		
<b>810-3201</b>	20" X 32" / 51 cm x 81 cm	1 / Box
<b>810-3204</b>	20" X 32" / 51 cm x 81 cm	4 / Box
<b>810-4801</b>	20" X 48" / 51 cm x 122 cm	1 / Box
<b>810-4804</b>	20" X 48" / 51 cm x 122 cm	4 / Box
<b>810-3601</b>	24" X 36" / 61 cm x 91 cm	1 / Box
<b>810-7201</b>	20" X 72" / 51 cm x 183 cm	1 / Box
<b>810-6001</b>	36" X 60" / 91 cm x 152 cm	1 / Box

# BATRIK GOLFF™ ANTI-FOG SOLUTION

Our GOLFF™ Anti-Fog Solution solves the fogging problem often encountered during Laparoscopic surgical procedures.

- Prevents the occurrence of fog on lenses
- Ensures clearer vision and avoids unnecessary interruptions during surgery
- Set includes an easy-squeeze 6 mL bottle and an adhesive radiopaque sponge for secure placement
- OEM (available)
- Sterile
- Latex free



## PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING
GOLFF	Anti-Fog Solution & Adhesive Sponge set	20 / Box

# BATRIK ENDOSCOPY BO2™ BITE BLOCKS


## WITH NASAL OXYGEN DELIVERY

Our BO2™ Bite Blocks significantly enhance patient safety by providing supplemental oxygenation during gastrointestinal endoscopy, preventing hypoxia in sedated patients, and ensuring superior patient comfort.

- Saves time and money by eliminating the need for separate nasal prongs and oxygen masks
- Designed to be free of sharp edges to ensure patient comfort and safety
- Available exclusively in Single-Use
- Optional retention strap is available to ensure a secure fit
- The Bite Blocks and optional retention straps are all latex free, lead free and DEHP free

**BO2**  
(FRONT SIDE)

**BO2**  
(BACK SIDE)



## PRODUCT SELECTION

ITEM #	TYPE	DESCRIPTION	PACKAGING
<b>BO2</b>	Single-Use	18mm / 54FR, includes strap	100 / Box
<b>BO2-N</b>	Single-Use	18mm / 54FR, without strap	100 / Box

BATRIK

# ENDOSCOPY BITESTRONG™ BITE BLOCKS

## ADULT & PEDIATRIC BITE BLOCKS

Our Bitestrong™ Bite Blocks provide a safe passage through the patient's mouth during TEE, protecting the patient's teeth and endoscope.

Our molding process ensures a smooth finish with no sharp edges (flashing) that can often irritate a patient's mouth. All Bite Blocks are individually wrapped in sealed, easy to open PE bags and are latex free, lead free, and DEHP free. Additional straps are available, if required

### SINGLE USE

- Available without dental retention

### REUSABLE

- Designed for repeated usage and will withstand commonly practiced chemical and steam sterilization cycles (134°C / 275°F).
- Available with or without an adjustable retention strap

## PRODUCT SELECTION

ITEM #	TYPE	DESCRIPTION	PACKAGING
<b>B1</b>	Single-Use <a href="#">Pediatric</a>	16 mm / 48FR, includes strap	100 / Box
<b>B1-25</b>	Single-Use <a href="#">Pediatric</a>	16 mm / 48FR, includes strap	25 / Box
<b>B1-N</b>	Single-Use <a href="#">Pediatric</a>	16 mm / 48FR, without strap	100 / Box
<b>B3</b>	Single-Use	20 mm / 60FR, includes strap	100 / Box
<b>B3-N</b>	Single-Use	20 mm / 60FR , without strap	100 / Box
<b>B3-SA</b>	Single-Use	20 mm / 60FR, with strap <u>attached</u> to one side	100 / Box
<b>B3-NDR</b>	Single-Use	20 mm / 60FR, no dental retention rim	100 / Box
<b>B3-NDR-N</b>	Single-Use	Same as B3-N, no dental retention rim, without strap.	100 / Box
<b>B4</b>	Single-Use	19 mm / 60FR	100 / Box
<b>B4T</b>	Single-Use	19 mm / 60FR, includes strap	100 / Box
<b>B4T-N</b>	Single-Use	19 mm / 60FR, without strap	100 / Box
<b>B6</b>	Reusable	19 mm / 60FR	100 / Box
<b>B6T</b>	Reusable	19 mm / 60FR, includes strap	100 / Box
<b>B6T-N</b>	Reusable	19 mm / 60FR, without strap	100 / Box
<b>B6-25</b>	Reusable	19 mm / 60FR	25 / Box
<b>B6T-25</b>	Reusable	19 mm / 60FR, includes strap	25 / Box
<b>B6T-N-25</b>	Reusable	19 mm / 60FR, without strap	25 / Box
<b>BS</b>	Single-Use	Bite-Block straps 21.7" x .79" / 55 cm x 2 cm	100 / Box



**BS**



**B1**



**B3**



**B4**



**B6**

# BATRIK COUVRETTE™ VALVE COVERS

Our Couvrette™ Valve Covers eliminate the need for manual cleaning and reprocessing of reusable valves. Our single use Valve Covers effectively reduce cross contamination between patients and also reduce fluid leakage from the parts during procedures.

- Easily and securely fit to all valves
- Reduce the potential for errors
- Single-Use
- Non-Sterile
- Latex free



16-0044

TIP: REUSABLE VALVE COVERS ARE  
NOTORIOUSLY DIFFICULT TO CLEAN.  
SAVE TIME AND REDUCE RISK  
WITH COUVRETTE™.



16-0022

## PRODUCT SELECTION

ITEM #	DESCRIPTION	COLOR	PACKAGING
<b>16-0022</b>	Biopsy Valve Cover - Olympus® & Fujinon®	Blue	100 / Box
<b>16-0044</b>	Biopsy Valve Cover - Pentax®	Pink	100 / Box

**NOTE:** Contact us to verify compatibility with your scope


Batrik

# DAYBRUSH™ CHANNEL BRUSHES

## REUSABLE

Our Reusable DayBrush™ Channel Brushes can be used for cleaning endoscopes for an entire day. In fact, it's been proven to effectively clean up to 10 scopes before the need to discard.\*

- Soft tip ensures easy manipulation within the channel with no risk of damage
- The back-end is sealed and softened to prevent fluid from getting trapped inside the tube. The softening process also prevents damage if the cleaning brush is used in reverse.
- Can be steam sterilized.
- Bristle volume is 20mm in length – nearly double most brushes on the market. A large brush provides more surface contact during the cleaning process.
- Significant per use savings when compared to other reusable or disposable channel cleaning brushes.



*(\*) Discard after cleaning 10 scopes or sooner if any brush degradation is noticed. After ten applications, the brush may have the appearance that it can continue to clean effectively, but we discourage prolonged use due to risk of cross contamination. Actual number of uses will vary based on handling and processing.*

## PRODUCT SELECTION

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
<b>30-18050-T</b>	70" / 180 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Single
<b>30-23050-T</b>	90" / 230 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Single
<b>30-23050-2T</b>	90" / 230 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Double

PACKAGING: 50 / Box

## OPTIONAL KIT ADD-ONS

ADD THE SUFFIX:

"F" to receive a kit with our Flat Sponge (Packaging: Box of 40). Applicable on certain models only, call for details.

"K" to receive a kit with the single use Valve / Control Cleaning Brush 5/11 Double Brushes (item # 27-57787). Applicable on all models.

"FK" to receive both the sponge and valve cleaning brushes.




# BATRIK ARC™ CHANNEL BRUSHES

## SINGLE-USE

BATRIK offers one of the largest selection of Endoscopy Channel Cleaning Brushes in the industry. Our brushes are available in a variety of lengths, materials, and styles to accommodate the widest range of endoscopic devices used in the healthcare industry today.

No matter what your requirements, you'll be able to find the brush that suits your specific needs. If you don't see the brush you require, simply call us. We have over 50 different brushes available and can custom fabricate to your precise specifications.



### OPTIONAL KIT ADD-ONS

ADD THE SUFFIX:

"F" to receive a kit with our Flat Sponge (Packaging: Box of 40). Applicable on certain models only, call for details.

"K" to receive a kit with the single use Valve / Control Cleaning Brush 5/11 Double Brushes (item # 27-57787). Applicable on all models.

"FK" to receive both the sponge and valve cleaning brushes.

## SERIES 333 & 311



### PRODUCT SELECTION

#### 5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.85 MM - PUNCHY PINK POLYPROPYLENE

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
<b>311-23050</b>	90" / 230 cm	0.94" / 24 mm	0.20" / 5 mm	0.073" / 1.85 mm	Single
<b>311-23050-2</b>	90" / 230 cm	0.94" / 24 mm	0.20" / 5 mm	0.073" / 1.85 mm	Double
<b>311-18050</b>	70" / 180 cm	0.94" / 24 mm	0.20" / 5 mm	0.073" / 1.85 mm	Single
<b>311-10050</b>	39" / 100 cm	0.94" / 24 mm	0.20" / 5 mm	0.073" / 1.85 mm	Single

PACKAGING: 50 / BOX

#### 5/6 MM TAPERED BRUSHES – SHEATH Ø 1.85 MM - PUNCHY PINK POLYPROPYLENE

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
<b>333-23056</b>	90" / 230 cm	0.79" / 20 mm	0.20", 0.24" / 5, 6 mm	0.073" / 1.85 mm	Single
<b>333-23056-2</b>	90" / 230 cm	0.79" / 20 mm	0.20", 0.24" / 5, 6 mm	0.073" / 1.85 mm	Double

PACKAGING: 50 / Box

## SERIES 300 & 322



### PRODUCT SELECTION

#### 5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.75 MM - METALLIC GREY DELRIN

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
<b>300-18050</b>	70" / 180 cm	0.79" / 20 mm	0.20" / 5 mm	0.070" / 1.75 mm	Single
<b>300-23050</b>	90" / 230 cm	0.79" / 20 mm	0.20" / 5 mm	0.070" / 1.75 mm	Single
<b>300-23050-2</b>	90" / 230 cm	0.79" / 20 mm	0.20" / 5 mm	0.070" / 1.75 mm	Double


PACKAGING: 50 / Box

#### 5/6 MM TAPERED BRUSHES – SHEATH Ø 1.75 MM - METALLIC GREY DELRIN

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
<b>322-23056</b>	90" / 230 cm	0.79" / 20 mm	0.20", 0.24" / 5, 6 mm	0.070" / 1.75 mm	Single
<b>322-23056-2</b>	90" / 230 cm	0.79" / 20 mm	0.20", 0.24" / 5, 6 mm	0.070" / 1.75 mm	Double

PACKAGING: 50 / Box

## SERIES 344




### PRODUCT SELECTION

#### 3 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.55 MM - METALLIC PURPLE DELRIN

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
<b>344-23030</b>	90" / 230 cm	0.79" / 20 mm	0.12" / 3 mm	0.060" / 1.55 mm	Single
<b>344-18030</b>	70" / 180 cm	0.79" / 20 mm	0.12" / 3 mm	0.060" / 1.55 mm	Single
<b>344-10030</b>	39" / 100 cm	0.79" / 20 mm	0.12" / 3 mm	0.060" / 1.55 mm	Single

PACKAGING: 50 / Box

## SERIES 35



### PRODUCT SELECTION

#### 2.5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.00 MM - YELLOW DELRIN

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
<b>35-23025</b>	90" / 230 cm	0.59" / 15 mm	0.098" / 2.5 mm	0.039" / 1.00 mm	Single
<b>35-12025</b>	47.2" / 120 cm	0.59" / 15 mm	0.098" / 2.5 mm	0.039" / 1.00 mm	Single

PACKAGING: 50 / Box

#### 2 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.00 MM - YELLOW DELRIN

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
<b>35-23020</b>	90" / 230 cm	0.59" / 15 mm	0.079" / 2.0 mm	0.039" / 1.00 mm	Single
<b>35-12020</b>	47.2" / 120 cm	0.59" / 15 mm	0.079" / 2.0 mm	0.039" / 1.00 mm	Single
<b>35-18020</b>	70.8" / 180 cm	0.59" / 15 mm	0.079" / 2.0 mm	0.039" / 1.00 mm	Single

PACKAGING: 50 / Box

BATRIK

# VALVE / CONTROL BRUSHES


27-05197



27-57787



40-18511



36-10010



27-15311



## PRODUCT SELECTION


ITEM #	DESCRIPTION	LENGTH	BRUSH 1 LENGTH	BRUSH 2 LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)
<b>27-05197</b>	Double Trumpet 5/7 White Color	5.00" / 13 cm	0.20" / 5 mm	1.00" / 25 mm	0.20" / 5 mm	0.28" / 7 mm
PACKAGING: Bag of 3 (single-use recommended). Volume Packaging: Bag of 50 use ITEM # 27-05197-50						
<b>27-57787</b>	Double. Universal 5/11 Yellow Color	6.00" / 15 cm	1.18" / 30 mm	0.79" / 20 mm	0.43" / 11 mm	0.20" / 5 mm
PACKAGING: Bag of 3 (single-use recommended). Volume Packaging: Bag of 50 use ITEM # 27-57787-50						
<b>40-18511</b>	Double Universal 5/11 Extra Long Blue Color	8.07" / 20.5 cm	1.57" / 40 mm	0.79" / 20 mm	0.43" / 11 mm	0.20" / 5 mm
PACKAGING: Bag of 3 (reusable and can be steam sterilized). Volume Packaging: Bag of 50 use ITEM # 40-18511-50						
<b>36-10010</b>	Single end large Red Color	4.3" / 11 cm	1.57" / 40 mm	N/A	0.39" / 10 mm	N/A
PACKAGING: Bag of 3 (single-use recommended). Volume Packaging: Bag of 50 use ITEM # 36-10010-50						
<b>27-15311</b>	Single end large Blue Color	4" / 10.3 cm	1.57" / 40 mm	N/A	0.39" / 10 mm	N/A
PACKAGING: Bag of 3 (single-use recommended). Volume Packaging: Bag of 50 use ITEM # 27-15311-50						

# BATRIK

## V-BRUSHES

### EXTERNAL CLEANING BRUSHES

Our External Cleaning Brushes for endoscopes / optical instruments provide excellent cleaning capabilities. With their unique U-shape bristles and innovative design, our brushes are the ultimate tool for cleaning the exterior of your endoscopic instruments. Free from sharp or metallic ends that can scratch or cause damage to valuable equipment, both single-use and reusable models are latex free and feature ergonomic handles.



## PRODUCT SELECTION



ITEM #	DESCRIPTION	BRUSH LENGTH / (Ø)	HANDLE LENGTH / (Ø)	TOTAL LENGTH	PACKAGING
<b>36-20005</b>	Single-Use External Cleaning Brush	70 mm / 25 mm	80 mm / 3 mm	150 mm	50 / Box
<b>36-20010</b>	Reusable External Cleaning Brush	90 mm / 25 mm	120 mm / 12 mm	210 mm	3 / Bag




# BATRIK SILICONE PIN MAT

Our Silicone Pin Mat cushions, protects, and organizes delicate surgical instruments for sterilization. The elevated nubs provide greater air-flow ventilation and drying, significantly reducing residual moisture collection.

- Reusable
- Autoclavable
- Latex free



**FRONT SIDE**  
Pins keep instruments securely in place



**BACK SIDE**  
Numerous holes and elevated nubs provide superior venting and reduce water collection

## PRODUCT SELECTION


ITEM #	DIMENSIONS	PACKAGING
<b>13-2001</b>	13" x 18" (30 cm x 45 cm)	1 unit
<b>13-2005</b>	13" x 18" (33 cm x 45 cm)	5 units

GOOD TO KNOW:  
CUSTOM SIZES  
ARE AVAILABLE  
UPON REQUEST


## BATRIK CORNER POCKETS™

Our Corner Pockets™ are the perfect tool to prevent tearing associated with sharp corners or feet on surgical trays. They fit all types of corners, including round, and the reusable model can withstand repeated sterilization up to 134°C. After 10 cycles, no product degradation was observed.

**18-CP07**  
(REUSABLE)



**18-CP02**  
(SINGLE USE)



- Latex Free
- No moisture retention
- Easy to use and position
- Soft, flexible plastic feels like silicone


### PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING
<b>18-CP07</b>	Corner Pockets™ Reusable	80 / Bag
<b>18-CORNER</b>	Corner Pockets™ Reusable	320 / Box (4 Bags of 80 pcs)
<b>18-EVAL7</b>	Corner Pockets™ Reusable	24 / Bag
<b>18-CP02*</b>	Corner Pockets™ Single-Use	240 / Bag
<b>18-EVAL2*</b>	Corner Pockets™ Single-Use	24 / Bag
<b>18-CPF-02</b>	Soft Pockets™ Single-Use Foam Corner Protectors	80 / Box (24 Boxes / Pallet)

\* Although the single use Corner Pockets™ appears similar to the reusable model it is blended differently and has only been validated as single-use. Re-using it may compromise sterility.

## SOFT POCKETS™

Our Soft Pockets™ Foam Corner Protectors feature foam that is 0.25" thicker than competitor products. They are especially useful with heavier instrument sets.





## BATRIK

# INSTRUMENT PROTECTORS


Our Instrument Protectors protect the tips and edges of instruments to prevent damage during handling and sterilization. Available in a wide array of color coded sizes for almost every application, including osteotome, retractor, double, and perforated.

- Steam and gas permeable
- Perforated protectors reduce vacuum seal and moisture build-up
- Validated for pre-vacuum sterilization
- Multiple color-coded sizes
- Latex Free
- Non-Sterile

### REGULAR / CLEAR



### PERFORATED / CLEAR



### REGULAR / TINTED



### PERFORATED / TINTED




### REGULAR / SOLID COLOR



### PERFORATED / SOLID COLOR



## PRODUCT SELECTION

	ITEM # REGULAR	ITEM # PERFORATED	COLOR	DIMENSIONS
<b>REGULAR &amp; PERFORATED / CLEAR</b>				
 RC-07  PC-07	RC-01	PC-01	Clear	0.062" x 0.75" (1.57 mm x 19 mm)
	RC-02	PC-02	Clear	0.075" x 0.75" (2 mm x 19 mm)
	RC-03	PC-03	Clear	0.109" x 0.75" (2.77 mm x 19 mm)
	RC-04	PC-04	Clear	0.125" x 1" (3.18 mm x 25.4 mm)
	RC-05	PC-05	Clear	0.187" x 1" (4.75 mm x 25.4 mm)
	RC-06	PC-06	Clear	0.406" x 0.75" (10.3 mm x 19 mm)
	RC-07	PC-07	Clear	0.062" x 0.375" x 1" (1.57 mm x 9.53 mm x 25.4 mm)
	RC-08	PC-08	Clear	0.078" x 0.625" x 1" (2 mm x 15.875 mm x 25.4 mm)
	RC-09	PC-09	Clear	0.125" x 1" x 1" (3.18 mm x 25.4 mm x 25.4 mm)
	RC-10	PC-10	Clear	Assorted
PACKAGING: Bag of 100				
<b>REGULAR &amp; PERFORATED / TINTED</b>				
 RT-07  PT-07	RT-01	PT-01	Tinted White	0.062" x 0.75" (1.57 mm x 19 mm)
	RT-02	PT-02	Tinted Blue	0.075" x 0.75" (2 mm x 19 mm)
	RT-03	PT-03	Tinted Green	0.109" x 0.75" (2.77 mm x 19 mm)
	RT-04	PT-04	Tinted Red	0.125" x 1" (3.18 mm x 25.4 mm)
	RT-05	PT-05	Tinted Yellow	0.187" x 1" (4.75 mm x 25.4 mm)
	RT-06	PT-06	Tinted Maroon	0.406" x 0.75" (10.3 mm x 19 mm)
	RT-07	PT-07	Tinted Orange	0.062" x 0.375" x 1" (1.57 mm x 9.53 mm x 25.4 mm)
	RT-08	PT-08	Tinted Brown	0.078" x 0.625" x 1" (2 mm x 15.875 mm x 25.4 mm)
	RT-09	PT-09	Tinted Black	0.125" x 1" x 1" (3.18 mm x 25.4 mm x 25.4 mm)
	RT-10	PT-10	Tinted Assorted	Assorted
PACKAGING: Bag of 100				
<b>REGULAR &amp; PERFORATED / SOLID</b>				
 RS-07  PS-07	RS-01	PS-01	Solid White	0.062" x 0.75" (1.57 mm x 19 mm)
	RS-02	PS-02	Solid Blue	0.075" x 0.75" (2 mm x 19 mm)
	RS-03	PS-03	Solid Green	0.109" x 0.75" (2.77 mm x 19 mm)
	RS-04	PS-04	Solid Red	0.125" x 1" (3.18 mm x 25.4 mm)
	RS-05	PS-05	Solid Yellow	0.187" x 1" (4.75 mm x 25.4 mm)
	RS-06	PS-06	Solid Maroon	0.406" x 0.75" (10.3 mm x 19 mm)
	RS-07	PS-07	Solid Orange	0.062" x 0.375" x 1" (1.57 mm x 9.53 mm x 25.4 mm)
	RS-08	PS-08	Solid Brown	0.078" x 0.625" x 1" (2 mm x 15.875 mm x 25.4 mm)
	RS-09	PS-09	Solid Black	0.125" x 1" x 1" (3.18 mm x 25.4 mm x 25.4 mm)
	RS-10	PS-10	Solid Assorted	Assorted
PACKAGING: Bag of 100				

ITEM #	DIMENSIONS
<b>OP-11</b>	0.25" x 0.75" (6.35 mm x 19.05 mm)
<b>OP-12</b>	0.375" x 1" (9.53 mm x 25.4 mm)
<b>OP-13</b>	0.5" x 1" (12.7 mm x 25.4 mm)
<b>OP-14</b>	0.75" x 1" (19.05 mm x 25.4 mm)
<b>OP-15</b>	1" x 1" (25.4 mm x 25.4 mm)
<b>OP-16</b>	1.25" x 1.25" (31.75 mm x 31.75 mm)
<b>OP-17</b>	1.5" x 1.5" (38.1 mm x 38.1 mm)
<b>OP-18</b>	1.75" x 1.25" (44.45 mm x 31.75 mm)
<b>OP-20</b>	Assorted

**PACKAGING:** Bag of 50

ITEM #	DIMENSIONS	PACKAGING
<b>RGS-01</b>	1.5" (38 mm)	55 / Bag
<b>RGS-02</b>	2" (51 mm)	45 / Bag
<b>RGS-03</b>	3" (76.2 mm)	40 / Bag


ITEM #	COLOR	SIZE
<b>RGM-01</b>	Pink	Small
<b>RGM-02</b>	Clear	Small
<b>RGM-03</b>	Green	Medium
<b>RGM-04</b>	Clear	Medium

**PACKAGING:** Bag of 50


ITEM #	COLOR	DIMENSIONS
<b>DG-01</b>	Clear	0.0625" x 0.75" (1.59 mm x 19.05 mm)
<b>DG-02</b>	Blue	0.078" x 0.75" (1.98 mm x 19.05 mm)
<b>DG-03</b>	Green	0.109" x 0.75" (2.78 mm x 19.05 mm)
<b>DG-04</b>	Red	0.125" x 1" (3.18 mm x 25.4 mm)
<b>DG-05</b>	Yellow	0.187" x 1" (4.76 mm x 25.4 mm)
<b>DG-07</b>	Orange	0.062" x 0.375" x 1" (1.59 mm x 9.53 mm x 25.4 mm)
<b>DG-10</b>	Assorted	Assorted

**PACKAGING:** Bag of 50

## OSTEOTOME PROTECTORS




## SINGLE-TIP RETRACTOR PROTECTORS



## WEITLANER RETRACTOR PROTECTORS



## PERFORATED DOUBLE PROTECTORS




# BATRIK S.T.I.C.S.

Our Simple Tracking Identification Color System, or S.T.I.C.S., is an easy to use, latex free, tool that allows you to quickly identify your instruments by department and/or procedure. The tape adheres securely and every tape sheet comes with a protective cover to ensure it lasts as long as you need.

Each 8.5" x 11" (216 mm x 279 mm) sheet contains the following adhesive strips:

- 1" (25.4 mm) x 220
- 1.5" (38.1 mm) x 44
- 2" (50.8 mm) x 44
- Total tape per sheet 374" (9500 mm)



## PRODUCT SELECTION

SOLID COLORS	
ITEM #	COLOR
T-850-2-01	White
T-850-2-02	Royal Blue
T-850-2-03	Kelly Green
T-850-2-04	Red
T-850-2-05	Yellow
T-850-2-06	Purple
T-850-2-07	Orange
T-850-2-08	Brown
T-850-2-09	Black
T-850-2-10	Gray
T-850-2-11	Lavender
T-850-2-12	Beige
T-850-2-13	Light Blue
T-850-2-14	Lime Green
T-850-2-15	Pink
T-850-2-17	Plum
T-850-2-18	Mint Green

BLACK STRIPES	
ITEM #	COLOR
T-850-3-01	Black/White
T-850-3-03	Black/Kelly Green
T-850-3-04	Black/Red
T-850-3-05	Black/Yellow
T-850-3-07	Black/Orange
T-850-3-10	Black/Gray
T-850-3-11	Black/Lavender
T-850-3-12	Black/Beige
T-850-3-13	Black/Light Blue
T-850-3-14	Black/Lime Green
T-850-3-15	Black/Pink
T-850-3-17	Black/Plum

TRI-COLOR	
ITEM #	COLOR
T-850-5-51	Red/White/Blue
T-850-5-54	Red/Black/Green

SPECIALTY	
ITEM #	COLOR
T-850-7-04	Yellow/Red


WHITE STRIPES	
ITEM #	COLOR
T-850-4-02	White/Royal Blue
T-850-4-03	White/Kelly Green
T-850-4-04	White/Red
T-850-4-05	White/Yellow
T-850-4-07	White/Orange
T-850-4-08	White/Brown
T-850-4-11	White/Lavender
T-850-4-13	White/Light Blue
T-850-4-14	White/Lime Green
T-850-4-15	White/Pink
T-850-4-17	White/Plum

FOUR-COLOR	
ITEM #	COLOR
T-850-2-31	Red/Yellow/ Purple/Orange
T-850-2-32	Blue/White/ Black/Green

## BATRIK ID TAPE



### AVAILABLE WIDTHS



A simple solution for marking surgical instruments with the aid of color-coding, our ID Tape is durable and will not chip, flake, or peel. It can be used in autoclaves, withstanding repeated exposure to sterilization. Independent validation tests demonstrated the efficiency of sterilization of marked instruments, even underneath the BATRIK ID Tape.



Quickly identifies instruments by department or procedure



Easy to apply



Withstands repeated exposure to sterilization



Has more per tape roll (300" / 7620 mm)



Available in 9 colors and 2 widths



Conveniently packaged in a clear box



Non-Sterile

### PRODUCT SELECTION

ITEM # 0.125" / 318 MM	ITEM# 0.25" / 635 MM	COLOR
T-125-01	T-250-01	White
T-125-02	T-250-02	Blue
T-125-03	T-250-03	Green
T-125-04	T-250-04	Red
T-125-05	T-250-05	Yellow
T-125-06	T-250-06	Purple
T-125-07	T-250-07	Orange
T-125-08	T-250-08	Brown
T-125-09	T-250-09	Black



**BATRIK**

# TAPE-B-GONE™ TAPE REMOVER

**Tape-B-Gone™** is a uniquely formulated and ready-to-use proprietary formula that easily removes adhesives from medical equipment and instruments. Unlike many other brands, Tape-B-Gone™ does not contain highly dangerous or caustic chemicals, is not flammable or combustible, and is safe for use on metal, glass, plastic, and rubber without any special equipment or ventilation.



## PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING
<b>T-08RMV</b>	8 oz (237 ml) bottle	1 Unit
<b>T-32RMV</b>	8 oz (237 ml) bottle	4 / Pack
<b>T-SAMP-RMV</b>	0.5 oz (15 ml) bottle	1 Unit

BATRIK


# PADLOCKS SECURITY TAGS

Our Padlock Security Tags for instrument containers and general hospital applications are specifically designed for the healthcare market. Our material selection, fabrication process, and engineering ensure consistent quality and high temperature compliance with the various steam sterilization temperatures in worldwide hospital applications.


Available with optional reactive ink indicator dots that provide visual sterility indication. Light blue dots turn brown/dark brown when steam sterilized.

## AVAILABLE IN THREE FORMATS


### 700-SERIES



### 800-SERIES



### 900-SERIES



## PRINTED & CUSTOM NUMBERED PADLOCKS

Available in 700, 800 & 900 SERIES

### PRODUCT SELECTION

COLOR	PRINTED WORD
<b>PRINTED PADLOCKS</b>	
WHITE	CLEAN
BLUE	STERILE
RED	REPAIR
YELLOW	SHARPEN
YELLOW	LOANER
BLACK	DISCARD
<b>CUSTOM NUMBERED</b>	
ALL COLORS	up to 8 numbers



**CUSTOM NUMBERED**  
Custom printed padlocks (up to 8 characters)



**PRINTED**  
6 different printed labels available

NOTE: Printed & Custom Numbered Padlocks are available in 700, 800, and 900 Series and can be ordered in 100, 200, and 500 unit quantities. Please contact us for ordering details.



Made of polypropylene



Latex and DEHP Free



Will not melt or soften with any established steam sterilization protocol



Fully open design (45°) facilitates insertion into instrument containers



Easy twist breaking design



Large base for easier handling



Available in different size



Custom Colors and OEM available upon request

INDICATOR  
DOTS

## REGULAR



## 700-SERIES PADLOCKS

The 700 series padlocks are perfect for "Aesculap" containers and those with a short, large opening.

## PRODUCT SELECTION

## REGULAR PADLOCKS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-701	15-711	15-710
BLUE	15-702	15-712	15-720
GREEN	15-703	15-713	15-730
RED	15-704	15-714	15-740
YELLOW	15-705	15-715	15-750
PINK	15-706	15-716	15-760
ORANGE	15-707	15-717	15-770
BLACK	15-709	15-719	15-790
LIGHT BLUE	15-722	15-732	15-742

## WITH INDICATORS DOTS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-701-IND	15-711-IND	15-710-IND
BLUE	15-702-IND	15-712-IND	15-720-IND
GREEN	15-703-IND	15-713-IND	15-730-IND
RED	15-704-IND	15-714-IND	15-740-IND
YELLOW	15-705-IND	15-715-IND	15-750-IND
PINK	15-706-IND	15-716-IND	15-760-IND
ORANGE	15-707-IND	15-717-IND	15-770-IND
BLACK	15-709-IND	15-719-IND	15-790-IND
LIGHT BLUE	15-722-IND	15-732-IND	15-742-IND

Double locking hooks can be inserted in a large guide hole




The **grooves** on the padlock facilitate handling



## INDICATOR DOTS

## REGULAR



# 800-SERIES PADLOCKS

The 800 series padlocks have a longer arm (32 mm) and have a single, thinner hook than the 700 series. They are perfect for almost any type of container, even those with long and thin channels.

## PRODUCT SELECTION

### REGULAR PADLOCKS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-801	15-811	15-810
BLUE	15-802	15-812	15-820
GREEN	15-803	15-813	15-830
RED	15-804	15-814	15-840
YELLOW	15-805	15-815	15-850
PINK	15-806	15-816	15-860
ORANGE	15-807	15-817	15-870
BLACK	15-809	15-819	15-890
LIGHT BLUE	15-822	15-832	15-842

### WITH INDICATORS DOTS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-801-IND	15-811-IND	15-810-IND
BLUE	15-802-IND	15-812-IND	15-820-IND
GREEN	15-803-IND	15-813-IND	15-830-IND
RED	15-804-IND	15-814-IND	15-840-IND
YELLOW	15-805-IND	15-815-IND	15-850-IND
PINK	15-806-IND	15-816-IND	15-860-IND
ORANGE	15-807-IND	15-817-IND	15-870-IND
BLACK	15-809-IND	15-819-IND	15-890-IND
LIGHT BLUE	15-822-IND	15-832-IND	15-842-IND

**Single locking hook**  
can be inserted  
in a large guide hole



The locking system  
allows for **easy closure**



**INDICATOR  
DOTS****REGULAR**

## 900-SERIES PADLOCKS

The 900 series padlocks are identical to the 800 series, but with a larger rectangular base which allows for the placement of a 20 mm x 30 mm tracking label.

### PRODUCT SELECTION

#### REGULAR PADLOCKS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-901	15-911	15-910
BLUE	15-902	15-912	15-920
GREEN	15-903	15-913	15-930
RED	15-904	15-914	15-940
YELLOW	15-905	15-915	15-950
PINK	15-906	15-916	15-960
ORANGE	15-907	15-917	15-970
BLACK	15-909	15-919	15-990
LIGHT BLUE	15-922	15-932	15-942

#### WITH INDICATORS DOTS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-901-IND	15-911-IND	15-910-IND
BLUE	15-902-IND	15-912-IND	15-920-IND
GREEN	15-903-IND	15-913-IND	15-930-IND
RED	15-904-IND	15-914-IND	15-940-IND
YELLOW	15-905-IND	15-915-IND	15-950-IND
PINK	15-906-IND	15-916-IND	15-960-IND
ORANGE	15-907-IND	15-917-IND	15-970-IND
BLACK	15-909-IND	15-919-IND	15-990-IND
LIGHT BLUE	15-922-IND	15-932-IND	15-942-IND

**Single locking hook**  
for use with smaller  
diameters



Ideal for adding labels/  
indicators indicators or  
markings to the large  
rectangular base







# BATRIK SAVEFACE™ FACE SHIELD

Our SaveFace™ Face Shield (Visor) is lightweight and allows for good peripheral vision while acting as a physical barrier against splashes. Can be used for extended periods due to comfort-focused design and superior quality materials.

Choose amongst two models:

- The classic **non-draped** model closes the gap between the shield and the chest bone when head is tilted forward for added protection.
- **Draped** offers full-length coverage extending from the top of the eyebrow to below the chest, covering the anterior neck area as per CDC recommendations.



ASSEMBLED IN CANADA

## PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	TOTAL LENGTH	PACKAGING
SF-100	SAVEFACE™ no drape	Shield Length: 9" / 23 cm	9" (23 cm)	100 / Box
SF-180	SAVEFACE™ with ANSI/AAMI PB70 Level 3 compliant protection drape	Shield Length: 9" / 23 cm Drape Length: 6" / 15 cm	15" (38 cm)	100 / Box

# BATRIK SHIELD™ FACE SHIELD

A cost-effective alternative to our premium SAVEFACE™ models, the BATRIK Shield™ provides excellent splash protection. An ideal safety addition to all pandemic kits.

- Self-shaping foam head piece molds comfortably to forehead
- Shields can be used for longer periods due to comfort and improved materials
- Optical shield with anti-fog and anti-glare provide a clearer viewing area
- Soft plastic fasteners adhere the elastic band to the shield for additional comfort and easy use
- Latex free



## PRODUCT SELECTION


ITEM #	DESCRIPTION	DIMENSIONS	TOTAL LENGTH	PACKAGING
<b>SF-096</b>	BATRIK Shield™ no drape	Shield Length: 9" / 23 cm*	9" (23 cm)*	96 / Box
<b>SF-896</b>	BATRIK Shield™ with ANSI/AAMI PB70 Level 3 compliant protection drape	Shield Length: 9" / 23 cm* Drape Length: 7" / 17.8 cm	15.9" / 40.3 cm*	96 / Box

\*Shield product tolerances: 8.9" - 9.5" / 22.5 cm - 23.2 cm

# BATRIK SCRUB PAD

Our Scrub Pad is ideal for pre-cleaning your surgical instruments without leaving any metal residue that often occurs when using steel wool or metal sponges. The abrasiveness is perfectly balanced to provide you with efficient cleaning without damaging your instruments. Its fine fibers allow the user to polish, clean, and debur all types of material surfaces.

- Abrasive Material: Aluminum Oxide
- Grit: 400
- Grade: Very Fine
- Style: Flexible
- Pad material: Non-Woven
- Color: Maroon
- Type: Hand Pad



## PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
<b>60-46050</b>	Scrub Pad	4" x 6" / 10 cm x 15 cm	150 / Box (50 bags of 3 units each)
<b>60-46003</b>	Scrub Pad	4" x 6" / 10 cm x 15 cm	3 / Bag

BATRIK

# PrAmLi Power<sup>3</sup> Enzymes™

## ENZYMATIC DETERGENTS

All **PrAmLi Power<sup>3</sup> Enzymes™** detergents are multi-enzyme pre-soak cleaners for medical equipment and instrumentation. Our proprietary unique blend of enzymes is effective in the removal of all proteins, blood, fatty tissue, and more. It will also remove stains from surgical scrubs, linens, and medical clothing.

Formulated with **PrAmLi Power<sup>3</sup> Enzymes™**, a unique blend of three powerful enzymes along with a surfactant to provide the ultimate in cleaning efficacy



### PRODUCT SELECTION

ITEM #	DILUTION DETERGENT / WATER	PACKAGING
<b>LOW FOAM DETERGENT</b>		
<b>SS-21-1001</b>	4 to 5 ml / L (0.5 to 1oz / Gallon)	1 x 3.8 L (1 Gallon)
<b>SS-21-1004</b>	4 to 5 ml / L (0.5 to 1oz / Gallon)	4 x 3.8 L (4 Gallons)
<b>SS-21-1025</b>	4 to 5 ml / L (0.5 to 1oz / Gallon)	2 x 9.5 L (2 x 2.5 Gallons)
<b>MILD FOAM DETERGENT</b>		
<b>SS-21-2001</b>	4 to 5 ml / L (0.5 to 1oz / Gallon)	1 x 3.8 L (1 Gallon)
<b>SS-21-2004</b>	4 to 5 ml / L (0.5 to 1oz / Gallon)	4 x 3.8 L (4 Gallons)
<b>SUPER CONCENTRATE</b>		
<b>SS-21-3001</b>	1 to 2 ml / L (0.125 to 0.25 oz / Gallon)	1 x 3.8 L (1 Gallon)
<b>SS-21-3004</b>	1 to 2 ml / L (0.125 to 0.25 oz / Gallon)	4 x 3.8 L (4 Gallons)
<b>SS-21-3025</b>	1 to 2 ml / L (0.125 to 0.25 oz / Gallon)	2 x 9.5 L (2 x 2.5 Gallons)

BATRIK

# PrAmLi Power<sup>3</sup> LUBE™ INSTRUMENT LUBRICANT

Specially formulated to lubricate surgical instruments, hinges, box locks, and general equipment during autoclaving. PrAmLi Power<sup>3</sup> LUBE™ is uniquely formulated with an all natural corrosion inhibitor.

- Forms a protective coating which can extend the instrument's, useful life
- Bacteriostatic
- Neutral pH 6.5 - 7.5
- Prevents corrosion



## PRODUCT SELECTION

ITEM #	DILUTION LUBRICANT / WATER	PACKAGING
<b>SS-21-6001</b>	1 to 4 ml / L (0.125 to 0.5 oz / Gallon)	1 x 3.8 L (1 Gallon)
<b>SS-21-6004</b>	1 to 4 ml / L (0.125 to 0.5 oz / Gallon)	4 x 3.8 L (4 Gallons)
<b>SS-21-6025</b>	1 to 4 ml / L (0.125 to 0.5 oz / Gallon)	2 x 9.5 L (2 x 2.5 Gallons)
<b>SS-21-6032</b>	No dilution required	1 x 946 ml (32 oz) Spray bottle
<b>SS-21-6012</b>	No dilution required	12 X 946 ml (32 oz) Spray Bottles



Your First ™

# ENVIRO-SACHET

## ENZYMATIC GEL

Our Enviro-Sachet Enzymatic Gel is a unique triple enzyme cleaning solution that breaks down all bioburden for effective surface and instrument pre-cleaning. It is a simple yet **effective, economical and eco-friendly system** consisting of a 4ml sachet of triple enzyme gel, sponge, and a reusable stainless steel bowl.



Formulated with **PrAmLi Power<sup>3</sup> Enzymes™**



## PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING
17-166	1 Dry Tubular Sponge & Enviro-Sachet kit	25 / Box 4 Boxes / Case
17-266	1 Dry Flat Sponge & Enviro-Sachet kit	50 / Box 4 Boxes / Case
17-177	1 Dry Tubular Sponge, Enviro-Sachet & Stainless Steel Enzymatic Bowl kit	1 Kit
17-277	1 Dry Flat Sponge, Enviro-Sachet & Stainless Steel Enzymatic Bowl kit	1 Kit
17-EZSAMP	SAMPLE KIT INCLUDES: 3 Dry Flat Sponges with Enviro-Sachet kits (17-266) 1 Dry Tubular Sponge with Enviro-Sachet kit (17-166) 1 Reusable Stainless Steel Enzymatic Bowl (H-EZ500)	1 Kit
H-EZ500	Reusable Stainless Steel Enzymatic Bowl 11 cm x 13.5 cm x 7.5 cm / 4.33" x 5.31" x 2.95"	1 Unit

NOTE: Wave Sponge Kits coming soon

Your First ™

# FOAMING SPRAY

## READY-TO-USE ENZYME FOAMING SPRAY

Our Foaming Spray is ideal for instrument pre-cleaning as the thick foam spray enzymatically breaks down all bioburden to allow for effective cleaning and sterilization.

Formulated with **PrAmLi Power<sup>3</sup> Enzymes™**, a unique blend of three powerful enzymes along with a surfactant to provide the ultimate in cleaning efficacy



## PRODUCT SELECTION

ITEM #	PACKAGING
<b>SS-21-4001</b>	1 x 946 ml (32 oz) Spray Bottle
<b>SS-21-4012</b>	12 x 946 ml (32 oz) Spray Bottles

Your First ™

# F1RST WIPES

Our cleanroom wipes are designed to do the first wipe-down of an endoscope after a procedure. They are made from hydroentangled non-woven material composed of 45% polyester and 55% cellulose. These wipes are manufactured and packaged in a cleanroom environment.



- Can be used for a cleanroom environment of ISO Class 5 through 8
- Production process guarantees extremely low particle generation and chemical extractables
- Compatible with virtually any solvent or solution for cleaning and wiping (such as **PrAmLi Power<sup>3</sup> Enzymes™**)
- Double-bag packaging
- Available in two versatile sizes




## PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
<b>17-A0909</b>	Non-woven Cleanroom Wipes	9" x 9" (23 cm x 23 cm)	300 / Bag 14 Bags / Case
<b>17-A1212</b>	Non-woven Cleanroom Wipes	12" x 12" (31 cm x 31 cm)	150 / Bag 10 Bags / Case

Your First ™

# NON-ENZYMATIC BEDSIDE SPONGES

Our sponges are specifically fabricated for use with all medical instrumentation. Unlike most competitor sponges, our sponges are white to better evidence the effectiveness of the cleaning process. Dark colored sponges easily hide contaminants, but with our white sponges you can visually evaluate the pre-clean efficacy.



## PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING
17-100	DRY TUBULAR SPONGE <a href="#">White Dispenser Box</a>	25 / Box (4 Boxes / Case)
17-101	DRY TUBULAR SPONGE <a href="#">Bulk</a>	560 / Case
17-200	DRY FLAT SPONGE <a href="#">White Dispenser Box</a>	50 / Box (4 Boxes / Case)
17-201	DRY FLAT SPONGE <a href="#">Bulk</a>	1,000 / Case

Your First ™

# ENZYMATIC BEDSIDE SPONGE KITS



Our Enzymatic Bedside Sponge Kits feature sponges all pre-treated with **PrAmLi Power<sup>3</sup> Enzymes™** to provide the most effective pre-clean or surface clean solution. This unique three-enzyme blend provides you with confident and effective cleaning. The tubular sponge is ideal for all endoscopes.



## PRODUCT SELECTION


ITEM #	DESCRIPTION	PACKAGING
<b>17-111</b>	DRY TUBULAR SPONGE Moisten in any 250-500 ml basin to activate the enzymatic concentrate.	25 / Box (4 Boxes / Case)
<b>17-211</b>	DRY FLAT SPONGE Moisten in any 250-500 ml basin to activate the enzymatic concentrate.	50 / Box (4 Boxes / Case)
<b>17-122</b>	DRY TUBULAR SPONGE Comes in a freestanding pouch to use as a basin. Simply add 500 ml of water.	15 / Box (4 Boxes / Case)
<b>17-222</b>	DRY FLAT SPONGE Comes in a freestanding pouch to use as a basin. Simply add 500 ml of water.	25 / Box (4 Boxes / Case)
<b>17-133</b>	READY TO USE TUBULAR SPONGE Comes bathing in 500 ml of enzymatic solution.	24 / Case
<b>17-233</b>	READY TO USE FLAT SPONGE Comes bathing in 500 ml of enzymatic solution.	24 / Case
<b>17-144</b>	READY TO USE DUAL PURPOSE TUBULAR SPONGE Comes in a freestanding pouch soaked in 100 ml of enzymatic solution.	12 / Box (4 Boxes / Case)
<b>17-244</b>	READY TO USE DUAL PURPOSE FLAT SPONGE Comes in a freestanding pouch soaked in 100 ml of enzymatic solution.	12 / Box (4 Boxes / Case)
<b>17-155</b>	READY TO USE TUBULAR SPONGE Comes bathing in 220 ml of enzymatic solution.	24 / Case
<b>17-255</b>	READY TO USE FLAT SPONGE Comes bathing in 220 ml of enzymatic solution.	24 / Case

COMING SOON: Wave style sponges for bedside enzymatic kits.



Your First ™

# NON-ENZYMATIC DETERGENT BEDSIDE SPONGE KITS

Our Bedside Sponge Kits with Non-Enzymatic Detergent are **Your First ™** in the medical instrument cleaning process.



## PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING
<b>717-122</b>	DRY TUBULAR SPONGE Comes in a freestanding pouch to use as a basin. Simply add 500 ml of water.	15 / Box (4 Boxes / Case)
<b>717-222</b>	DRY FLAT SPONGE Comes in a freestanding pouch to use as a basin. Simply add 500 ml of water.	25 / Box (4 Boxes / Case)
<b>717-133</b>	READY TO USE TUBULAR SPONGE Comes bathing in 500 ml of detergent solution.	24 / Case
<b>717-233</b>	READY TO USE FLAT SPONGE Comes bathing in 500 ml of detergent solution.	24 / Case
<b>717-144</b>	READY TO USE DUAL PURPOSE TUBULAR SPONGE Comes in a freestanding pouch soaked in 100 ml of detergent solution.	12 / Box (4 Boxes / Case)
<b>717-244</b>	READY TO USE DUAL PURPOSE FLAT SPONGE Comes in a freestanding pouch soaked in 100 ml of detergent solution.	12 / Box (4 Boxes / Case)
<b>717-155</b>	READY TO USE TUBULAR SPONGE Comes bathing in 220 ml of detergent solution.	24 / Case
<b>717-255</b>	READY TO USE FLAT SPONGE Comes bathing in 220 ml of detergent solution.	24 / Case


COMING SOON: Wave style sponges for bedside non-enzymatic detergent kits

Your First ™

# TRANSPORT BAGS FOR ENDOSCOPE AND MEDICAL INSTRUMENTS

Thoughtfully designed to eliminate touch-points during endoscope transport from point-of-care to the reprocessing area, and back, focused on safety of patients and medical staff.

- Made of recyclable LDPE material known to be quite flexible, highly resilient, and tough
- Latex and DEHP-free
- Excellent resistance to dilute and concentrated acids, alcohols, bases and esters\*
- Good resistance to aldehydes, ketones, and vegetable-based oils\*
- Can withstand high temperatures of up to 90 °C (194°F)
- Tamper-evident seal for added safety, also aids in keeping the instrument moist to avoid drying of bio-debris
- Clearly marked bags with color-coded visuals facilitates immediate identification of clean vs. used instruments
- Users can write important information directly on the bag in spaces marked with visual icons
- Available in two convenient sizes

GOOD TO KNOW: USING **Your First ™** SEALABLE TRANSPORT BAGS REDUCES USER ERROR, PREVENTS CROSS-CONTAMINATION, AND PROMOTES PATIENT AND MEDICAL STAFF SAFETY.

## CLEAN INSTRUMENTS BAG



## SOILED INSTRUMENTS BAG

MADE IN CANADA 

## PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
<b>17-HAP24</b>	Clean Instruments Bag	24" x 22" / 61.0 cm x 56.0 cm	200 / Case
<b>17-SAD24</b>	Soiled Instruments Bag	24" x 22" / 61.0 cm x 56.0 cm	200 / Case
<b>17-HAP30</b>	Clean Instruments Bag	30.5" x 26" / 77.0 cm x 66.0 cm	200 / Case
<b>17-SAD30</b>	Soiled Instruments Bag	30.5" x 26" / 77.0 cm x 66.0 cm	200 / Case

\*Reference on file

BATRIK

# ARC™ WIRE STYLE BRUSHES

## SINGLE-USE AND REUSABLE

You've got options with our ARC™ Wire Brushes. Our single-use and reusable brushes are available in both flexible and rigid varieties and come in a multitude of sizes. ARC™ Wire Brushes are guaranteed to last longer than twisted wire and preserve the integrity of your valuable instruments.

### BRUSH SELECTION TIPS:

#### INSTRUMENT SPECIFICATIONS

**What is the maximum length of the instrument(s) to be cleaned?**

Add at least 2" (5cm) to that length, so you have an adequate length to work with.

**What is the inner diameter of the instrument?**

Add at most 0.125" (3mm) to determine the diameter of the bristle to ensure effective cleaning. The bristle diameter has to slightly exceed the inner tube diameter. Bristles that

are too large will bend too much and not clean effectively.

#### BRUSH SELECTION

Simply find the section in our catalogue that best corresponds to the overall length and diameter of your instrument. As a rule, if your present brush is 8" (20.32cm) in length and we only offer 10" (25.4cm) length with your required diameter, that is entirely acceptable.

In many cases, our ARC™ Wire Brushes will have more bristles on the brush. As an example, many 5mm (0.19") diameter brushes are only 1" (2.5cm) in length. At BATRIK, we strived to provide the greatest amount of bristles to ensure the most effective cleaning.

#### CUSTOMIZABLE WIRE BRUSHES


A unique aspect of single-use ARC™ Wire Brushes is that they can be customized to any length:

- Up to 2000mm for Wire Brushes on plastic semi-rigid sheaths, rigid brushes from 5 to 20mm in diameter.
- No limit in length for Wire Brushes on flexible plastic sheaths, flexible brushes from 2-5mm in diameter

## RIGID



## FLEXIBLE



**12" / 30 CM SINGLE-END WIRE BRUSH**

	ITEM # SINGLE-USE	ITEM # REUSABLE	COLOR	TOTAL LENGTH	BRUSH (Ø)	BRUSH LENGTH	SHEATH (Ø)
<b>FLEXIBLE</b> BRUSH & SHEATH	<b>41-31502-S</b>	<b>41-31502</b>	White	12.401" / 31.5 cm	0.078" / 2 mm	0.59" / 15 mm	0.39" / 1 mm
	<b>41-32503-S</b>	<b>41-32503</b>	Yellow	12.795" / 32.5 cm	0.118" / 3 mm	0.98" / 25 mm	0.39" / 1 mm
	<b>411-32504-S</b>	<b>411-32504</b>	Violet	12.795" / 32.5 cm	0.157" / 4 mm	0.98" / 25 mm	0.61" / 1.55 mm
	<b>411-32505-S</b>	<b>411-32505</b>	Grey	12.795" / 32.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75 mm
<b>RIGID</b> BRUSH & SHEATH	<b>41-35005-S</b>	<b>41-35005</b>	White	13.779" / 35 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
	<b>41-35006-S</b>	<b>41-35006</b>	Green	13.779" / 35 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
	<b>41-35008-S</b>	<b>41-35008</b>	Blue	13.779" / 35 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-35010-S</b>	<b>41-35010</b>	Red	13.779" / 35 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-35012-S</b>	<b>41-35012</b>	Yellow	13.779" / 35 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-35015-S</b>	<b>41-35015</b>	Green	13.779" / 35 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
	<b>41-35020-S</b>	<b>41-35020</b>	Blue	13.779" / 35 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING Single-Use: Bag of 20 / Reusable: Bag of 3

**18" / 45 CM SINGLE-END WIRE BRUSH**

	ITEM # SINGLE-USE	ITEM # REUSABLE	COLOR	TOTAL LENGTH	BRUSH (Ø)	BRUSH LENGTH	SHEATH (Ø)
<b>FLEXIBLE</b> BRUSH & SHEATH	<b>41-46502-S</b>	<b>41-46502</b>	White	18.307" / 46.5 cm	0.078" / 2 mm	0.59" / 15 mm	0.39" / 1 mm
	<b>41-47503-S</b>	<b>41-47503</b>	Yellow	18.7" / 47.5 cm	0.118" / 3 mm	0.98" / 25 mm	0.39" / 1 mm
	<b>411-47504-S</b>	<b>411-47504</b>	Violet	18.7" / 47.5 cm	0.157" / 4 mm	0.98" / 25 mm	0.61" / 1.55 mm
	<b>411-47505-S</b>	<b>411-47505</b>	Grey	18.7" / 47.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75 mm
<b>RIGID</b> BRUSH & SHEATH	<b>41-50005-S</b>	<b>41-50005</b>	White	19.685" / 50 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
	<b>41-50006-S</b>	<b>41-50006</b>	Green	19.685" / 50 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
	<b>41-50008-S</b>	<b>41-50008</b>	Blue	19.685" / 50 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-50010-S</b>	<b>41-50010</b>	Red	19.685" / 50 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-50012-S</b>	<b>41-50012</b>	Yellow	19.685" / 50 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-50015-S</b>	<b>41-50015</b>	Green	19.685" / 50 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
	<b>41-50020-S</b>	<b>41-50020</b>	Blue	19.685" / 50 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Single -Use: Bag of 20 / Reusable: Bag of 3

## 24" / 60 CM SINGLE-END WIRE BRUSH

	ITEM # SINGLE-USE	ITEM # REUSABLE	COLOR	TOTAL LENGTH	BRUSH (Ø)	BRUSH LENGTH	SHEATH (Ø)
<b>FLEXIBLE</b> BRUSH & SHEATH	<b>41-61502-S</b>	<b>41-61502</b>	White	24.212" / 61.5 cm	0.078" / 2 mm	0.59" / 15 mm	0.39" / 1mm
	<b>41-62503-S</b>	<b>41-62503</b>	Yellow	24.6" / 62.5 cm	0.118" / 3 mm	0.98" / 25 mm	0.39" / 1mm
	<b>411-62504-S</b>	<b>411-62504</b>	Violet	24.6" / 62.5 cm	0.157" / 4 mm	0.98" / 25 mm	0.61" / 1.55mm
	<b>411-62505-S</b>	<b>411-62505</b>	Grey	24.6" / 62.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75mm
<b>RIGID</b> BRUSH & SHEATH	<b>41-65005-S</b>	<b>41-65005</b>	White	25.59" / 65 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3mm
	<b>41-65006-S</b>	<b>41-65006</b>	Green	25.59" / 65 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3mm
	<b>41-65008-S</b>	<b>41-65008</b>	Blue	25.59" / 65 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-65010-S</b>	<b>41-65010</b>	Red	25.59" / 65 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-65012-S</b>	<b>41-65012</b>	Yellow	25.59" / 65 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-65015-S</b>	<b>41-65015</b>	Green	25.59" / 65 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
	<b>41-65020-S</b>	<b>41-65020</b>	Blue	25.59" / 65 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Single-Use: Bag of 20 / Reusable: Bag of 3

## 30" / 75 CM SINGLE-END WIRE BRUSH

	ITEM # SINGLE-USE	ITEM # REUSABLE	COLOR	TOTAL LENGTH	BRUSH (Ø)	BRUSH LENGTH	SHEATH (Ø)
<b>FLEXIBLE</b> BRUSH & SHEATH	<b>41-76502-S</b>	<b>41-76502</b>	White	30.118" / 76.5 cm	0.078" / 2 mm	0.59" / 15 mm	0.39" / 1 mm
	<b>41-77503-S</b>	<b>41-77503</b>	Yellow	30.511" / 77.5 cm	0.118" / 3 mm	0.98" / 25 mm	0.39" / 1 mm
	<b>411-77504-S</b>	<b>411-77504</b>	Violet	30.511" / 77.5 cm	0.157" / 4 mm	0.98" / 25 mm	0.61" / 1.55 mm
	<b>411-77505-S</b>	<b>411-77505</b>	Grey	30.511" / 77.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75 mm
<b>RIGID</b> BRUSH & SHEATH	<b>41-80005-S</b>	<b>41-80005</b>	White	31.496" / 80 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
	<b>41-80006-S</b>	<b>41-80006</b>	Green	31.496" / 80 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
	<b>41-80008-S</b>	<b>41-80008</b>	Blue	31.496" / 80 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-80010-S</b>	<b>41-80010</b>	Red	31.496" / 80 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-80012-S</b>	<b>41-80012</b>	Yellow	31.496" / 80 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-80015-S</b>	<b>41-80015</b>	Green	31.496" / 80 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
	<b>41-80020-S</b>	<b>41-80020</b>	Blue	31.496" / 80 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Single-Use: Bag of 20 / Reusable: Bag of 3



FLEXIBLE 2MM




FLEXIBLE 3MM



FLEXIBLE 4MM



FLEXIBLE 5MM



## FLEXIBLE

Flexible brush diameter (Ø) ranging from 2 mm to 5 mm

RIGID 20MM



RIGID 15MM



RIGID 12MM



RIGID 10MM



RIGID 8MM



RIGID 6MM



RIGID 5MM



## RIGID

Rigid brush diameter (Ø) ranging from 5 mm to 20 mm

### BRUSH SELECTION SYSTEM (Actual Size)

2mm  
0.078"



3 mm  
0.118"



4 mm  
0.157"



5 mm  
0.196"



6 mm  
0.236"



8 mm  
0.314"



10 mm  
0.394"



12 mm  
0.472"



15 mm  
0.591"



20 mm  
0.787"



# BATRIK REUSABLE DOUBLE ARC™ WIRE BRUSHES ASYMMETRIC

## 14" / 35 CM ASYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-35056	White	13.779" / 35 cm	0.19" / 5 mm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-35810	Blue	13.779" / 35 cm	0.31" / 8 mm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-35112	Yellow	13.779" / 35 cm	0.39" / 10 mm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-35152	Green	13.779" / 35 cm	0.59" / 15 mm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

## 20" / 50 CM ASYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-50056	White	19.685" / 50 cm	0.19" / 5 mm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-50810	Blue	19.685" / 50 cm	0.31" / 8 mm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-50112	Yellow	19.685" / 50 cm	0.39" / 10 mm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-50152	Green	19.685" / 50 cm	0.59" / 15 mm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

## 26" / 65 CM ASYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-65056	White	25.59" / 65 cm	0.19" / 5 mm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-65810	Blue	25.59" / 65 cm	0.31" / 8 mm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-65112	Yellow	25.59" / 65 cm	0.39" / 10 mm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-65152	Green	25.59" / 65 cm	0.59" / 15 mm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

## 31" / 80 CM ASYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-80056	White	31.496" / 80 cm	0.19" / 5 mm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-80810	Blue	31.496" / 80 cm	0.31" / 8 mm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-80112	Yellow	31.496" / 80 cm	0.39" / 10 mm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-80152	Green	31.496" / 80 cm	0.59" / 15 mm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

# BATRIK REUSABLE DOUBLE ARC™ WIRE BRUSHES SYMMETRIC

## 14" / 35 CM SYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-35005	White	13.779" / 35 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
40-35006	Green	13.779" / 35 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-35008	Blue	13.779" / 35 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
40-35010	Red	13.779" / 35 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-35012	Yellow	13.779" / 35 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-35015	Green	13.779" / 35 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-35020	Blue	13.779" / 35 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

## 20" / 50 CM SYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-50005	White	19.685" / 50 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
40-50006	Green	19.685" / 50 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-50008	Blue	19.685" / 50 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
40-50010	Red	19.685" / 50 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-50012	Yellow	19.685" / 50 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-50015	Green	19.685" / 50 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-50020	Blue	19.685" / 50 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

## 26" / 65 CM SYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-65005	White	25.59" / 65 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
40-65006	Green	25.59" / 65 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-65008	Blue	25.59" / 65 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
40-65010	Red	25.59" / 65 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-65012	Yellow	25.59" / 65 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-65015	Green	25.59" / 65 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-65020	Blue	25.59" / 65 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

## 31" / 80 CM SYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-80005	White	31.496" / 80 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
40-80006	Green	31.496" / 80 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-80008	Blue	31.496" / 80 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
40-80010	Red	31.496" / 80 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-65112	Yellow	31.496" / 80 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-80015	Green	31.496" / 80 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-80020	Blue	31.496" / 80 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

## ASYMMETRIC BRUSHES

5 MM

6 MM



8 MM

10 MM



10 MM

12.5 MM



15 MM

20 MM




BATRIK

# REUSABLE IN-BRUSH™ ANTIMICROBIAL TWISTED WIRE BRUSHES

Our In-Brush™ Antimicrobial Reusable Twisted Wire Brushes are designed for use with Frazier/Baron suction tubes. The twisted wires have a fine/medium diameter yet are rigid enough to ensure excellent cleaning capabilities.

Our wide assortment of Nylon Twisted Wire Brushes are ideal for large bore, uniquely shaped, and endoscopic instruments.

GOOD TO KNOW: OUR TWISTED  
WIRE BRUSHES FEATURE  
A PROTECTIVE BLUE TIP  
FOR ADDED SAFETY.



## 10" / 250 MM FRAZIER /BARON SUCTION TUBE BRUSHES

ITEM #	SIZE	BRUSH LENGTH	BRUSH (Ø)
<b>21-10039</b>	3 FR	0.5" / 12.7 mm	0.039" / 1 mm
<b>21-10067</b>	5 FR	0.5" / 12.7 mm	0.065" / 1.67 mm
<b>20-10059</b>	5 FR	1" / 25.4 mm	0.059" / 1.50 mm
<b>20-10079</b>	6 FR	1" / 25.4 mm	0.079" / 2 mm
<b>20-10090</b>	7 FR	1" / 25.4 mm	0.090" / 2.3 mm
<b>20-10106</b>	8 FR	1" / 25.4 mm	0.106" / 2.7 mm
<b>20-10118</b>	9 FR	1" / 25.4 mm	0.118" / 3 mm
<b>20-10129</b>	10 FR	1.5" / 38 mm	0.129" / 3.3 mm
<b>20-10146</b>	11 FR	1.5" / 38 mm	0.0146" / 3.7 mm
<b>20-10157</b>	12 FR	1.5" / 38 mm	0.157" / 4 mm

PACKAGING: Bag of 3

## 12" / 305 MM NYLON TWISTED WIRE BRUSHES

ITEM #	BRUSH LENGTH	BRUSH (Ø)
<b>22-12078</b>	2" / 50.8 mm	0.078" / 2 mm
<b>22-12118</b>	2" / 50.8 mm	0.118" / 3 mm
<b>22-12157</b>	3" / 76 mm	0.157" / 4 mm
<b>22-12197</b>	3" / 76 mm	0.197" / 5 mm
<b>22-12236</b>	3" / 76 mm	0.236" / 6 mm
<b>22-12393</b>	3" / 76 mm	0.394" / 10 mm

PACKAGING: Bag of 3



## 16" / 406 MM NYLON TWISTED WIRE BRUSHES

ITEM #	BRUSH LENGTH	BRUSH (Ø)
<b>23-16197</b>	2.5" / 63.5 mm	0.197" / 5 mm
<b>23-16236</b>	2.5" / 63.5 mm	0.236" / 6 mm
<b>23-16276</b>	3" / 76 mm	0.276" / 7 mm
<b>23-16394</b>	3" / 76 mm	0.394" / 10 mm
<b>23-16472</b>	3" / 76 mm	0.472" / 12 mm
<b>23-16591</b>	3" / 76 mm	0.591" / 15 mm
<b>23-16787</b>	3" / 76 mm	0.787" / 20 mm

PACKAGING: Bag of 3

## 18" / 457 MM NYLON TWISTED WIRE BRUSHES

ITEM #	BRUSH LENGTH	BRUSH (Ø)
<b>19-18010</b>	0.47" / 12 mm	0.039" / 1 mm
<b>19-18015</b>	0.5" / 12 mm	0.06" / 1.6 mm

PACKAGING: Bag of 3

## 24" / 610 MM NYLON TWISTED WIRE BRUSHES

ITEM #	BRUSH LENGTH	BRUSH (Ø)
<b>24-24098</b>	2.5" / 63.5 mm	0.098" / 2.5 mm
<b>24-24197</b>	2.5" / 63.5 mm	0.197" / 5 mm
<b>24-24276</b>	3" / 76 mm	0.276" / 7 mm
<b>24-24394</b>	3" / 76 mm	0.394" / 10 mm
<b>24-24472</b>	3" / 76 mm	0.472" / 12 mm

PACKAGING: Bag of 3



## BRONCHO BRUSHES

ITEM #	DESCRIPTION	BRUS LENGTH	BRUSH (Ø) INCHES / MM
<b>26-06175</b>	Small	6" / 152 mm	1.75" / 44.5 mm
<b>26-08325</b>	Large	8" / 203 mm	3.25" / 82.6 mm

PACKAGING: Bag of 3

## TRACHEOSTOMY BRUSHES

ITEM #	DESCRIPTION	TOTAL LENGTH	BRUSH LENGTH	BRUSH (Ø)
<b>25-08500</b>	Large 6-12.3 3 / 4"	8" / 203 mm	3.75" / 95 mm	0.5" / 12.7 mm
<b>25-37550</b>	Small 000-5.1 7 / 8"	4" / 93 mm	1.88" / 48 mm	0.5" / 12.7 mm

PACKAGING: Bag of 3

BATRIK

# STAINLESS STEEL ARC™ MICRO BRUSHES

Our Stainless Steel ARC™ Micro Brushes are some of the smallest brushes in the industry. They are designed to clean channels with an inner diameter ranging from 0.5 mm to 1.40 mm.

- Flexible stainless steel sheath will not bend like the regular metal twisted wire brushes
- Polyamide bristles prevents or slows bacterial growth between sterilizations
- Reusable\*
- Available in two lengths

TIP: USING THE PULL-THROUGH  
TECHNIQUE IS HIGHLY  
RECOMMENDED WITH THESE  
BRUSHES. SIMPLY INSERT THE  
BLUNT, NON-BRUSH END FIRST.



## PRODUCT SELECTION

ITEM #	TOTAL LENGTH	BRUSH SIZE (Ø)	BRUSH LENGTH	SHEATH (Ø)	PACKAGING
<b>35-03015</b>	12" / 30 cm	0.06" / 1.5 mm	0.6" / 15 mm	0.02" / 0.40 mm	6 / Tube
<b>35-04515</b>	18" / 45 cm	0.06" / 1.5 mm	0.6" / 15 mm	0.02" / 0.40 mm	6 / Tube

\*Although all materials used in the fabrication of these brushes can be re-processed repeatedly, due to the delicate nature of these brushes re-use is conditional on their handling and application.

## BATRIK INSTRUMENT CLEANING FLAT BRUSHES

Our Instrument Cleaning Flat Brushes are specifically designed for hospital use in the pre-clean and preparation of instruments.

- Long ergonomic handle permits ease-of-use
- Extra long length of brush head (60 [±5] mm in Single-Use and 75 [±5] mm in Reusable)
- Transparent brush fibers allow for quick verification of its cleanliness
- The plastic handle will not release any toxic components
- Available in 6 models (4 reusable/autoclavable and 2 single-use)



### FLEXIBLE

This brush is designed for large instruments that require a more gentle cleaning application such as optics, plastics, or hard to reach spot.



### EXTRA FLEXIBLE

The very long, flexible bristles eliminate any stress during sanitation, and makes this the ideal brush for fragile instruments as well as both pliable and rigid endoscopes. The length of the bristles assist in cleaning areas that would otherwise be difficult to reach, such as the accessory channel or tip deflection of pliable endoscopes.



### RIGID

The rigidity of this brush, along with the length of the bristles and handle, improves and simplifies instrument cleaning.



### EXTRA RIGID

Instrument safe polyamide fibers offer steel-like rigidity and allows deep cleaning of soiled instruments, particularly in grooved or uneven surfaces.



50-40919



50-40922



50-40920



50-40918



50-40917



50-40921

## PRODUCT SELECTION


ITEM #	BRUSH TYPE	DESCRIPTION	COLOR	BRUSH LENGTH	BRUSH SIZE	PACKAGING
50-40917	Extra Rigid	Reusable Autoclavable Flat Brush	Blue	220 mm	75 mm x 20 mm x 10 mm	Bag of 3
50-40918	Rigid	Reusable Autoclavable Flat Brush	Green	220 mm	75 mm x 20 mm x 15 mm	Bag of 3
50-40919	Flexible	Reusable Autoclavable Flat Brush	Yellow	220 mm	75 mm x 20 mm x 15 mm	Bag of 3
50-40920	Extra Flexible	Reusable Autoclavable Flat Brush	Red	220 mm	75 mm x 20 mm x 25 mm	Bag of 3
50-40921	Rigid	Single-Use Flat Brush	Light Blue	220 mm	60 mm x 20 mm x 15 mm	Bag of 20
50-40922	Flexible	Single-Use Flat Brush	Orange	220 mm	60 mm x 20 mm x 15 mm	Bag of 20

**BATRIK**

# REUSABLE ORTHOPEDIC RASPS CLEANING BRUSH

Our Reusable Orthopedic Rasps Cleaning Brush is used to clean orthopedic rasps and other similar instruments used in orthopedic surgery.

- Long lasting and can withstand many sterilization cycles
- Ergonomic steamable thermoplastic handle
- Antibacterial polyamide bristles
- Reusable
- Latex free



## PRODUCT DESCRIPTION

ITEM #	TOTAL LENGTH	HANDLE LENGTH	BRUSH SIZE	PACKAGING
<b>50-RASP</b>	9" / 23 cm	5.5" / 14 cm	0.78" x 2.7" x 3.5" 2 cm x 7 cm x 9 cm	3 / Bag


BATRIK

# REUSABLE PRECISION BRUSH

## INSTRUMENT CLEANING

Our Reusable Precision Instrument Cleaning Brush is the ideal brush for the quick removal of difficult debris. It cleans externally as well as deep into instrument jaws. Although the firm and durable bristles feel like stainless steel, there is absolutely no risk of instrument damage or personal injury.

- Bristles have an extra firm density and superior durability
- Polyamide (TYNEX®) bristles
- Steamable thermoplastic ergonomic handle
- Steam / chemical sterilizable
- Latex free



GOOD TO KNOW: ELIMINATES  
THE NEED FOR THE STAINLESS  
STEEL TOOTHBRUSH-STYLE BRUSH,  
WHICH CAN DAMAGE VALUABLE  
INSTRUMENTS!

## PRODUCT SELECTION

ITEM #	DESCRIPTION	TOTAL LENGTH
<b>50-40800</b>	Reusable Precision External Cleaning Brush	6.29" / 160 mm
PACKAGING: Bag of 3		






**BATRIK**

# REUSABLE INSTRUMENT CLEANING BRUSHES

Designed specifically for hospital applications, our Reusable Instrument Cleaning Brushes feature long and flexible bristles to maximize the effectiveness of large-surface cleaning.

- High quality polyamide clear bristles
- Easy to identify debris
- Autoclavable
- Latex free

**REUSABLE JAR/BOTTLE BRUSH****STERILIZER DRAIN BRUSH**

## PRODUCT SELECTION

ITEM #	DESCRIPTION	HANDLE LENGTH	BRUSH LENGTH	BRUSH (Ø)	TOTAL LENGTH
<b>50-10295</b>	Reusable Jar / Bottle Brush	120 mm ± 2 mm	120 mm ± 5 mm	7 mm ± 2 mm	13.19" / 355 mm
<b>50-10224</b>	Sterilizer Drain Cleaning Brush	N/A	4" / 10.1 cm	2" / 5.08 cm	24" / 609 mm
PACKAGING: Bag of 3					

BATRIK

# ACETABULAR REAMER BRUSH

Our double-ended cleaning brush is manufactured exclusively for cleaning acetabular bone reamers. The “Pig Tail” rounded style circular head brush is designed to effectively clean inside the reamer, while the rigid straight brush on the reverse side is used for cleaning individual cutting holes. Special fusion of the brush heads to the hand piece along with its high quality antimicrobial medical grade bristles ensures a durable, reusable, and highly effective cleaning tool.



STRAIGHT BRUSH



ROUND BRUSH



## PRODUCT SELECTION

ITEM #	DESCRIPTION	HANDLE LENGTH	STRAIGHT BRUSH	ROUND BRUSH HEAD
<b>50-ACETAB</b>	Acetabular Reamer Cleaning Brush	4.72\" / 120 mm Ø 0.47\" / 12 mm	1.18\" / 30 mm Ø 0.3125\" / 8 mm	2\" / 50 mm Ø 1.75\" / 45 mm

PACKAGING: Bag of 3

NOTE: Specifications subject to small variances

**BATRIK****ALL PURPOSE REUSABLE BRUSHES****Plastic Handle Brushes**

We offer a wide selection of brushes for instrument cleaning applications. These brushes are ideal for large bore, uniquely shaped, and endoscopic instruments.

**In-Brush™ Plastic Handle**

Available in regular toothbrush style (nylon or stainless) and double-ended brushes.

**28-72550****28-69394****29-72550****29-69394****29-70050****50-40911****PRODUCT SELECTION**

ITEM #	DESCRIPTION	LENGTH TOTAL	BRUSH LENGTH	BRUSH HEIGHT
<b>28-72550</b>	Plastic Single End Nylon Reusable Handle Brush	7.25" / 184 mm	1.5" / 38 mm	0.5" / 12.7 mm
<b>29-72550</b>	Plastic Single End Stainless Steel Reusable Handle Brush	7.25" / 184 mm	1.5" / 38 mm	0.5" / 12.7 mm
<b>29-70050</b>	Plastic Single End Brass Reusable Handle Brush	7.25" / 184 mm	1.5" / 38 mm	0.5" / 12.7 mm
<b>28-69394</b>	Plastic Double End Nylon Reusable Handle Brush	6.9" / 175 mm	1.38" / 35 mm 1" / 25 mm	0.39" / 10 mm 0.2" / 5 mm
<b>29-69394</b>	Plastic Double End Stainless Steel Reusable Handle Brush	6.9" / 175 mm	1.38" / 35 mm 1" / 25 mm	0.39" / 10 mm 0.2" / 5 mm
<b>50-40911</b>	Hand Scrub Nylon Brush Reusable autoclavable	4.09" / 104 mm	3.74" / 95 mm	1.1" / 28 mm

**PACKAGING:** Bag of 3 - For a bag of 50, add "-50" following the item code.

BATRIK

# AUTOCLAVE CLEANING BRUSH KITS



## ROUND BRUSH

- Features 3 holes for adjustable pole placement
- Adjustable brush placement (parallel or perpendicular) and longer brush length provides easy access to hard-to-reach corners
- Made with FDA approved polyester bristles

## HAND BRUSH

- Ergonomic handle for better hand comfort
- Longer, FDA approved polyester bristles
- Easy to clean (no cavity)

## PRODUCT SELECTION

### REGULAR AUTOCLAVE CLEANING BRUSH KIT

ITEM # (KIT)	DESCRIPTION	ITEM # (UNIT)	DESCRIPTION	DIMENSIONS
<b>50-CLAVEREG</b>	Autoclave REGULAR Brush Cleaning Kit	<b>50-35900</b>	1 Hand Brush	Length: 3.5" x 8" / 8.89 cm x 20.32 cm
		<b>50-55800</b>	1 Rounded Corner Brush (Head Brush)	Head size: 5.5 x 10" / 14 x 25.4 cm
		<b>50-EXTPOLE</b>	1 Anodized Aluminum Extension Pole	2-4 ft / 60-122 cm

### SUPERIOR AUTOCLAVE CLEANING BRUSH KIT

ITEM # (KIT)	DESCRIPTION	ITEM # (UNIT)	DESCRIPTION	DIMENSIONS
<b>50-CLAVESUP</b>	Autoclave SUPERIOR Brush Cleaning Kit	<b>50-35900</b>	1 Hand Brush	Length: 3.5" x 8" / 11.5 cm x 20.32cm
		<b>50-55800</b>	1 Rounded Corner Brush (Head Brush)	Head size: 5.5 x 10" / 14 x 25.4 cm
		<b>50-EXTPOLE</b>	1 Anodized Aluminum Extension Pole	2-4 ft / 60-122 cm
		<b>SF-100</b>	2 SAVEFACE™ Face Shields	9" / 23 cm

## REPLACEMENT COMPONENTS



**HAND BRUSH**  
**50-35900**

Packaging: 1



**BRUSH HEAD**  
**50-55800**

Packaging: 1



**EXTENSION POLE**  
**50-EXTPOLE**

Packaging: 1



**SAVEFACE™**  
**FACE SHIELDS SF-100**

Packaging: 100 / Box

**BATRIK**

# HYSTEROSONOGRAPHY CATHETER

Hystero-sonography (also known as sonohysterography) can now be performed safely, simply and atraumatically with the gentle infusion of saline into the uterus with the use of our H/S Catheter.


- Available in two sizes: 5FR and 7FR
- Minimally invasive and ideal for office use
- Gentle placement sheath included (decreasing the need for a tenaculum) to minimize patient discomfort
- Flexible and easy-to-use
- Total length of 410mm with a working length of 280 mm
- Latex Free

GOOD TO KNOW: ACHIEVE  
IMPROVED IMAGERY  
AND DIAGNOSIS WITH  
SONOHYSTEROGRAMS

“We depend upon the reliability, precision and functionality of the Batrik H/S Catheter in our clinical and research operations.

Having performed more than 100,000 sonohysterography procedures, we can unequivocally say that this is the best HSG catheter that we have ever used

”  
-Dr. Alex Hartman



## PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING
<b>51-5005</b>	5 FR H/S Catheter	10 / box
<b>51-5007</b>	7 FR H/S Catheter	10 / box

NOTES




Lined area for notes, consisting of 20 horizontal lines.



100

This image shows a full page of white paper with horizontal blue or grey ruling lines. The lines are evenly spaced and run across the width of the page, typical of notebook paper. There are no margins, text, or other markings on the page.



SURGMED GROUP

850 Halpern Ave., (Dorval) Montreal QC H9P 1G6  
Phone: +1 514-631-7988 | Toll Free: 800-661-5432 | Fax.: +1 514-631-9083  
Email: [info@surgmed.com](mailto:info@surgmed.com) | Website: [www.batrik.com](http://www.batrik.com)