

PROHS®

Horizontal Steam Sterilizer PJ PREMIUM VERSION

CE 1155 PED

PROHS Equipamento Hospitalar
e Serviços Associados, S.A.

Rua do Castanhal 316
Zona Industrial Maia I, sector II
4475-122 Maia, Portugal

| T | +351 229 059 170
| F | +351 229 015 900
| E | prohs@prohs.pt
| W | www.prohs.pt




Horizontal Steam Sterilizer PJ

PREMIUM VERSION

The **PREMIUM Version** of PROHS Horizontal Steam Sterilizer PJ brings, **elegance** through its **innovative** and **futuristic** design. This version enhances the **user interaction** with the machine, including, as standard, a SD card with **PC connection software** and a **second touch screen** on the unloading side, taking the equipment interface to the next level.

The interior was also the subject of the utmost attention, the **steam pipes** were upgraded to AISI 316L **stainless steel** while the doors and chamber are delivered with a **finer polishing**, reducing the roughness to 0.8 micron, leading to a splendid finish.

This version is always manufactured with partial jacket chamber and a new exclusive frontal design to offer the user a truly **PREMIUM** experience.



Model	Useful Dimensions	External Dimensions	Energy		Water Consumption	Weight	Useful Space
	/cm	/cm	/with Steam Generator	/without Steam Generator	/litres (per cycle)	/Kg (aproximate, with steam generator)	/STU
PROHS PJ 190L	50x50x75	180x102x100	21kW	2,1kW	120	495	1
PROHS PJ 250L	50x50x100	180x102x130	32kW	2,1kW	130	570	1
PROHS PJ 325L	50x50x130	180x102x160	32kW	2,1kW	145	685	2
PROHS PJ 350L	68x68x75	195x122x105	32kW	2,1kW	170	750	4
PROHS PJ 360L	60x60x100	180x112x134	32kW	2,1kW	170	700	1
PROHS PJ 460L	68x68x100	195x122x134	42kW	2,1kW	170	800	6
PROHS PJ 600L	68x68x130	195x122x164	53kW	2,8kW	200	830	8
PROHS PJ 740L	68x68x160	195x122x194	53kW	2,8kW	225	1000	10
PROHS PJ 880L	68x68x190	195x122x224	68kW	2,8kW	250	1150	12



Load View



Side View



Unload View

Due to product improvement and development, PROHS reserves the right to make changes without notice.