

## NICU | PICU




CLARAVUE is the best-in-class, hygienic solution for simple, accurate and fully radiolucent patient monitoring. CLARAVUE includes pre-wired and pre-positioned disposable electrodes and patented reusable bioactive trunk cables.

For the past ten years, CLARAVUE has been the proven choice for superior trace quality and infection prevention in prominent European hospitals.

### Lightweight, disposable, pre-wired electrodes for the NICU | PICU

#### FEATURES

- 5-day wear time means fewer disruptions ensuring the body is focused on essential growth and repair
- Lightweight, disposable leadwires for use on even the tiniest of patients, reducing the risk of cross contamination
- Soft and gentle hydrogel provides comfortable contact without the use of harsh acrylic adhesives
- Pre-wired electrodes helps to eliminate motion artifact



#### PACKAGING

Electrodes are packaged in a single-use pouch with an easy-to-open tear

##### Box




##### Single-Use Kit



FRONT

BACK



## CLARAVUE RADIOLUCENT DISPOSABLE PRE-WIRED ELECTRODES

Product #	Description	Per Kit	Kits/Box	Kits/Case	Solid Gel	Duration of Use
4009832C	Pediatric Patient Kit - 3 Electrodes	3	40	480		

## COMPATIBLE MONITOR CONNECTIONS

Compatible Cables include adaptors for Bard, Biosense, Mindray, Datascope, GE Marquette, Philips, Datex, Mennen, Mortara, Nihon Kohden, GE, Getemed, Spacelabs, Siemens, Drager, Schiller, Zoll and more.



Mindray



Philips



Spacelabs  
TruLink



GE



Physio



Nihon  
Kohden

## CABLE CONNECTIONS

Electrode → Patented Electrode Cable Connector → Cable → Monitor

