

WarwickSASCo
Moulding The Future of Medical Plastics



5 reasons to choose Warwick SASCo

Recognised as the market leader in the design and manufacture of medical polyware

Wide range of high quality re-usable and single-use polypropylene products for all areas of healthcare

Unrivalled reputation for providing quality products, delivery and service to healthcare professionals worldwide

Re-usable products are guaranteed for a minimum of 1000 wash cycles – delivering excellent value for money

Certification since 1997 demonstrates our continued commitment to maintaining industry standards. Our quality management system complies to ISO13485:2016. Products comply with (EU)2017/745 Medical Device Regulation Class 1 and show CE where applicable.



Warwick SASCo

For over 40 years Warwick SASCo has been part of evolving supply of decontamination and nursing products worldwide to hospitals, healthcare and patient nursing.

Recognised as a market leader, we have a reputation for high quality products designed for hospital decontamination and sterile services, operating and surgical applications, patient care and home nursing - where infection management and good user application is paramount.

Warwick SASCo's low cost reusable and single use product ranges are chosen by healthcare professionals and end users worldwide.

Our growing family owned business is proud to be represented by leading distributor partners in over 70 countries.

All our products are manufactured within the UK with Quality Assured Suppliers. Our customer support and distribution centre is in Warwick, England.


The team at Warwick SASCo is immensely proud of delivering high quality products and a professional friendly service to you - our customer. You will find our experienced team are knowledgeable, caring and quick to respond to your requirements.

A handwritten signature in blue ink, appearing to read 'Darby Booth'.

Darby Booth,
Managing Director – Warwick SASCo

Products Overview


Our wide range of reusable and disposable medical polyware products are chosen by healthcare professionals for their high quality and low cost.



Reusable Medical Polyware
WARWICK RANGE

The Warwick Range for use in Decontamination Services, Operating Rooms, Surgeries, Clinics and healthcare services. Products are designed for repeated washing, disinfection and sterilization.

The Warwick Range of products can be processed by steam autoclaving at 135°C, Ethylene Oxide (EO) and Hydrogen Peroxide (H₂O₂). Products comply with (EU)2017/745 Medical Device Regulation Class 1 and show CE where applicable. Material used is polypropylene unless otherwise stated.



Single Use Medical Polyware
STERILITE RANGE

Sterilizable lightweight polyware for single-use and disposable sets. Sterilite Range is designed for Custom Pack Manufacturers and Decontamination Services.

The Sterilite Range of products can be processed by steam autoclaving at 135°C, Ethylene Oxide (EO) and Hydrogen Peroxide (H₂O₂). Products comply with (EU)2017/745 Medical Device Regulation Class 1 and show CE where applicable. Material used is polypropylene unless otherwise stated.



Nursing & Care Daily Living Aids & Accessories
SPA RANGE

Spa Range of personal living aids, patient care, toileting and bedside products.

Products comply with BS2588 For Portable Sanitary Products, BS3215 for Reusable Urinals, (EU)2017/745 Medical Device Regulation Class 1 and show CE where applicable.

Material used is polypropylene unless otherwise stated.



How to Order

To place an order, request a quotation or for general enquiries:

CALL us on +44 (0) 1926 422 427
 FAX +44 (0) 1926 315 405
 EMAIL enquiries@sasco.co.uk


We work with a wide range of industry sectors including:



STERILE SERVICES & DECONTAMINATION
Operating Room, Instrument Packing




PERSONAL USE & HOME CARE
Daily Living Aids, Community Nurse, Incontinent Care



HOSPITALS
Surgery, Clinics, Wards, Veterinary, Dentistry



NURSING HOMES
Patient Care, Rehabilitation



CUSTOM PACK MANUFACTURERS
Single Use Packs, Contract Instrument Processing

WARWICK - reusable range

Compartment Trays	page	6
Disinfectant & Soaking Containers	page	6
Instrument Trays	page	7-8
DIN Style Baskets	page	8
Instrument Tray Lids	page	9
Gallipots	page	9
Lotion & Sponge Bowls	page	10
Jugs	page	10
Kidney Dishes	page	11
Theatre Bowls	page	11
Wash Bowls	page	11
Medical Boxes, Instrument Boxes & Storage Containers	page	12
Scope Containers	page	13
Silicone Protection	page	13
Tray Tags	page	14

STERILITE - single use range

Kidney Dishes	page	16
Compartment Trays	page	17
Gallipots	page	17
Quivers	page	18
Sponge Bowls	page	18
Theatre Bowls	page	18

SPA - care and living aid range

Bedpans	page	20-21
Bedpan Holders	page	21
Commode Pans	page	21
Bidets	page	21
Urinal Bottles and Urinal Pans	page	22
Urinal Accessories	page	22
Drinking Beakers and Cups	page	23
Drinking Cup Lids	page	23
Denture Cups	page	23
Jugs and Jug Sets	page	24
Medicine Measures	page	24
Patient Trays	page	25
Kidney Dishes	page	25
Vomit Bowls	page	25



To find out more about our service, product ranges, product decontamination cycle and industry news:

visit www.sasco.co.uk

WARWICK RANGE

Reusable Medical Polyware

The Warwick Range is a wide range of sterilizable long-life plastic polyware products.

Designed as an alternative to stainless steel products, they are ideal for use in Decontamination Services, Operating Rooms, Clinics, Surgeries and other surgical or healthcare services.

Products are designed for repeated washing, disinfection and sterilization. The Warwick Range of products can be processed by steam autoclaving at 135°C, Ethylene Oxide (EO) and Hydrogen Peroxide (H₂O₂). Products comply with (EU)2017/745 Medical Device Regulation Class 1 and show CE where applicable.

Material used is polypropylene unless otherwise stated.

Warwick Benefits

- Low cost.
- Robust and hard wearing - guaranteed for 1000 cycles.
- Easy to clean - blood protein is easily removed from the surface, leaving no residue.
- 50% lighter than stainless steel, assisting users with storage, lifting and handling.
- Can be sterilized by steam autoclaving at 135°C, Ethylene Oxide (EO) and Hydrogen Peroxide (H₂O₂).
- Can be used as standalone items or as a support component in surgical sets.
- Widely used in Sterile Services for packing instrument sets and procedure kits for the autoclave.
- Disinfection sets can be used for manual sterilization in Sterile Services, Clinics, etc.
- Used non-sterile in hospital wards and clinics.
- Smooth corners and surfaces, reducing risk of paper rips and the need for re-processing.
- Autoclave sterilization temperature is 134.6°C + 2°C (generally referred to as 135°C).

Compartment Trays, Disinfectant and Soaking Containers



Compartment Tray - 3 compartments

Code: CIT1814

Dimensions: 181x134x54mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 38



Compartment Tray - 3 compartments

Code: CIT2013

Dimensions: 205x130x30mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 120



Compartment Tray - 4 compartments

Code: CIT2718

Dimensions: 270x180x41mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 25



Compartment Tray - 4 compartments

Code: LCT2215

Dimensions: 220x150x24mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 160



Disinfection Set - Soaking container, strainer and hinged lid
Code: DSSW1000

Dimensions: 260x130x65mm
Volume: 1litre
Material: Polypropylene
Colour: White
Case Size: 10



Disinfection Set - Soaking container, strainer and lid
Code: DSS2711

Dimensions: 257x180x128mm
Volume: 3litre
Material: Polypropylene
Colour: Blue
Case Size: 5



Disinfection Set - Soaking container, strainer and lid
Code: DSSW2711

Dimensions: 257x180x128mm
Volume: 3litre
Material: Polypropylene
Colour: White
Case Size: 5



Disinfection Set - Soaking container, strainer and vented lid
Code: DSSV2711

Dimensions: 257x180x128mm
Volume: 3litre
Material: Polypropylene
Colour: Blue
Case Size: 5




Disinfection Set - Soaking container, strainer, drip tray and lid
Code: DSS4215

Dimensions: 426x156x138mm
Volume: 4litre
Material: Polypropylene
Colour: Blue
Case Size: 1



Disinfection Set - Soaking container, strainer, drip tray and lid
Code: DSS8515

Dimensions: 890x200x145mm
Volume: 12litre
Material: Polypropylene
Colour: Blue
Case Size: 1



Disinfection Set - 2 x soaking containers, 1 x strainer, 2 x lids and 1 x drip tray
Code: DSS8515TP

Dimensions: 890x200x145mm
Volume: 12litre
Material: Polypropylene
Colour: Blue
Case Size: 1



R Reusable

W Washer Decontamination

S Suitable for Steam Sterilization

E Suitable for Ethylene Oxide Sterilization

HP Suitable for Hydrogen Peroxide Sterilization

LF Latex Free

Instrument Trays and DIN Style Baskets



Instrument Tray -
Solid base

Code: IT1813

Dimensions: 181x134x54mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 40

R W S E HP LF



Instrument Tray -
Solid base

Code: IT1913

Dimensions: 195x135x23mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 115

R W S E HP LF



Instrument Tray -
Solid base

Code: T2015

Dimensions: 200x150x51mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 60

R W S E HP LF



Instrument Tray -
Solid ribbed base

Code: IT2015

Dimensions: 200x150x51mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 60

R W S E HP LF



Instrument Tray -
Solid base

Code: T3025

Dimensions: 300x250x52mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 25

R W S E HP LF



Instrument Tray -
Solid base

Code: IT2715

Dimensions: 270x150x30mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 60

R W S E HP LF

How to order

To place an order, request
a quotation or for general
enquiries:

CALL us on
+44 (0) 1926 422 427

FAX
+44 (0) 1926 315 405

EMAIL
enquiries@sasco.co.uk



Instrument Tray -
Solid ribbed base

Code: IT3025

Dimensions: 300x250x52mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 25

R W S E HP LF



Instrument Tray -
Solid ribbed base

Code: IT4030

Dimensions: 424x305x75mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 12

R W S E HP LF

Instrument Trays and DIN Style Baskets



Instrument Tray -
Mesh base and sides

Code: MIT1515

Dimensions: 200x180x30mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 50

R W S E HP LF



Instrument Tray -
Mesh base

Code: MIT2015

Dimensions: 200x150x51mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 60

R W S E HP LF



Instrument Tray -
Mesh base and sides

Code: MIT2709

Dimensions: 275x110x30mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 60

R W S E HP LF



Instrument Tray -
Mesh base

Code: MIT2715

Dimensions: 270x150x30mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 60

R W S E HP LF



Instrument Tray -
Mesh base and sides

Code: MIT2817

Dimensions: 280x170x55mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 25

R W S E HP LF



Instrument Tray -
Mesh base and sides

Code: MIT3025

Dimensions: 300x250x52mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 24

R W S E HP LF



Instrument Tray -
Mesh base

Code: MIT4030

Dimensions: 424x305x75mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 12

R W S E HP LF



Instrument Tray -
Perforated base

Code: PIT3025

Dimensions: 300x250x52mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 25

R W S E HP LF



Instrument Tray -
Perforated base

Code: PIT4030

Dimensions: 424x305x75mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 12

R W S E HP LF



DIN Style Tray -
Mesh base and sides

Code: DINT2525

Dimensions: 255x255x50mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 20

R W S E HP LF



DIN Style Tray -
Mesh base and sides

Code: DINT4825

Dimensions: 480x255x50mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 25

R W S E HP LF

R Reusable W Washer Decontamination S Suitable for Steam Sterilization
E Suitable for Ethylene Oxide Sterilization HP Suitable for Hydrogen Peroxide Sterilization LF Latex Free

R Reusable W Washer Decontamination S Suitable for Steam Sterilization
E Suitable for Ethylene Oxide Sterilization HP Suitable for Hydrogen Peroxide Sterilization LF Latex Free

Instrument Tray Lids



Instrument Tray Lid -
For trays IT1813/CIT1814

Code: ITL181

Dimensions: N/A
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 10

R W S E HP LF



Instrument Tray Lid -
For trays IT2015/MIT2015/T2015
Code: ITL201

Dimensions: N/A
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 10

R W S E HP LF



Instrument Tray Lid -
For trays IT2715/MIT2715
Code: ITL271

Dimensions: N/A
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 10

R W S E HP LF



Instrument Tray Lid -
For tray CIT2718
Code: ITL281

Dimensions: N/A
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 10

R W S E HP LF



Instrument Tray Lid -
For tray MIT2817
Code: ITL287

Dimensions: N/A
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 10

R W S E HP LF



Instrument Tray Lid -
For trays IT3025/MIT3025/
PIT3025/T3025
Code: ITL302

Dimensions: N/A
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 10

R W S E HP LF

Gallipots



Gallipot

Code: GP40

Dimensions: 40x30mm
Volume: 35ml
Material: Polypropylene
Colour: Blue
Case Size: 100

R W S E HP LF



Gallipot

Code: GP60

Dimensions: 60x40mm
Volume: 95ml
Material: Polypropylene
Colour: Blue
Case Size: 300

R W S E HP LF



Gallipot

Code: GP80

Dimensions: 80x60mm
Volume: 280ml
Material: Polypropylene
Colour: Blue
Case Size: 80

R W S E HP LF



Gallipot - Flat base
Graduated

Code: GPLB60

Dimensions: 60x40mm
Volume: 50ml
Material: Polypropylene
Colour: Blue
Case Size: 400

R W S E HP LF



Gallipot - Flat base
Graduated

Code: GPLB80

Dimensions: 80x60mm
Volume: 200ml
Material: Polypropylene
Colour: Blue
Case Size: 200

R W S E HP LF

Lotion and Sponge Bowls

How to order

To place an order, request
a quotation or for general
enquiries:

CALL us on
+44 (0)1926 422 427

FAX
+44 (0)1926 315 405

EMAIL
enquiries@sasco.co.uk



Lotion Bowl - Flat base
Graduated

Code: GPLB100

Dimensions: 100x45mm
Volume: 300ml
Material: Polypropylene
Colour: Blue
Case Size: 135

R W S E HP LF



Lotion Bowl - Flat base
Graduated

Code: GPLB150

Dimensions: 150x70mm
Volume: 500ml
Material: Polypropylene
Colour: Blue
Case Size: 60

R W S E HP LF



Sponge Bowl - Flat base
Graduated

Code: SB140

Dimensions: 140x65mm
Volume: 300ml
Material: Polypropylene
Colour: Blue
Case Size: 160

R W S E HP LF



Lotion Bowl
Graduated

Code: LB100

Dimensions: 100x45mm
Volume: 280ml
Material: Polypropylene
Colour: Blue
Case Size: 200

R W S E HP LF



Lotion Bowl
Graduated

Code: LB150

Dimensions: 150x70mm
Volume: 900ml
Material: Polypropylene
Colour: Blue
Case Size: 60

R W S E HP LF



Lotion Bowl
Graduated

Code: LB200

Dimensions: 200x90mm
Volume: 2000ml
Material: Polypropylene
Colour: Blue
Case Size: 30

R W S E HP LF



Lotion Bowl
Graduated

Code: LB250

Dimensions: 250x105mm
Volume: 3000ml
Material: Polypropylene
Colour: Blue
Case Size: 20

R W S E HP LF

Jugs



Jug
Graduated

Code: JG500

Dimensions: 95mm x 95mm
Volume: 500ml
Material: Polypropylene
Colour: Natural
Case Size: 50

R W S E HP LF



Jug
Graduated

Code: JG1000

Dimensions: 120mm x 125mm
Volume: 1000ml
Material: Polypropylene
Colour: Natural
Case Size: 45

R W S E HP LF



Jug
Graduated

Code: JG2000

Dimensions: 130mm x 200mm
Volume: 2000ml
Material: Polypropylene
Colour: Natural
Case Size: 20

R W S E HP LF



Kidney Dishes, Theatre Bowls and Wash Bowls



Kidney Dish
Graduated

Code: KD150

Dimensions: 150x35mm
Volume: 150ml
Material: Polypropylene
Colour: Blue
Case Size: 210

R W S E HP LF



Kidney Dish
Graduated

Code: KD200

Dimensions: 200x45mm
Volume: 500ml
Material: Polypropylene
Colour: Blue
Case Size: 150

R W S E HP LF



Kidney Dish
Graduated

Code: KD250

Dimensions: 250x55mm
Volume: 750ml
Material: Polypropylene
Colour: Blue
Case Size: 60

R W S E HP LF



Kidney Dish
Graduated

Code: KDY250

Dimensions: 250x55mm
Volume: 750ml
Material: Polypropylene
Colour: Yellow
Case Size: 60

R W S E HP LF



Kidney Dish
Graduated

Code: KD300

Dimensions: 300x65mm
Volume: 1500ml
Material: Polypropylene
Colour: Blue
Case Size: 25

R W S E HP LF



Kidney Dish Tray -
Shallow kidney shaped tray

Code: KDT185

Dimensions: 185x92mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 70

R W S E HP LF



Wash Bowl
Graduated

Code: WBN3112

Dimensions: 310x105mm
Volume: 3000ml
Material: Polypropylene
Colour: Natural
Case Size: 25

R W S E HP LF



Wash Bowl
Graduated

Code: WBN3513

Dimensions: 350x120mm
Volume: 5000ml
Material: Polypropylene
Colour: Natural
Case Size: 25

R W S E HP LF



Theatre Bowl
Graduated

Code: WB3112

Dimensions: 310x105mm
Volume: 3000ml
Material: Polypropylene
Colour: Blue
Case Size: 25

R W S E HP LF



Theatre Bowl
Graduated

Code: WB3513

Dimensions: 350x120mm
Volume: 5000ml
Material: Polypropylene
Colour: Blue
Case Size: 25

R W S E HP LF

R Reusable

W Washer Decontamination

S Suitable for Steam Sterilization

E Suitable for Ethylene Oxide Sterilization

HP Suitable for Hydrogen Peroxide Sterilization

LF Latex Free

Medical Boxes, Instrument Boxes & Storage Containers



Instrument Box -
Storage box with lid

Code: IBX2718

Dimensions: 257x180x128mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 10

R W S E HP LF



Instrument box -
Storage box with vented lid

Code: IBXV2718

Dimensions: 257x180x128mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 10

R W S E HP LF



Instrument Box -
Storage box with lid

Code: IBXN290

Dimensions: 300x224x140mm
Volume: N/A
Material: Polypropylene
Colour: Natural
Case Size: 12

R W S E HP LF

How to order

To place an order, request
a quotation or for general
enquiries:

CALL us on
+44 (0) 1926 422 427

FAX
+44 (0) 1926 315 405

EMAIL
enquiries@sasco.co.uk



Transportation Box -
With hinged lid

Code: MB2318

Dimensions: 230x170x200mm
Volume: N/A
Material: ABS
Colour: Orange
Other colours available
Case Size: 10

R W S E HP LF

R Reusable

W Washer Decontamination


S Suitable for Steam Sterilization

E Suitable for Ethylene Oxide Sterilization

HP Suitable for Hydrogen Peroxide Sterilization

LF Latex Free

Scope Containers



Scope Container -
With silicone dividers and lid

Code: SMT5510

Dimensions: 660x110x66mm
Volume: N/A
Material: Polypropylene
Colour: Lilac
Case Size: 1

R W S E HP LF




Silicone Plug -
For SMT5510

Code: SG5510

Dimensions: N/A
Volume: N/A
Material: Silicone
Colour: Pink
Case Size: 10

R W S E HP LF


Silicone Protection



Silicone Pad -
For trays IT/MIT2715 and
small metal DIN tray
Code: SP2313

Dimensions: 230x132mm
Volume: N/A
Material: Silicone
Colour: Pink
Case Size: 10

R W S E HP LF




Silicone Pad -
For trays IT/MIT/PIT/T3025

Code: SP2823

Dimensions: 280x230mm
Volume: N/A
Material: Silicone
Colour: Pink
Case Size: 10

R W S E HP LF




Silicone Pad -
Fits ½ tray IT/MIT/PIT4030

Code: SP2719

Dimensions: 270x190mm
Volume: N/A
Material: Silicone
Colour: Pink
Case Size: 10

R W S E HP LF



Silicone Pad -
Fits ½ tray DINT4825

Code: SP4622

Dimensions: 230x220mm
Volume: N/A
Material: Silicone
Colour: Pink
Case Size: 10

R W S E HP LF



Silicone Pad -
For medium metal DIN tray

Code: SP2324

Dimensions: 235x245mm
Volume: N/A
Material: Silicone
Colour: Pink
Case Size: 10

R W S E HP LF




Silicone Pad -
For large metal DIN tray

Code: SP4724

Dimensions: 470x245mm
Volume: N/A
Material: Silicone
Colour: Pink
Case Size: 10

R W S E HP LF



Silicone Roll -
1 metre length with 10
divisible sections
Code: SR10

Dimensions: 100x245mm
Volume: N/A
Material: Silicone
Colour: Pink
Case Size: 1

R W S E HP LF

Tray Tags



Tray Tag

Code: TAG-010

Dimensions: 105x40mm
Volume: N/A
Material: Nylon
Colour: White
Case Size: 100

R W S E HP LF



Tray Tag

Code: TAG-015

Dimensions: 105x40mm
Volume: N/A
Material: Nylon
Colour: Yellow
Case Size: 100

R W S E HP LF



Tray Tag

Code: TAG-020

Dimensions: 105x40mm
Volume: N/A
Material: Nylon
Colour: Blue
Case Size: 100

R W S E HP LF



Tray Tag

Code: TAG-022

Dimensions: 105x40mm
Volume: N/A
Material: Nylon
Colour: Light Green
Case Size: 100

R W S E HP LF



Tray Tag

Code: TAG-025

Dimensions: 105x40mm
Volume: N/A
Material: Nylon
Colour: Orange
Case Size: 100

R W S E HP LF



Tray Tag

Code: TAG-030

Dimensions: 105x40mm
Volume: N/A
Material: Nylon
Colour: Red
Case Size: 100

R W S E HP LF



Tray Tag

Code: TAG-050

Dimensions: 105x40mm
Volume: N/A
Material: Nylon
Colour: Green
Case Size: 100

R W S E HP LF



Tray Tag

Code: TAG-055

Dimensions: 105x40mm
Volume: N/A
Material: Nylon
Colour: Grey
Case Size: 100

R W S E HP LF



Tray Tag

Code: TAG-060

Dimensions: 105x40mm
Volume: N/A
Material: Nylon
Colour: Lavender
Case Size: 100

R W S E HP LF



Tray Tag

Code: TAG-065

Dimensions: 105x40mm
Volume: N/A
Material: Nylon
Colour: Tan
Case Size: 100

R W S E HP LF



Tray Tag

Code: TAG-070

Dimensions: 105x40mm
Volume: N/A
Material: Nylon
Colour: Black
Case Size: 100

R W S E HP LF



Stainless Steel Attachment Ring

Code: TR30

Dimensions: 30mm
Volume: N/A
Material: Steel
Colour: Silver
Case Size: 100

R W S E HP LF

R Reusable

W Washer Decontamination

S Suitable for Steam Sterilization

E Suitable for Ethylene Oxide Sterilization

HP Suitable for Hydrogen Peroxide Sterilization

LF Latex Free

STERILITE RANGE

Single Use Medical Polyware

The Sterilite Range is a wide range of single use plastic polyware products.

Designed to speed up packing and sterilizing procedures. Sterilite is ideal for use in Sterile Pack Manufacturing, Decontamination Services, Operating Rooms and Clinics due to its low cost, product design and performance in hospitals.

The Sterilite Range of products can be processed by steam autoclaving at 135°C, Ethylene Oxide (EO) and Hydrogen Peroxide (H₂O₂).

Products comply with (EU)2017/745 Medical Device Regulation Class 1 and show CE where applicable. Material used is polypropylene unless otherwise stated.

Sterilite Benefits

- Low cost and disposable.
- No need to invest time and money into retrieval and decontamination after use.
- Clean and ready for sterilization.
- Can be sterilized by steam autoclaving at 135°C, Ethylene Oxide (EO) and Hydrogen Peroxide (H₂O₂).
- Can be used as stand-alone items or as a support (supplementary) component in operating rooms and surgical procedures.
- Widely used in Sterile Services for packing instrument sets and procedure kits.
- Used non-sterile in hospital wards and clinics.



Kidney Dishes



Kidney Dish
Graduated

Code: KB250

Dimensions: 252x100x48mm
Volume: 750ml
Material: Polypropylene
Colour: Natural
Case Size: 450



Kidney Dish
Graduated

Code: KB300

Dimensions: 300x140x57mm
Volume: 1500ml
Material: Polypropylene
Colour: Natural
Case Size: 135



Single Use

Washer Decontamination


Suitable for Steam Sterilization

Suitable for Ethylene Oxide Sterilization

Suitable for Hydrogen Peroxide Sterilization

Latex Free

Compartment Tray



Packing Tray -
4 compartments

Code: LCT2215S

Dimensions: 220x150x24mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Case Size: 160



Gallipots



Gallipot
Graduated

Code: GPS60

Dimensions: 62x30mm
Volume: 60ml
Material: Polypropylene
Colour: Natural
Case Size: 2000



Gallipot
Graduated

Code: GPS120

Dimensions: 85x30mm
Volume: 120ml
Material: Polypropylene
Colour: Natural
Case Size: 1000



Gallipot
Graduated

Code: GP100

Dimensions: 100x53mm
Volume: 250ml
Material: Polypropylene
Colour: Natural
Case Size: 1050



Gallipot
Graduated

Code: GP130

Dimensions: 130x66mm
Volume: 500ml
Material: Polypropylene
Colour: Natural
Case Size: 735



Single Use

Washer Decontamination

Suitable for Steam Sterilization

Suitable for Ethylene Oxide Sterilization

Suitable for Hydrogen Peroxide Sterilization

Latex Free

Quiver



Quiver

Code: QV30

Dimensions: 190x65mm
Volume: N/A
Material: Polypropylene
Colour: Natural
Case Size: 75



Quiver Extension -
For QV30

Code: QV47

Dimensions: 75x244mm
12mm(loop)
Volume: N/A
Material: Polypropylene
Colour: Pink
Case Size: 60



Sponge Bowls and Theatre Bowl

Sponge Bowl
Graduated

Code: SB170

Dimensions: 170x84mm
Volume: 1000ml
Material: Polypropylene
Colour: Natural
Case Size: 340



Sponge Bowl
Graduated

Code: SB200

Dimensions: 210x90mm
Volume: 2000ml
Material: Polypropylene
Colour: Natural
Case Size: 220



Theatre Bowl
Graduated

Code: TB350

Dimensions: 352x116mm
Volume: 5000ml
Material: Polypropylene
Colour: Blue
Case Size: 50



How to order

To place an order, request
a quotation or for general
enquiries:

CALL us on

+44 (0) 1926 422 427

FAX

+44 (0) 1926 315 405

EMAIL

enquiries@sasco.co.uk

Single Use

Washer Decontamination

Suitable for Steam Sterilization

Suitable for Ethylene Oxide Sterilization

Suitable for Hydrogen Peroxide Sterilization

Latex Free

SPA RANGE

Nursing & Care, Daily Living Aids and Accessories

The Spa range includes reusable products for personal living, patient care, toileting and bedside comfort.

The Spa Range is designed for use in Hospital Wards, Rehabilitation, Nursing Homes, Community Care and for Personal Use in the Home.

Products comply with BS2588 For Portable Sanitary Products, BS3215 for Reusable Urinals, (EU)2017/745 Medical Device Regulation Class 1 and show CE where applicable.

Material used is polypropylene unless otherwise stated.

Spa Benefits

- Toileting products that preserve patient dignity and independence by empowering patients to toilet themselves with little or no assistance.
 - Urinal Bottles specifically designed to the male and female anatomy for a comfortable and safe fit during use.
 - Product design reduces risk of cross infection, unpleasant odours and sights.
 - Easy to use, clean and disinfect,
- meaning both patient and carer can be confident that products are safe to re-use.
 - Beakers and Drinking Cups are wide, easy to grip, non-scratch, drop-proof and have a nondrip rim.
 - Handled beakers are ideal for patients with unsteady or infirm hands.
- Medicine Measures are sturdy, ideal for long-term use and regularly checked by Trading Standards for measurement accuracy and compliance.
 - All Spa products can be cleaned in automatic washer disinfectors. Also suitable for hand washing using detergents and disinfectants.

Bedpans

Bedpan - Adult

Code: BPN

Dimensions: 370x310x90mm
Volume: 2500ml
Material: Polypropylene
Colour: Green
Case Size: 30

R W S E HP LF



How to order

To place an order, request a quotation or for general enquiries:

CALL us on
+44 (0) 1926 422 427
FAX
+44 (0) 1926 315 405

EMAIL
enquiries@sasco.co.uk



Slipper Pan - Bedpan with handle and lid

Code: FPN

Dimensions: 449x309x99mm
Volume: 1500ml
Material: Polypropylene
Colour: Green
Case Size: 12

R W S E HP LF



Slipper Pan - Bedpan with handle and lid

Code: FPNW

Dimensions: 449x309x99mm
Volume: 1500ml
Material: Polypropylene
Colour: White
Case Size: 12

R W S E HP LF



Bedpan - Hospital style bedpan

Code: VP

Dimensions: 355x308x108mm
Volume: 2500ml
Material: Polypropylene
Colour: Green
Case Size: 18

R W S E HP LF



Bedpan - Hospital bedpan with lid

Code: VPN

Dimensions: 355x308x108mm
Volume: 2500ml
Material: Polypropylene
Colour: Green
Case Size: 15

R W S E HP LF



Bedpan - Hospital bedpan with lid

Code: VPNW

Dimensions: 355x308x108mm
Volume: 2500ml
Material: Polypropylene
Colour: White
Case Size: 15

R W S E HP LF



Bedpan - European bedpan with handle and lid

Code: DPNF

Dimensions: 420x315x100mm
Volume: 1500ml
Material: Polypropylene
Colour: Green
Case Size: 15

R W S E HP LF

R Reusable

W Washer Decontamination

S Suitable for Steam Sterilization

E Suitable for Ethylene Oxide Sterilization

HP Suitable for Hydrogen Peroxide Sterilization

LF Latex Free

Bedpans and Bedpan Holders

Slipper Pan -
Female urinal and
infants bedpan
Code: SPN

Dimensions: 335x245mm
Volume: 1500ml
Material: Polypropylene
Colour: Green
Case Size: 25

R W S E HP LF



Bedpan Support -
For pulp disposable midi
slipper pan 1.3 litres
Code: PSPS

Dimensions: 470x290x75mm
Volume: 1500ml
Material: Polypropylene
Colour: Cream
Case Size: 30

R W S E HP LF



Bedpan Support -
For pulp disposable bedpan
2 litres
Code: PBPS

Dimensions: 370x310x90mm
Volume: 2000ml
Material: Polypropylene
Colour: Cream
Case Size: 30

R W S E HP LF

Commode Pans and Bidets



Commode Pan / Bidet

Code: CMP

Dimensions: 359x338x130mm
Volume: 3000ml
Material: Polypropylene
Colour: Green
Case Size: 20

R W S E HP LF



Commode Pan / Bidet

Code: CMPW

Dimensions: 359x338x130mm
Volume: 3000ml
Material: Polypropylene
Colour: White
Case Size: 20

R W S E HP LF



Commode Pan -
Commode pan with lid

Code: CBL

Dimensions: 259x259x144mm
Volume: 2500ml
Material: Polypropylene
Colour: White
Case Size: 25

R W S E HP LF

How to order

To place an order, request
a quotation or for general
enquiries:

CALL us on

+44 (0)1926 422 427

FAX

+44 (0)1926 315 405

EMAIL

enquiries@sasco.co.uk

Urinal Bottles and Urinal Pans

Urinal Bottle -
Male urinal with handle
Graduated
Code: UML1000

Dimensions: 180x100x85mm
Volume: 1000ml
Material: Polyethylene
Colour: Natural
Case Size: 50

R W S E HP LF



Urinal Bottle -
Male urinal, hospital style
Graduated
Code: UM1000

Dimensions: 170x90x 80mm
Volume: 1000ml
Material: Polypropylene
Colour: Natural
Pack Size: 25

R W S E HP LF



Urinal Bottle -
Male urinal with handle
Graduated
Code: UMA1000

Dimensions: 190x100x80mm
Volume: 1000ml
Material: Polyethylene
Colour: Natural
Case Size: 40

R W S E HP LF



Slipper Pan -
Female urinal
Code: SPN

Dimensions: 335x245mm
Volume: 1500ml
Material: Polypropylene
Colour: Green
Case Size: 25

R W S E HP LF



Urinal Bottle -
Female urinal with handle
Code: UF1000

Dimensions: 170x100x80mm
Volume: 1000ml
Material: Polyethylene
Colour: Natural
Case Size: 30

R W S E HP LF



Urinal Pan -
Female Urinal
Code: PUF1000

Dimensions: 165x85x30mm
Volume: 1000ml
Material: Polyethylene
Colour: White
Case Size: 20

R W S E HP LF



Non Leak Valve -
For use with UMA1000 &
UML1000
Code: UMT

Dimensions: 50x90mm
Volume: N/A
Material: Polypropylene
Colour: White
Case Size: 25

R W S E HP LF



Urinal Bottle with UMT -Male urinal
bottle with non leak valve
Graduated
Code: USP1200

Dimensions: 180x100x85mm
Volume: 1000ml
Material: Polypropylene and
Polyethylene
Colour: N/A
Case Size: 10

R W S E HP LF



Urinal Bottle Holder -
For UBS1000, UML1000 and
UM1000
Code: UBH2002

Dimensions: N/A
Volume: 1200ml
Material: Plastic Coated Wire
Colour: White
Case Size: 10

R W S E HP LF

2 Single Use

W Washer Decontamination

S Suitable for Steam Sterilization

E Suitable for Ethylene Oxide Sterilization

HP Suitable for Hydrogen Peroxide Sterilization

LF Latex Free

R Reusable

W Washer Decontamination

S Suitable for Steam Sterilization

E Suitable for Ethylene Oxide Sterilization

HP Suitable for Hydrogen Peroxide Sterilization

LF Latex Free

Drinking Beakers, Cups, Lids and Denture Cups



Beaker - Feeder/drinking cup
Graduated

Code: BFC250

Dimensions: 62x75mm
Volume: 250ml
Material: Polypropylene
Colour: Natural
Case Size: 300

R W S E HP LF



Beaker - Feeder/drinking cup with
handles
Graduated

Code: BFH250

Dimensions: 62x80mm
Volume: 250ml
Material: Polypropylene
Colour: Natural
Case Size: 140

R W S E HP LF



Lid -
Narrow Spout for BFC250/
BFH250
Code: BNS

Dimensions: 4x4mm spout
Volume: N/A
Material: Polypropylene
Colour: Natural
Case Size: 300

R W S E HP LF



Lid -
Medium Spout for BFC250/
BFH250
Code: BMS

Dimensions: 7x4mm spout
Volume: N/A
Material: Polypropylene
Colour: Natural
Case Size: 300

R W S E HP LF



Lid -
Wide spout for BFC250
& BFH250
Code: BWS

Dimensions: 13x8mm spout
Volume: N/A
Material: Polypropylene
Colour: Natural
Case Size: 300

R W S E HP LF



Lid -
Rim feeder for beakers BFC250 &
BFH250
Code: BRF

Dimensions: N/A
Volume: N/A
Material: Polypropylene
Colour: Natural
Case Size: 200

R W S E HP LF



Drinking Cup -
With handles, spout lid and feeder
lid
Code: GPMC6250

Dimensions: 73x50mm
Volume: 250ml
Material: Polypropylene
Colour: Cream
Case Size: 90

R W S E HP LF



Teapot Style Beaker

Code: TFB

Dimensions: 88x55mm
Volume: 250ml
Material: Polypropylene
Colour: White
Case Size: 45

R W S E HP LF



Beaker Set - Feeder/drinking cup
with handle, spout lid and feeder lid
Graduated

Code: HDMS

Dimensions: 74x90mm
Volume: 300ml
Material: SAN
Colour: Clear
Case Size: 40

R W S E HP LF



Clear Glass Style Beaker

Code: TBC

Dimensions: 70x70mm
Volume: 280ml
Material: SAN
Colour: Clear
Case Size: 120

R W S E HP LF



Denture Cup -
With strainer and lid

Code: DNT

Dimensions: 90x80x60mm
Volume: N/A
Material: Polypropylene
Colour: Green
Case Size: 75

R W S E HP LF

R Reusable

W Washer Decontamination

S Suitable for Steam Sterilization

E Suitable for Ethylene Oxide Sterilization

HP Suitable for Hydrogen Peroxide Sterilization

LF Latex Free

Jugs and Jug Sets



Jug
Graduated

Code: JG500

Dimensions: N/A
Volume: 500ml
Material: Polypropylene
Colour: Natural
Case Size: 50

R W S E HP LF



Jug
Graduated

Code: JG1000

Dimensions: N/A
Volume: 1000ml
Material: Polypropylene
Colour: Natural
Case Size: 45

R W S E HP LF



Jug with lid
Graduated

Code: JGL1000

Dimensions: 120x125mm
Volume: 1000ml
Material: Polypropylene
Colour: Natural
Case Size: 40

R W S E HP LF



Jug Set - Includes tray, carafe with
lid and beaker
Graduated

Code: LST

Dimensions: N/A
Volume: 1000ml
Material: Polypropylene
Colour: N/A
Case Size: 20

R W S E HP LF

Medicine Measures



Medicine Measure -
Graduated

Code: MM60

Dimensions: 48x45mm
Volume: 60ml
Material: Polypropylene
Colour: Natural
Case Size: 1500

R W S E HP LF



Medicine Measure -
Graduated black lines

Code: MM60BL

Dimensions: 48x45mm
Volume: 60ml
Material: Polypropylene
Colour: Natural
Case Size: 1500

R W S E HP LF

How to order

To place an order, request
a quotation or for general
enquiries:

CALL us on
+44 (0) 1926 422 427
FAX
+44 (0) 1926 315 405

EMAIL
enquiries@sasco.co.uk

Patient Trays, Kidney Dishes and Vomit Bowls



Dispensing/Dressing Tray

Code: ITC2715

Dimensions: 270x150x30mm
 Volume: N/A
 Material: Polypropylene
 Colour: Green
 Case Size: 60

- R
- W
- S
- E
- HP
- LF



Kidney Dish
Graduated

Code: KDC250

Dimensions: 250x120x50mm
 Volume: 750ml
 Material: Polypropylene
 Colour: Green
 Case Size: 60

- R
- W
- S
- E
- HP
- LF



Vomit Bowl - With handle
Graduated

Code: VB300

Dimensions: 300x145x60mm
 Volume: 1500ml
 Material: Polypropylene
 Colour: Green
 Case Size: 25

- R
- W
- S
- E
- HP
- LF



- R Reusable
- W Washer Decontamination
- S Suitable for Steam Sterilization
- E Suitable for Ethylene Oxide Sterilization
- HP Suitable for Hydrogen Peroxide Sterilization
- LF Latex Free



Warwick House, Heathcote Way, Warwick, CV34 6TE
T: +00 44 (0) 1926 422 427
F: +00 44 (0) 1926 315 405
E: enquiries@sasco.co.uk

www.sasco.co.uk