



Endoscopy product catalog

Clinical and general practice supplies,
surgical instruments

WWW.BODY-PRODUCTS.DE



Protection and hygiene

Endoscopic accessories



Biopsy valves

Protection of patients
and nursing staff

Reliably leakproof

Securely attached

OLYMPUS, FUJIFILM, AND PENTAX

Biopsy valves

Protection of patients and nursing staff



Sterile

- Cap and seal the biopsy channel opening
- Prevents bodily fluids from leaking out
- Secure attachment to the biopsy channel opening
- Optimal leak tightness
- Easy to open and close
- Individually sterile packed



Color	Product description	Pieces per box	Product code
For OLYMPUS endoscopes			
Blue	For single use, sterile	200	BV-OX-02-EO
For PENTAX endoscopes			
Red	For single use, sterile	200	BV-PX-02-EO
For FUJIFILM endoscopes			
Green	For single use, sterile	200	BV-FX-02-EO

Biopsy valves



Protection of patients and nursing staff

Not sterile

- Cap and seal the biopsy channel opening
- Prevents bodily fluids from leaking out
- Secure attachment to the biopsy channel opening
- Optimal leak tightness
- Easy to open and close
- Individually packed



Color	Product description	Pieces per box	Product code
For OLYMPUS endoscopes			
Blue	For single use, not sterile	200	BV-OX-01
For PENTAX endoscopes			
Red	For single use, not sterile	200	BV-PX-01
For FUJIFILM endoscopes			
Green	For single use, not sterile	200	BV-FX-01

Endoscopic accessories



Valve sets

Valves and adapters in practical single-use sets

OLYMPUS, FUJIFILM, AND PENTAX

High safety

Constant performance

Valve sets for endoscopes

Valves and adapters in practical single-use sets



OLYMPUS

Biopsy valve, air/water valve, suction valve

- Suitable for OLYMPUS endoscopes from the 140 series
- Minimizes cross-contamination risks

Jet channel adapter

- With integrated check valve
- Individually sterile packed



Jet channel adapter, single-use



Biopsy valve



Jet channel adapter, 24h/use



Air/water valve



Air/water cleaning valve



Suction valve

Set contents	Pieces per box	Product code
Water jet channel adapter, single-use	100	ADA-02-JC-OUS
Water jet channel adapter, 24-hour use	100	ADA-JC-OUS
Biopsy valve, air/water valve, suction valve and water jet channel adapter	25	AWV-KIT-01-OLY
Biopsy valve, air/water valve, suction valve	25	AWV-KIT-02-OLY
Air/water valve, suction valve	25	AWV-KIT-03-OLY
Biopsy valve, air/water valve, suction valve, air/water cleaning valve	25	AWV-KIT-04-OLY
Biopsy valve, air/water valve, suction valve, air/water cleaning valve and water jet channel adapter	25	AWV-KIT-05-OLY
Air/water valve, suction valve, air/water cleaning valve, and water jet channel adapter	25	AWV-KIT-06-OLY
Air/water cleaning valve (not sterile) for OLYMPUS endoscopes	100	AWV-KIT-07-OLY

Valve sets for endoscopes

Valves and adapters in practical single-use sets



FUJIFILM

- Suitable for FUJIFILM endoscopes in the 500/600 and 700 series
- Minimizes cross-contamination risks

Jet channel adapter

- With integrated check valve
- Individually sterile packed



Jet channel adapter, single-use



Biopsy valve



Suction valve, FUJIFILM 500/600 series



Air/water cleaning valve, FUJIFILM 700 series



Air/water valve, FUJIFILM 500/600 series



Suction valve, FUJIFILM 700 series



Air/water valve, FUJIFILM 700 series

Set contents	Pieces per box	Product code
Suitable for FUJIFILM 700 series		
Water jet channel adapter, single-use	100	ADA-02-JC-FJ
Water jet channel adapter, 24-hour use	100	ADA-JC-FJ
Biopsy valve, air/water valve, and suction valve	25	AWV-KIT-01-FJ
Air/water valve, suction valve	25	AWV-KIT-02-FJ
Biopsy valve, air/water valve, suction valve, and water jet channel adapter	25	AWV-KIT-05-FJ
Biopsy valve, air/water valve, suction valve, and air/water cleaning valve	25	AWV-KIT-06-FJ
Air/water cleaning valve (not sterile)	100	AWV-KIT-07-FJ
Biopsy valve, air/water valve, suction valve, air/water cleaning valve, and water jet channel adapter	25	AWV-KIT-08-FJ
Suitable for FUJIFILM 500/600 series		
Biopsy valve, air/water valve, and suction valve	25	AWV-KIT-03-FJ
Air/water valve, suction valve	25	AWV-KIT-04-FJ
Biopsy valve, air/water valve, suction valve, and water jet channel adapter (from 580 series)	25	AWV-KIT-09-FJ

Valve sets for endoscopes

Valves and adapters in practical single-use sets



PENTAX

- Suitable for PENTAX endoscopes from 90/i10 series
- Minimizes cross-contamination risks

Jet channel adapter

- With integrated check valve
- Individually sterile packed



Water jet channel adapter, single-use



Water jet channel adapter, 24-hour use



Air/water valve



Biopsy valve



Suction valve

Set contents	Pieces per box	Product code
Suitable for PENTAX from 90/i10 series		
Water jet channel adapter, single-use	100	ADA-02-JC-PE
Water jet channel adapter, 24-hour use	100	ADA-JC-PE
Biopsy valve, air/water valve, suction valve and water jet channel adapter	25	AWV-KIT-01-PX
Biopsy valve, air/water valve, suction valve	25	AWV-KIT-02-PX
Air/water valve, suction valve	25	AWV-KIT-03-PX

Endoscopic accessories



Bite blocks

Hygiene and protection for patients and endoscopes

Ergonomically shaped

Soft edges

With and without strap

Bite blocks



Hygiene and protection for patients and endoscopes

- Soft edges increase safety and patient protection
- Optional holding straps increase patient comfort
- Single-use articles reduce the risk of cross-contamination

Adult

- Adapts to the facial contours thanks to ergonomic design



Color	Product description	Type of strap	Opening dimensions	Pieces per box	Product code
Silver	Without strap	-	25 x 22 mm	100	BB-05-30-30
Silver	With textile strap	Textile	25 x 22 mm	100	BB-05-30-30-T
Silver	With rubber strap	Rubber	25 x 22 mm	100	BB-05-30-30-R

Bite blocks



Hygiene and protection for patients and endoscopes

- Soft edges increase safety and patient protection
- Optional holding straps increase patient comfort
- Single-use articles reduce the risk of cross-contamination

Adult

- Adapts to the facial contours thanks to ergonomic design



Color	Product description	Type of strap	Opening dimensions	Pieces per box	Product code
Yellow	Without strap	-	27 x 23 mm	100	BB-01-85-50
Yellow	With white textile strap	Textile	27 x 23 mm	100	BB-01-85-50-T
Yellow	With rubber strap	Rubber	27 x 23 mm	100	BB-01-85-50-R

Pediatric

- Suitable for children and people with small mouths
- Adapts to the facial contours thanks to ergonomic design



Color	Product description	Type of strap	Opening dimensions	Pieces per box	Product code
Orange	Without strap	-	19 x 16 mm	100	BB-03-30-33
Orange	With textile strap	Textile	19 x 16 mm	100	BB-03-30-33-T
Orange	With rubber strap	Rubber	19 x 16 mm	100	BB-03-30-33-R



Cleaning

Endoscopic accessories



Cleaning brushes

Safe and gentle cleaning of
the endoscope channels

Rounded brush tips
Resistant to bending
Latex-free

Cleaning brushes




Safe and gentle cleaning of the endoscope channels

- Ball heads on the brush ends and rounded brush tips for gentle cleaning
- Brushes made of soft, durable nylon
- Ideal brush head size with optimal density
- Brush head diameter allows for maximum mechanical brushing of the endoscope channel
- Bend-resistant thanks to top-quality catheter rigidity

One-/two-sided cleaning brushes

- Efficient cleaning thanks to two brushes



Color	Catheter		Brush 1		Brush 2		Pieces per box	Product code
	∅	Length	∅	Length	∅	Length		
For endoscopes with working channels from 2.8 – 4.2 mm								
Yellow	1.7 mm	230 cm	6 mm	20 mm	-	-	100	SCBS-23-20-06-230
Blue	1.7 mm	230 cm	6 mm	20 mm	6 mm	20 mm	100	DCBS-23-20-06-230
Yellow	1.8 mm	180 cm	5 mm	20 mm	5 mm	20 mm	100	DCBS-23-20-05-180
Yellow	1.8 mm	230 cm	5 mm	20 mm	5 mm	20 mm	100	DCBS-23-20-05-230
Yellow	1.8 mm	240 cm	6 mm	20 mm	6 mm	20 mm	100	DCBS-23-20-06-240
For endoscopes with working channels from 2.8 – 5.2 mm								
Red	1.8 mm	230 cm	7 mm	20 mm	7 mm	20 mm	100	DCBS-23-20-07-230
Red	1.8 mm	180 cm	7 mm	20 mm	7 mm	20 mm	100	DCBS-23-20-07-180
Blue	1.5 mm	120 cm	6 mm	15 mm	6 mm	15 mm	100	DCBS-15-15-06-120
For endoscopes with working channels from 2.0 – 2.8 mm								
Green	1.5 mm	120 cm	3 mm	20 mm	3 mm	20 mm	100	DCBS-18-20-03-120
Green	1.5 mm	180 cm	3 mm	20 mm	3 mm	20 mm	100	DCBS-18-20-03-180
Green	1.5 mm	230 cm	3 mm	20 mm	3 mm	20 mm	100	DCBS-18-20-03-230

Cleaning brushes



Safe and gentle cleaning of the endoscope channels

- Ball heads on the brush ends and rounded brush tips for gentle cleaning
- Brushes made of soft, durable nylon
- Ideal brush head size with optimal density
- Brush head diameter allows for maximum mechanical brushing of the endoscope channel
- Bend-resistant thanks to top-quality catheter rigidity

Color	Catheter		Brush 1		Brush 2		Pieces per box	Product code
	∅	Length	∅	Length	∅	Length		
For endoscopes with working channels from 1.0 – 1.5 mm								
Blue	0.8 mm	120 cm	1.5 mm	10 mm	1.5 mm	10 mm	100	DCBS-08-10-15-120
Blue	0.8 mm	180 cm	1.5 mm	10 mm	1.5 mm	10 mm	100	DCBS-08-10-15-180
Blue	0.8 mm	230 cm	1.5 mm	10 mm	1.5 mm	10 mm	100	DCBS-08-10-15-230
For endoscopes with working channels from 1.2 – 1.8 mm								
Yellow	1.0 mm	120 cm	2 mm	15 mm	2 mm	15 mm	100	DCBS-15-15-02-120
Yellow	1.0 mm	230 cm	2 mm	15 mm	2 mm	15 mm	100	DCBS-15-15-02-230
For enteroscopes with working channels from von 2.8 – 4.2 mm								
Blue	1.7 mm	300 cm	6 mm	20 mm	6 mm	20 mm	100	DCBS-23-20-06-300

Cleaning brushes




Safe and gentle cleaning of the endoscope channels

- Ball heads on the brush ends and rounded brush tips for gentle cleaning
- Brushes made of soft, durable nylon
- Ideal brush head size with optimal density
- Brush head diameter allows for maximum mechanical brushing of the endoscope channel
- Bend-resistant thanks to top-quality catheter rigidity


Multipoint cleaning brush

- Combination brushes for working/suction channels
- With handle



Color	Catheter		Brush 1		Brush 2		Brush 3		Pieces per box	Product code
	∅	Length	∅	Length	∅	Length	∅	Length		
For endoscopes with working channels from 2.8 – 4.2 mm										
Blue	1.7 mm	230 cm	6 mm	20 mm	6 mm	20 mm	6 mm	20 mm	100	MCBS-23-20-06-230
Blue	1.7 mm	260 cm	6 mm	20 mm	6 mm	20 mm	6 mm	20 mm	100	MCBS-23-20-06-260

Endoscopic accessories



Combination cleaning brushes

Safe and gentle cleaning
of endoscope channels
and valves

High cleaning efficiency
Rounded brush tips
Latex-free

Combination cleaning brushes



Safe and gentle cleaning of endoscope channels and valves

- Ball heads on the brush ends and rounded brush tips for gentle cleaning
- Brushes made of soft, durable nylon
- Ideal brush head size with optimal density
- Brush head diameter allows for maximum mechanical brushing of the endoscope channel
- Bend-resistant thanks to top-quality catheter rigidity

With handle – Two-sided


- Combination brushes for the valve plane, valves, and working/suction channels
- With handle



Color	Catheter		Valve brush		Channel brush		Pieces per box	Product code
	∅	Length	∅	Length	∅	Length		
For endoscopes with working channels from 2.8 – 4.2 mm								
Blue	1.7 mm	230 cm	11 mm	40 mm	6 mm	20 mm	100	DVBS-23-40-11-20-06-230
Blue	1.7 mm	230 cm	11 mm	40 mm	5 mm	20 mm	100	DVBS-23-40-11-20-05-230
White	1.7 mm	230 cm	10 mm	40 mm	5 mm	15 mm	100	DVEC-23-40-10-15-05-230

Without handle – Two-sided

- Combination brushes for the valve plane, valves, and working/suction channels
- Without handle



Color	Catheter		Valve brush		Channel brush		Pieces per box	Product code
	∅	Length	∅	Length	∅	Length		
For endoscopes with working channels from 2.8 – 4.2 mm								
Yellow	1.7 mm	230 cm	10 mm	40 mm	6 mm	20 mm	100	DVBS-23-40-10-20-06-230

Combination cleaning brushes




Safe and gentle cleaning of endoscope channels and valves

- Ball heads on the brush ends and rounded brush tips for gentle cleaning
- Brushes made of soft, durable nylon
- Ideal brush head size and density
- Brush head diameter allows for maximum mechanical brushing of the endoscope channel
- Bend-resistant thanks to top-quality catheter rigidity

With sweeper

- Sweeper for greater cleaning efficiency



Color	Catheter		Brush		Sweeper		Pieces per box	Product code
	∅	Length	∅	Length	∅	Length		
For endoscopes with working channels from 2.8 – 4.2 mm								
Yellow	1.7 mm	230 cm	5 mm	20 mm	5 mm	300 mm	100	SCBS-23-20-5-300-5-230
Blue	1.7 mm	160 cm	5 mm	20 mm	5 mm	300 mm	100	SCBS-23-20-5-300-5-160
For endoscopes with working channels from 1.8 – 2.8 mm								
Orange	1.5 mm	160 cm	3 mm	20 mm	3 mm	300 mm	100	SCBS-18-20-3-300-3-160
Green	1.5 mm	120 cm	3 mm	20 mm	3 mm	300 mm	100	SCBS-18-20-3-300-3-120

Two-sided (Christmas tree shape)

- Combination brushes for working/ suction channels
- Conical brush heads for optimized cleaning




Color	Catheter		Channel brush 1		Channel brush 2		Pieces per box	Product code
	∅	Length	∅	Length	∅	Length		
For endoscopes with working channels from 2.8 – 4.2 mm								
Blue	1.7 mm	230 cm	4.5 – 7 mm	10 mm	3 – 6 mm	20 mm	100	DCBS-23-20-07/4.5-20-05-230

Combination cleaning brushes



Safe and gentle cleaning of endoscope channels and valves

- Ball heads on the brush ends and rounded brush tips for gentle cleaning
- Brushes made of soft, durable nylon
- Ideal brush head size and density
- Brush head diameter allows for maximum mechanical brushing of the endoscope channel
- Bend-resistant thanks to top-quality catheter rigidity



Two-sided (Christmas tree shape) with sweeper

- Combination brushes for working/ suction channels
- Conical brush heads for optimized cleaning
- Sweeper for greater cleaning efficiency

Color	Catheter		Channel brush 1		Sweeper		Channel brush 2		Pieces per box	Product code
	∅	Length	∅	Length	∅	Length	∅	Length		
For endoscopes with working channels from 2.8 – 4.2 mm										
Blue	1.7 mm	260 cm	4.5 – 7 mm	10 mm	5 mm	100 mm	3 – 6 mm	20 mm	100	DVBS-23-20-07/4.5-20-05-260



Valve/channel cleaning brushes

Easy cleaning of working channels and valves

High cleaning efficiency
Rounded brush tips
Maximum brush density


Valve/channel cleaning brushes



Easy cleaning of working channels and valves


- Rounded brush tips for gentle cleaning
- Brushes made of soft, durable nylon
- Rounded handle

Two-sided valve/ channel cleaning brush




Color	Length	Valve brush		Channel brush		Pieces per box	Product code
		∅	Length	∅	Length		
Blue	15 cm	11 mm	35 mm	5 mm	20 mm	200	DVS-30-35-10-20-05-15

Valve cleaning brush with handle



Color	Length	Valve brush		Pieces per box	Product code
		∅	Length		
Green	9.5 cm	11 mm	35 mm	200	SVS-30-35-11

Cleaning brush for Albarrán lever on duoendoscopes



Color	Length	Brush		Pieces per box	Product code
		∅	Length		
Red	8 cm	2 mm	15 mm	200	SES-30-2-15

Endoscopic accessories



Cleaning brush sets

Practical sets for efficient cleaning of working channels and valves

High cleaning efficiency
Rounded brush tips
Ergonomic handle

Cleaning brush sets



Practical sets for efficient cleaning of working channels and valves


Two-sided cleaning brush and two-sided valve cleaning brush

Two-sided cleaning brush

- Efficient cleaning thanks to two brushes

Two-sided valve and channel cleaning brush

- For the valve plane, valves, and working/suction channels
- Rounded brush tips for gentle cleaning
- Brushes made of soft, durable nylon



Two-sided cleaning brushes							Two-sided valve cleaning brushes				Pieces per box	Product code
Color	Catheter		Brush 1 and 2				Valve brush		Channel brush			
	∅	Length	∅	Length	∅	Length	∅	Length	∅	Length		
For endoscopes with working channels from 2.8 – 4.2 mm												
Yellow	1.8 mm	180 cm	5 mm	20 mm	5 mm	20 mm	11 mm	35 mm	5 mm	20 mm	100	DCBS-SET-180-20-05-35-10
Yellow	1.8 mm	230 cm	5 mm	20 mm	5 mm	20 mm	11 mm	35 mm	5 mm	20 mm	100	DCBS-SET-20-05-35-10
Blue	1.7 mm	230 cm	4.5 – 7 mm	10 mm	6 – 3 mm	20 mm	11 mm	35 mm	5 mm	20 mm	100	DBC-SET-23-20-07/4.5-20-05-230
Pink	1.65 mm	240 cm	5 mm	20 mm	5 mm	20 mm	11 mm	35 mm	5 mm	20 mm	100	DCBS-SET-23-20-05-240
Yellow	1.7 mm	230 cm	5 mm	20 mm	-	-	11 mm	35 mm	5 mm	20 mm	100	SCBS-SET-20-05-300-5-230

Cleaning brush sets



Practical sets for efficient cleaning of working channels and valves


Multipoint cleaning brush and two-sided valve cleaning brush

Multipoint cleaning brush

- Efficient cleaning thanks to three brushes

Two-sided valve and channel cleaning brush

- For the valve plane, valves, and working/suction channels
- Rounded brush tips for gentle cleaning
- Brushes made of soft, durable nylon



Multipoint cleaning brush					Two-sided valve cleaning brush				Pieces per box	Product code
Color	Catheter		Brush 1/2/3		Valve brush		Channel brush			
	∅	Length	∅	Length	∅	Length	∅	Length		
For endoscopes with working channels from 2.8 – 4.2 mm										
Blue	1.7 mm	230 cm	6 mm	20 mm	11 mm	35 mm	5 mm	20 mm	100	MCBS-SET-23-20-06-230
Blue	1.7 mm	260 cm	6 mm	20 mm	11 mm	35 mm	5 mm	20 mm	100	MCBS-SET-23-20-06-260

Endoscopic accessories



Single-use EcoFlushCap

Efficient lens cleaning

Can be used for 24 hours

With and without CO₂ connection

OLYMPUS/FUJIFILM/PENTAX

EcoFlushCap HybriX single-use system



Higher patient protection combined with lower costs

The innovative, 24-hour single-use **EcoFlushCap HybriX** system supports the high demands on patient safety and significantly reduces the risk of cross-contamination. Suitable for use with any sterile water bottle system, (500 ml and 1000 ml), the HybriX system ensures the effective flushing of the working channel/jet channel and reliably cleans the endoscope lens. The intelligent Twin Tubing system makes the HybriX system ideal for use with all common flushing pumps. Thanks to the optional single-use connectors for various CO₂ insufflators, the **EcoFlushCap HybriX system** also supports CO₂ insufflation.

- Jet channel connector including check valve
- For OLYMPUS/FUJIFILM/PENTAX endoscopes



Lens cleaning

Ensures a clear view and optimal visibility



External flushing

Additional flushing option



CO₂ insufflation

Reduces pressure and tension in patients

Twin Tubing for compatibility with all common flushing pump units



Pieces per box	Product code			
	For OLYMPUS	For FUJIFILM 500/600 series	For FUJIFILM 700 series	For PENTAX
With CO ₂ connection tube				
10	WB-HYB-OL4	WB-HYB-FJ4	WB-HYB-FJ5	WB-HYB-PX5

Single-use EcoFlushCap

Efficient lens cleaning



- Can be used for 24 h
- With CO₂ tube system
- Compatible with all common sterile water bottles (500 ml and 1000 ml)



EcoFlushCap for sterile water bottles with CO₂ connection

Pieces per box	Product code			
	For OLYMPUS	For FUJIFILM 500/600 series	For FUJIFILM 700 series	For PENTAX
With CO ₂ connection				
10	WB-OL3-WC2	WB-FJ3-WC2	WB-FJ6-WC2	WB-PX6-WC2

EcoFlushCap for sterile water bottles with 110 cm CO₂ connection tube

Pieces per box	Product code			
	For OLYMPUS	For FUJIFILM 500/600 series	For FUJIFILM 700 series	For PENTAX
With 110 cm CO ₂ connection tube				
15	WB-OL4-WC2	WB-FJ4-WC2	WB-FJ5-WC2	WB-PX5-WC2

Single-use EcoFlushCap

Efficient lens cleaning



- Can be used for 24 h
- Compatible with all common sterile water bottles (500 ml and 1000 ml)



EcoFlushCap for sterile water bottles without CO₂ connection



Pieces per box	Product code		
	For OLYMPUS	For FUJIFILM 500/600 series	For PENTAX
Without CO ₂ connection			
10	WB-OL2-NC2	WB-FJ2-NC2	WB-PX4-WC2

Twin EcoFlushCap

- Compatible with all common sterile water bottles (500 ml and 1000 ml)
- Can be used for 24 h
- Suitable for FUJIFILM endoscopes in the 500/600/700 series



FUJIFILM 700 series



FUJIFILM 500/600 series



Pieces per box	For FUJIFILM 500/600/700 series	Product code
15	with CO ₂ connection	WB-FJ7-WC2

Connection, suction, and flushing tubes



Suction and flushing tube with cap for flushing pumps

- Can be used for 24 h
- Compatible with all common sterile water bottles (500 ml and 1000 ml)
- Jet channel connector including check valve



Product description	Pieces per box	Product code
For all common flushing pumps	15	IT-OUS-210

Suction and flushing tube with 2 male Luer-Lock connections and check valve

- Can be used for 24 h



Product description	Pieces per box	Product code
For various flushing pumps	15	IT-ETS-210

Connection, suction, and flushing tubes



Flushing tube for Olympus endoscopes

- Individually packed



Product description	Pieces per box	Product code
Manual Cleaning tube adapter	50	OLY-CT-001

Adapters



Jet channel adapter

- With integrated check valve
- Individually sterile packed
- Can be used for 24 h



For OLYMPUS



For FUJIFILM



For PENTAX from the 90 series

Product description	Pieces per box	Product code
Jet channel adapter for OLYMPUS	100	ADA-JC-OUS
Jet channel adapter for FUJIFILM	100	ADA-JC-FJ
Jet channel adapter for PENTAX from the 90 series	100	ADA-JC-PE

Adapter CO₂ insufflator



For FUJIFILM



For OLYMPUS

Product description	Pieces per box	Product code
Single-use adapter for OLYMPUS CO ₂ insufflator (UCR)	100	ADA-CO2-OUS
Single-use adapter for FUJIFILM CO ₂ insufflator (GW-100)	100	ADA-CO2-FJ

Jet channel adapter

- Individually sterile packed
- Single-use
- With integrated check valve



For PENTAX from the 90 series



For FUJIFILM



For OLYMPUS

Product description	Pieces per box	Product code
Jet channel adapter for OLYMPUS	100	ADA-02-JC-OUS
Jet channel adapter for FUJIFILM	100	ADA-02-JC-FJ
Jet channel adapter for PENTAX from the 90 series	100	ADA-02-JC-PE

Bottle holder for sterile water bottles



Plastic

Compatible with all common sterile water bottles (500 ml and 1000 ml)

- 1 hook for rail mounting
- With bracket for rear attachment
- 118 x 113 x 113 mm (W x H x D)



Product description	Pieces per box	Product code
Bottle holder	1	BT-ABS-001



Protection and cleaning



Endoscope distal protection meshes

Protection at the distal end of the endoscope

Optimal protection
Reliable support
Flexible dimensions


Endoscope distal protection meshes

Protection at the distal end of the endoscope



Distal protection mesh

- Flexible net to protect the distal end of endoscopes during transport and storage
- 2 colors for marking different sizes
- Individually sterile packed



Color	Suitable for endoscopes with a diameter of	Length	Pieces per box	Product code
Blue	7 – 15 mm	50 mm	100	DTP-01
Red	2 – 10 mm	50 mm	100	DTP-02

Endoscope distal protection sponges

Protection and optimal support at the distal end of the endoscope



- Optimal protection of the distal end of endoscopes during transport and storage
- Reliable support for endoscopes during storage
- 2 colors for marking different sizes
- Individually sterile packed



Color	Diameter	Length	Outer diameter	Inner diameter	Pieces per box	Product code
Blue	> 5 mm	140 mm	37 mm	9 – 24 mm	50	EP-14-01-EO
Red	< 5 mm	140 mm	37 mm	4 – 24 mm	50	EP-05-01-EO

3-in-1 distal protection and cleaning sponge


Protection, support, and cleaning at the distal end of the endoscope



- 3-in-1 product
- Optimal protection of the distal end of endoscopes during transport and storage
- Reliable support for endoscopes during storage
- Reliable cleaning of endoscopes
- Individually sterile packed



Color	Diameter	Length	Outer diameter	Inner diameter	Pieces per box	Product code
White	>/< 5 mm	126 mm	46 mm	-	50	EP-05-14-CL-EO



Cleaning

Accessories



Single-use colonoscopy dignity shorts

Maintain dignity during the colonoscopy

Comfortable to wear
Lightweight non-woven fabric
Examination slit in the back

Medi-Inn® colonoscopy dignity shorts



Maintain patient dignity during the colonoscopy

The single-use Medi-Inn unisex dignity shorts are the cost-effective and hygienic alternative to reusable shorts. The dark blue shorts in Bermuda design are made of lightweight PP non-woven fabric (polypropylene) and are available in 2 sizes.

Reliable protection

The highly durable, lightweight PP non-woven fabric protects against moisture and dirt

Comfortable to wear

The elasticated waist ensures a perfect and secure fit

Easy access

The rear, slightly overlapping examination slit provides easy access for rectal examinations

- Unisex examination shorts made of lightweight PP non-woven fabric (polypropylene)
- Comfortable to wear
- Extremely durable and moisture-repellant
- Elastic rubber band
- Examination slit in the back

Size	Color	Packaging unit	Product code
XL – 120 x 65 cm	Dark blue	100	DS-L-ECO
XXL – 130 x 65 cm	Dark blue	100	DS-L-ECO-XXL



Your contact:

BODY PRODUCTS RELAX
Pharma und Kosmetik GmbH

Alfred-Nobel-Straße 1-3
50226 Frechen

Phone: +49 2234 91165-0
Fax: +49 2234 91165-23
Email: info@body-products.de

WWW.BODY-PRODUCTS.DE

REV. 00 | 09/22