



# **TRINITY** MEDICAL

PATIENT CENTRIC TECHNOLOGIES

## **Decontamination & IPC Product Catalogue**

**2025-2026**

[www.trinitymed.ie](http://www.trinitymed.ie)


# Your partner for....



CSSD & Surgical



Endo & SIUP



IPC



TECH SUPPORT



**CSSD**

**Non Linting Wipes**



Code	Descript.	Pack	Notes
HNL400C	Non Linting Wipes 30x38cm	400sheets/bx	MOQ 4 boxes

**Long Cuff Nitrile Glove**



Code	Descript.	Pack	Notes
LCG-S, LCG-M, LCG-L, LCG-XL	40cm Long Cuff Glove	500/case	-

**Water Repellent Gown**




Code	Descript.	Pack	Notes
810008	Water Repellent Gown	200/bx	-

**Disposable Face Shields**



Code	Descript.	Pack	Notes
PS2000	Disposable Face Shield	200/bx	-





### Trifecta™ Brush – 3-in-1 Instrument Cleaning Tool

The Trifecta™ Brush is a versatile cleaning solution designed for effective instrument maintenance. Available in nylon or a combination of nylon and stainless steel bristles, it features a precision tip for hard-to-reach areas and firm bristles for tough residues. Ideal for removing buildup while protecting delicate surfaces, it ensures thorough and safe cleaning.



Code	Descript.	Pack	Notes
28-23012	Triple Nylon Brush	3/pk	For a bag of 50, add "-50" following the item code.
29-23012	Triple Nylon Brush with Stainless Steel	3/pk	For a bag of 50, add "-50" following the item code.

**Reamer Cleaning Brush**

Code	Descript.	Pack	Notes
50-ACETAB	Acetabular Reamer Cleaning Brush	3/pk	Reusable

**General Cleaning Brush**

Code	Descript.	Pack	Notes
28-72550-50	Cleaning Brush Nylon	50/pk	-

**General Cleaning Brush**

Code	Descript.	Pack	Notes
29-72550-50	Cleaning Brush Steel Bristle	50/pk	-

**Precision Cleaning Brush**

Code	Descript.	Pack	Notes
50-40800	Precision Cleaning Brush	3/pk	Autoclavable

**Stainless Steel ARC™ Micro Brushes**

Code	Descript.	Pack	Notes
35-03015	Micro Brush 1.5mm	6/pk	30cm Long
35-04515	Micro Brush 1.5mm	6/pk	45cm Long

**General Cleaning Brush - Brass**

Code	Descript.	Pack	Notes
29-70050	Cleaning Brush Brass	50/pk	n/a



**Flat Cleaning Brush**



Code	Descript.	Pack	Notes
50-40917	Flat Brush Blue	3/pk	Extra Rigid
50-40918	Flat Brush Green	3/pk	Rigid
50-40919	Flat Brush Yellow	3/pk	Flexible
50-40920	Flat Brush Red	3/pk	Extra Flexible
50-40921	Single-Use Flat Brush Orange	20/pk	Flexible
50-40922	Single-Use Flat Brush Light Blue	20/pk	Rigid

**V- Brush**



Code	Descript.	Pack	Notes
36-20005	Single Use V-Brush	50/pk	Reusable also available

**Reusable Orthopedic Rasps Cleaning Brush**



Code	Descript.	Pack	Notes
50-RASP	Cleaning Brush for Rasps	3/pk	n/a

**Twisted Wire Cleaning Brush**



Scan for full brochure



Code	Descript.	Pack	Notes
n/a	Various options available	3/pk & 50/pk	Contact Trinity Med for more info.

**Arc Plastic Handle Brushes**



Scan for full brochure




Code	Descript.	Pack	Notes
n/a	Plastic Handle Cleaning Brush - Various Sizes - Single & Double	3/pk & 50pk	Contact Trinity Med for more info.

**Endo Valve Cleaning Brushes**



Scan for full brochure



Code	Descript.	Pack	Notes
Scan or click on QR code for full brochure	Valve Cleaning Brush - Various Sizes	3/pk	Contact Trinity Med for more info.

**Bottle & Drain Brushes**



Code	Descript.	Pack	Notes
50-10295	Reusable Jar / Bottle Brush	3/pk	Autoclavable
50-10224	Sterilizer Drain Cleaning Brush	3/pk	Autoclavable

**Channel Brushes - Individual & Sets**



Scan for full brochure

Code	Descript.	Pack	Notes
Scan or click on QR code for full brochure	Channel Cleaning Brushes - Various Sizes	various	Contact Trinity Med for more info.

**Double Ended Brushes**



Code	Descript.	Pack	Notes
28-69394 / 28-69394-50	Double Ended Brush - Nylon	3/pk or 50/pk	n/a
29-69394 / 29-69394-50	Double Ended Brush - Steel	3/pk or 50/pk	n/a

**Pipe Cleaners**




Code	Descript.	Pack	Notes
SOFTEM004	SoftClean Pipe-Cleaner 9mm x 30cm	Pk/50	More sizes available
SOFTEM033	SoftClean Pipe-Cleaner 3mm x 30cm	Pk/50	More sizes available



Padlock Security Tags - Complete Range




Scan for full brochure



Code	Descript.	Pack	Notes
Various	Full Range Available	-	Contact Trinity Med for more info.

**Protective Corner Pockets**



Scan for full brochure

Code	Descript.	Pack	Notes
Various	Variety of disposable and reusable protective corner pockets	n/a	Foam and Silicon variations available.

**Instrument Tip Protectors - Double**



Scan for full brochure

Code	Descript.	Pack	Notes
Various	Double Tip Protectors	50/pk	Contact Trinity Medical for information

**Instrument Tip Protectors - Single**



Code	Descript.	Pack	Notes
Various	Single Tip Protectors	100/pk	Contact Trinity Med for more info.

**Misc Instrument Tip Protectors**



Code	Descript.	Pack	Notes
Various	Miscellaneous Tip Protectors	-	Contact Trinity Med for more info.



**Scrub Pad for Instruments**

Code	Descript.	Pack	Notes
60-46003	Scrub Pad for Instruments	3/pk & 50pk	n/a

**Silicon Pin Mat**

Code	Descript.	Pack	Notes
13-2001	Silicon Pin Mat	1pk & 5pk	n/a

**Anti-Fog Solution**

Code	Descript.	Pack	Notes
G.O.L.F.F.	Anti-Fog Solution & Adhesive Sponge set	20/bx	Contact Trinity Med for more info.

**Clamp Jaw Covers**

Code	Descript.	Pack	Notes
Various	Clamp Jam Covers	Various	Sterile & Non Sterile available. Contact Trinity Medical for more info.



**Pin Guards**



Code	Descript.	Pack	Notes
Various	Pin Guards for K Wires	n/a	Contact Trinity Med for more info.

**Anti-Fatigue Mats**



Code	Descript.	Pack	Notes
800-2401	Anti-Fatigue Mat	1 / Box	18" X 24" / 46 cm x 61 cm
800-2405	Anti-Fatigue Mat	5 / Box	18" X 24" / 46 cm x 61 cm
800-3201	Anti-Fatigue Mat	1 / Box	20" X 32" / 51 cm x 81 cm
800-3203	Anti-Fatigue Mat	3 / Box	20" X 32" / 51 cm x 81 cm
800-4801	Anti-Fatigue Mat	1 / Box	20" X 48" / 51 cm x 122 cm
800-4802	Anti-Fatigue Mat	2 / Box	20" X 48" / 51 cm x 122 cm
800-7201	Anti-Fatigue Mat	1 / Box	30" X 72" / 76 cm x 183 cm

**Cautery Tip Cleaner**



Code	Descript.	Pack	Notes
60-10000	Cautery Tip Cleaner Non-Sterile	2,200 /Bx	4.5 cm x 4.5 cm



## Instrument Organisation &amp; Protection - Foam



Foam pouches prevent costly damage to instrument tips caused by shifting during handling and transport. Pouches neither retain nor absorb moisture, are 100% porous allowing for full penetration of sterilant gases, and validated for use in steam and EtO sterilisers.

Code	Descript.	Pack	Notes
FP-1001	FOAM POUCH 3.75 X 10.06	CS-50	MOQs apply
FP-1002	FOAM TIP PROTECTOR	CS-50	MOQs apply
FP-1003	FOAM PINHOLE HOLDER	CS-50	MOQs apply
FP-1004	FOAM CORNER PROTECTORS (SET OF 4)	CS-288	MOQs apply
FP-1005	12 INSTRUMENT FOAM PROTECTOR	CS-50	MOQs apply
FP-4600	FOAM PROTECTION PRODUCTS - ROLL 12W X 108L	CS-3	MOQs apply
FP-4800	FOAM PROTECTION PRODUCTS - PRECUT SHEETS 12W X 10L	CS-50	MOQs apply
FP-4850	FOAM PROTECTION PRODUCTS - PRECUT SHEETS 12W X 16L	CS-25	MOQs apply
FP-4880	FOAM PROTECTION PRODUCTS - PRECUT SHEETS 12W X 20L	CS-25	MOQs apply
FP-5200	FOAM PROTECTION PRODUCTS - MULTI-POUCH WITH BELTS - 11W X 15L	CS-10	MOQs apply
FP-5500	FOAM PROTECTION PRODUCTS - UNIVERSAL MULTI-POUCH - 13W X 17L	CS-10	MOQs apply
FP-6200	FOAM PROTECTION PRODUCTS - ENDOSCOPIC POUCH - 2W X 20L	CS-50	MOQs apply
FP-6210	FOAM PROTECTION PRODUCTS - LARGE POUCH - 4W X 7L	CS-50	MOQs apply
FP-6250	FOAM PROTECTION PRODUCTS - SMALL SINGLE POUCH - 2.1W X 7.5L	CS-50	MOQs apply
FP-6270	CAMERA/FIBER-OPTIC CABLE 9.5W X 10L	CS-10	MOQs apply

**Magnifying Lamp**



Code	Descript.	Pack	Notes
D25110-UK	5 LED Magnifier Lamp	each	Clamp fixing or table stand available.

**LED 4X MAGNIFIER**



Code	Descript.	Pack	Notes
MAG-002	LED 4x Magnifier	each	n/a

**Paper Style Tip Protector**



Code	Descript.	Pack	Notes
HE1-0006-W	Paper Style Tip Protector	100/bx	168 x 90mm

**Instrument Stringers**



Code	Descript.	Pack	Notes
Various	Multiple Options Available	-	Contact Trinity Medical for more info.

**Autoclave Cleaning Brush Kit**



Code	Descript.	Pack	Notes
50-CLAVESUP	Autoclave SUPERIOR Brush Cleaning Kit	each	n/a

**Spray Gun - Air & Water**



Code	Descript.	Pack	Notes
Various	Air / Water Spray Gun - Multiple options available	each	Contact Trinity Medical for more info.

**Identification Labels & Label Tags**



Code	Descript.	Pack	Notes
Various	Various colours available, including variants with text such as 'Bio-hazard' and 'Repair'.	-	Contact Trinity Medical for more info.

**Tape-B-Gone™ Tape Remover**



Code	Descript.	Pack	Notes
T-32RMV	Tape-Be-Gone Tape Remover	4 x 237ml	n/a



**Electrosurgical Insulation Defect Detector**



Electrosurgical insulation defect detector is a compact, hand held, battery operated unit that is used to detect and locate defects such as pinholes, cracks and bare spots in protective coatings applied over conductive instrument's surfaces.

Code	Descript.	Pack	Notes
n/a	Insulation Tester	each	Contact Trinity Medical for more info.

**Flexible Inspection Scope for Channels of Hollow Instruments & Flexible Endoscopes**



Visually inspect endoscopes and small lumens with remarkable clarity

- Ergonomic design
- Compact control box with flexible inspection scope in diameter of 2.1mm.
- This scope has a working length of 160cm and attaches to a control box which allows for light level adjustments.
- The scope is designed to inspect internal channels of 2.1mm in diameter or larger.

Code	Descript.	Pack	Notes
n/a	Sherlock Flexible Inspection Scope	each	Contact Trinity Medical for more info.

**Magic Touch Magnifier**



Code	Descript.	Pack	Notes
RNDMAG-1	Magic Touch Magnifier	each	Ø 9.53 x 4.57h

**Residual Protein Testing - Qualitative Result**



Code	Descript.	Pack	Notes
26933300	30 Tests with short swabs.	each	High sensitivity, with control test.
26933500	Endo - 30 Tests with long swabs 1.9mm diameter	each	High sensitivity, with control test.
26933400	Endo - 30 Tests with long swabs 2.8mm diameter	each	High sensitivity, with control test.
26933600	Endo - 30 Tests with long swabs 3.7mm diameter	each	High sensitivity, with control test.

**Cleaning Efficacy Indicators - CEIs**



Code	Descript.	Pack	Notes
RI-52002	Cleaning Efficacy Indicators CEIs, Bx/50	50/bx	Available for lumens also.

**Dispenser for Autoclave Tape**



Code	Descript.	Pack	Notes
TAPEDISDOUB	50mm Heavy Duty Bench Tape Dispenser	each	n/a
TAPEDISTRIP	75mm Heavy Duty Bench Tape Dispenser	each	n/a

**Terragene Ultrasonic Cavitation Test**



Code	Descript.	Pack	Notes
CDWU	Ultrasonic Cavitation Test	30/BX	n/a

**Terragene SPLAT - Washing Indicator**



Scan for to watch video



Code	Descript.	Pack	Notes
CDWA4 - High	Splat - Cleaning Efficacy Indicators	200/pk	Compatible with Trazanto
CDWA3 - MED	Splat - Cleaning Efficacy Indicators	200/PK	Compatible with Trazanto

**Terragene Bowie & Dick Test Pack**



Scan for to watch video



Code	Descript.	Pack	Notes
BD125X-2	Bowie & Dick Test Pack	20/bx	High Contrast Colour Change. Compatible with Trazanto.

**Terragene Chemical Integrator for Steam**



Code	Descript.	Pack	Notes
IT26-C	Chemical Integrator for Steam	250/bx	Compatible with Trazanto.
IT26-1YS	Chemical Integrator for Steam - Colour Change	200/bx	Compatible with Trazanto.



**Push Dye Test for Heat Seal Reels & Pouches**



Code	Descript.	Pack	Notes
44100	For Paper Packing	30/bx	n/a
44110	For Tyvek Packing	30/bx	n/a

**Seal Sure**



Code	Descript.	Pack	Notes
PS2100	Routine monitor For Sealers	250/pk	n/a

**Terragene High Sensitivity Pen for Quantifiable Protein Detection**



The HSE recommends the use of technologies to allow for the quantification of proteins.



Code	Descript.	Pack	Notes
PRO1MICRO	Swab & Incubate for Quantifiable Protein Detection	100/bx	Use with MiniPro Incubator

**Terragene Trazanto Analyser**



Code	Descript.	Pack	Notes
TRAZANTO	Analyser for Various Indicators	each	See video for more info



**Photon - Biological Indicator for Steam - Instant Readout 7 Seconds**



Code	Descript.	Pack	Notes
BT225	7 Second Instant Readout for Steam Sterilisation	50/bx	Geobacillus stearothermophilus

**Hyper - Biological Indicator for Hydrogen Peroxide - 5min**



Code	Descript.	Pack	Notes
BT98	5min Biological Indicator for Hydrogen Peroxide	50/BX	Geobacillus stearothermophilus

**Hyper & Photon Auto-Readers for Rapid Read-out of Biological Indicators**



Code	Descript.	Pack	Notes
BPH	Photon Autoreader	each	n/a
BHY	Hyper Autoreader	each	n/a

**Lumenia Six Flow Cleaning Challenge - For AERs**



Scan for to watch video



Code	Descript.	Pack	Notes
LSF1	6 Channel AER Challenge Device	each	Compatible with Trazanto



**EnziMed Pre-Cleaner**

Code	Descript.	Pack	Notes
OL20304	Multi-enzymatic foam spray for surgical instruments	6/bx	n/a

**EnziMed 4 - Standard Enzymatic Detergent for Manual Cleaning**

Code	Descript.	Pack	Notes
OL20506	Low Foaming Enzymatic Detergent 5ltr.	4/bx	For manual and ultrasonic cleaning applications

**EnziMed 5 - Enzymatic Detergent for Instrument Washer Disinfectors**


Code	Descript.	Pack	Notes
OL20206	EnziMed Instrument Washer 5ltr.	4/bx	Non foaming for automated washer disinfectors.

**OneLife - Detect 2**

Code	Descript.	Pack	Notes
OL22329	Detect 2 Reveals Organic Contamination Real-time 950ml	6/bx	Contact Trinity Med for more info




**Specialist Baskets for Rigid Scopes & Robotic Surgery (DaVinci)**



Code	Descript.	Pack	Notes
n/a	Various models & sizes available	each	Contact Trinity Med for more info.

**Instrument Trays DIN**



Code	Descript.	Pack	Notes
750046	1/2 DIN 240x250x70mm	each	n/a
750048	Full DIN 480x250x70mm	each	n/a
750049	Full DIN Deep 480x250x100mm	each	n/a
7500911	Lid for Full DIN	each	n/a
7500912	Lid for 1/2 DIN	each	n/a

**Fine Mesh Baskets**



Download the complete Hupfer brochure

Code	Descript.	Pack	Notes
7501282	Fine Mesh Basket with Lid & Clasp	each	80 x 40 x 30mm
7501220	Fine Mesh Basket with Lid & Clasp	each	80 x 49 x 42mm
7501132	Fine Mesh Basket with Lid & Clasp	each	80 x 80 x 40mm
7501204	Fine Mesh Basket with Lid & Clasp	each	200 x 100 x 50mm
7501133	Fine Mesh Basket with Lid & Clasp	each	220 x 140 x 50mm
7501135	Fine Mesh Basket with Lid & Clasp	each	274 x 174 x 36mm
7501161	Fine Mesh Basket with Lid & Clasp	each	210 x 150 x 40mm
7501162	Fine Mesh Basket with Lid & Clasp	each	310 x 210 x 40mm



**SMS Heavy Grade Wrap - Interleaved**



Code	Descript.	Pack	Notes
T4007575	Blue / Green SMS Interleaved Wrap for Steam and VH2O2	126 sheets per box	75x75cm
T4009090	Blue / Green SMS Interleaved Wrap for Steam and VH2O2	100 sheets per box	90X90cm
T500100100	Blue / Green SMS Interleaved Wrap for Steam and VH2O2	100 sheets per box	100x100cm
T500120120	Blue / Green SMS Interleaved Wrap for Steam and VH2O2	76 sheets per box	120x120cm

**Reels & Pouches**



Code	Descript.	Pack	Notes
Various	Reels and Pouches	n/a	Contact Trinity Med for more info.

**Tyvek Reels & Pouches**



Code	Descript.	Pack	Notes
Various	Flat Reels and Pouches for low temperature sterilisation VH202	n/a	Contact Trinity Med for more info.



**Gandus Rotary Heat Sealer - Pura**



Code	Descript.	Pack	Notes
PURA	Validatable Medical Rotary Heat Sealer. Touch Screen.	each	Available with front support plate. Printer included.

**Reel Holder with Cutting Blade**



Code	Descript.	Pack	Notes
n/a	Stainless steel Reel Holder, fitted with a cutting blade.	each	Available in various widths.

**Ultrasonic Baths**



Code	Descript.	Pack	Notes
n/a	Various models & sizes available	each	Contact Trinity Med for more info.

**Wrap Holder**



Code	Descript.	Pack	Notes
7500371	Holder for Sterile Barrier Sheets	each	Available with lower retaining grid.



**Packing Stations**



Code	Descript.	Pack	Notes
n/a	Customisable Height Adjustable Packing Stations	each	Contact Trinity Medical for more info.

**Ensure Touch ATP Hygiene Monitoring**



Code	Descript.	Pack	Notes
ETOUCH	Ensure Touch Cleaning Verification System	each	Contact Trinity Medical for more info.
US2020	Hygiena UltraSnap™ Surface Swap	pk/100	n/a
ASY0355	Endoswab 3mm	pk/100	n/a
ASY0469	Endoswab 2mm	pk/100	n/a
ASY0468	Endoswab 1mm	pk/100	n/a

**Sychem Control Plus Wipes**



Code	Descript.	Pack	Notes
SYC-SWC-W200	Control Plus Surface Wipes, Tub, Box of 6	Box of 6 Tubs	MOQ 3 Boxes

**Sychem Control Plus - HLD Spray RTU**



Code	Descript.	Pack	Notes
SYC-SWC-RTU	Control Plus Ready to Use Spray, High Level Disinfection	750mm x 12	-
SYC-020152-C-PLUS	Sychem Control Concentrate PLUS	5ltr x 2	-

**Warwick SASco - Moulded Medical Plastics**



Download the complete Warwick Sasco brochure

Code	Descript.	Pack	Notes
Various	Comprehensive range of medical plastics.	-	Contact Trinity Medical for more info.





**Ultrasonic Irrigation Wash Station - Ultrasonic MED**



Code	Descript.	Pack	Notes
n/a	Twin Bowl High Efficiency Ultrasonic Wash Station. Height adjustable, with in-built traceability software. Compatible with DaVinci Handpieces.	each	Contact Trinity Medical for more info.

**Custom Stainless Steel Fabrication**

Trinity Medical offers a design and custom fabrication service to help customers bring their ideas to fruition. Whether you require a bespoke wash-station, specialised instrument basket or customised transport system, our sales & engineering personnel will work with you to design a solution to fulfill your requirements. Contact us at [info@trinitymed.ie](mailto:info@trinitymed.ie).



**Batrik Soft Pockets**

Code	Descript.	Pack	Notes
18-CPF-02	Foam Corner Protectors - Designed to protect the corners of surgical trays from tearing wrap	80/bx	Single use - blue

**ChemEco Stainless Steel Cleaning Powder**

Code	Descript.	Pack	Notes
CECO0030	4 x 1Kg Tub	1 x 4	-

**Protective Transport Bags**


Code	Descript.	Pack	Notes
PPU150250	Transport Bag 150X250mm	500	-
PPU200300	Transport Bag 200X300mm	500	-
PPU300550	Transport Bag 300x550mm	500	-
PPU400550	Transport Bag 400x550mm	250	-
PPU400750	Transport Bag 400x750mm	200	-
PPU600750	Transport Bag 600x750mm	200	-
PPU600900	Transport Bag 600x900mm	150	-



**CSSD & Surgical Consumables**

Scan the QR code below to view the full range of CSSD & Surgical Consumable products from our partner Batrik.





Endo

## EndoPax 30 Day Endoscope Packing & Storage System



**The iM-MED EndoPax™** system provides clinical teams with a solution for the aseptic storage and safe transportation of endoscopes. This unique process safely contains the endoscopes for up to 30 days.

Following decontamination and compliant drying in an endoscope drying cabinet, thereby ensuring all residual bacteria are fully desiccated, an endoscope is packed using the space-saving, free standing iM EndoPax™ system.

With a rapid, user-friendly process, the system delivers a process ensuring the internal and external lumens are fully conditioned using an inert gas mix, to inhibit aerobic and anaerobic bacterial growth.

Without the addition of chemistry to the terminal packing process, potential risks are reduced to the endoscope, patient and user.

The medical device bag is heat sealed ensuring a reliable seal integrity until the user requires the endoscope.

The iM EndoPax™ system has an independent monitor to provide a fully validated and traceable process. Prolong the aseptic storage and safe transportation of endoscopes for up to 30 days.

iM EndoPax™ enables endoscopes to be prepared and be conveniently ready for routine lists, emergency situations, out of hours and off-site endoscopy work.

Endoscopes can be transported in a controlled and tamper evident environment using transportation cases or lockable trolleys.

Code	Descript.	Pack	Notes
600020	Endopax Solution with IMS	each	Contact Trinity Medical for more info.

**Ultrasonic Irrigation Wash Station - UltraZonic ENDO**



Code	Descript.	Pack	Notes
ENDO	Twin Endoscope Manual Wash System for automated pre-cleaning of Endoscopes. Height-adjustable.	each	Contact Trinity Medical for more info.

**Ultrasonic Irrigation Wash Station - UltraZonic ENDO Mini.**



Code	Descript.	Pack	Notes
ENDO MINI	MINI is ENDO's little brother, especially designed for the pre-cleaning of ENT scopes, cystoscopes and bronchoscopes. Like our acclaimed ENDO device, the ENDO MINI also covers all the steps in the pre-cleaning semi-automatically.	each	Contact Trinity Medical for more info.

## PAA Neutraliser



Code	Descript.	Pack	Notes
J94A	Neutralising Solution for Residual Peracetic Acid.	Box / 24	24 x 50ml bottles

## ChemDAQ Steri-Trac 2+ Hazardous Vapour Monitoring Solutions






Code	Descript.	Pack	Notes
n/a	Steri-Trac 2 Fixed, continuous gas monitoring system for measuring occupational exposure to Ethylene Oxide (EtO), Hydrogen Peroxide (H <sub>2</sub> O <sub>2</sub> ) and Peracetic Acid (PAA).	each	n/a

## Monmouth 'Circulaire' Self-Ventilating Chemical Storage Cabinet



Code	Descript.	Pack	Notes
n/a	Removes PAA vapours from ambient air. Available in a variety of sizes.	each	n/a

**Endo Tray Liners - iM CleanEst**

-  **Decontaminated Cover**
-  **Contaminated Cover**
-  **Clear Tray Liner**



Green tray cover for decontaminated endoscopes ready to use.




Red tray cover for contaminated endoscopes to be cleaned.



Code	Descript.	Pack	Notes
125010	iM CleanEST Liners. Sterile. Red, Green and Clear tray liner.	Box/200	n/a

**EndoTube - Disposable Tube Sets for Endoscopes**



Code	Descript.	Pack	Notes
Various	Available for Air/Water Bottle, Insufflation Line and Irrigation.	pack/20	Compatible with Olympus, Fuji and Pentax.

**Endoscopy Transport Trolley Systems**



Code	Descript.	Pack	Notes
Various	Available in various formats, open, closed, multi-level 6 & 10	each	Contact Trinity Medical for more info.



**Disposable Valve Sets**



Code	Descript.	Pack	Notes
Various	Disposable Valve Sets for <b>Olympus, Pentax &amp; Fujifilm</b> Endoscopes.	n/a	Contact Trinity Medical for more info.

**Biopsy Valve with Irrigation Tube**



Code	Descript.	Pack	Notes
830041	For Fujinon & Olympus	pack/400	Contact Trinity Medical for more info.

**Biopsy Valve - Various**



Code	Descript.	Pack	Notes
Various	Available for Olympus, Pentax & Fuji, Sterile & Non Sterile	each	Contact Trinity Medical for more info.

**Polyp Trap**



Code	Descript.	Pack	Notes
PT-01-05	Numbered collection chambers in a transparent transport container for accurate sample identification	10	5 Chambers



### Prefilled Bedside Cleaning Kit - 500ml Ready to Use



Code	Descript.	Pack	Notes
1033	Flexi-Clenz Prefilled Tub with Cylindrical Sponge	27	Contact Trinity Medical for more info.

### Endoscope Flushing & Leak Testing Device



Code	Descript.	Pack	Notes
SEB1000	Smart Endo Box - Flushing and Leak Testing Device	each	With digital interface. No internal tubing.

### Endoscope Flushing Device



Code	Descript.	Pack	Notes
Q88	Compact Device for Flushing of Endoscopes	each	With digital interface.



## Bite Block



Code	Descript.	Pack	Notes
BB-01-85-50	Without strap - Yellow	100	Opening: 27 x 23 mm
BB-01-85-50-T	With white textile strap - Yellow	100	Opening: 27 x 23 mm
BB-01-85-50-R	With rubber strap - Yellow	100	Opening: 27 x 23 mm

## Combination Cleaning Brushes



Code	Descript.	Pack	Notes
Various	Various Diameters Available	100	Contact Trinity Medical for more info.

## Cleaning Brush Sets



Code	Descript.	Pack	Notes
Various	Various Diameters Available	100	Contact Trinity Medical for more info.

## Cleaning Brushed



Code	Descript.	Pack	Notes
Various	Various Diameters Available	100	Contact Trinity Medical for more info.



## Endoscope Distal Protection



Code	Descript.	Pack	Notes
DTP-01	7 – 15 mm	100	50mm Length
DTP-02	2 – 10 mm	100	50mm Length

## 3-in-1 Distal Protection and Cleaning Sponge



Code	Descript.	Pack	Notes
EP-05-14-CL-EO	White - 126mm	50	Contact Trinity Medical for more info.

## Colonoscopy Dignity Shorts



Code	Descript.	Pack	Notes
DS-L-ECO	Dark Blue	100	Contact Trinity Medical for more info.
DS-L-ECO-XXL	Dark Blue	100	Contact Trinity Medical for more info.



**Manual Leak Test Device for Endoscopes**



Code	Descript.	Pack	Notes
200.00040.03	Manual Leak Tester	each	-

**Scope-Clenz Pouches**



Code	Descript.	Pack	Notes
1002P	Flushing Solution Specifically Formulated for use in the FlushTec / ScopeTech Flushing Devices only.	30x300ml	-

**Endo-Stack Equipment**



Code	Descript.	Pack	Notes
AQM	AQUA:MASTER Irrigation Device	each	Contact Trinity Medical for more info.
INSF	Insufflator 1800/16 C	each	Contact Trinity Medical for more info.
LM	Electronic Leak Test Device	each	Contact Trinity Medical for more info.



**Endotherapy Instrumentation**

Scan the QR code below to explore our full range of endotherapy products from our partner Body Products GMBH.





**SIUP**

ZAPARAY HLD3



**With the HLD-3 you're not just upgrading your disinfection method - you're transforming it!**

Replacing toxic biocidal manual procedures with a sustainable alternative. 100% validated by unique UV-C sensors and traceability features. The disinfection effectiveness is confirmed with each cycle. The device produces no heat and requires no toxic consumables. Being a compact device, the HLD-3 can reside next to your ultrasound workstation, avoiding costs with logistical damage.

The HLD-3 is a strategic investment for a safer future due to labour reduction, low energy consumption and long lifetime.

Using UV-C LED technology ensures effective bacteriological and viral reduction of minimal >99,99%. ZAPARAY is a respected market leader in UV-C LED disinfection technology. We develop genuine sustainable devices using scientific and academic partnerships.

With our automated disinfection, you reduce operating costs and staff load. Join us in creating safer, healthier environments with our unique UV-C LED disinfection for heat sensitive objects !

**USP's for our UV-C LED technology**



**ECO-FRIENDLY**

No chemicals, no water, a mercury-free solution.



**TIME EFFICIENT**

Easy to use. Guaranteed disinfection in less than 5 minutes.



**COST EFFICIENT**

One single radiation dose consumes only 10W of electricity.



**HEAT FREE ASSURANCE**

UV-C LED's work at room temperature.



**IMPROVED & DECENTRALIZED**

Manual disinfection is error prone, instrument logistics are a cross contamination risk.



**VALIDATED & AUTOMATED**

Prevent cross-contamination. Consistent, reproducible object disinfection.

### Customer Enquiry Form

Scan the QR code to complete a quick enquiry form, which will instantly notify the Trinity Med customer support team.





# **TRINITY** MEDICAL

PATIENT CENTRIC TECHNOLOGIES

[www.trinitymed.ie](http://www.trinitymed.ie)